

Product datasheet for **SC337155**

ACSBG2 (NM_001289177) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ACSBG2 (NM_001289177) Human Untagged Clone
Tag:	Tag Free
Symbol:	ACSBG2
Synonyms:	BGR; BRGL; PRTD-NY3; PRTDNY3
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC337155 representing NM_001289177.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGACTGGAACCCCAAGACTCAAGAAGGAGCTAAAGATCTTGAAGTAGACATGAATAAACAGAAAGTT
ACTCCCAGGCTGTGGACCACCTGTCGAGATGGAGAAGTCCTTCTGAGGCTATCCAAACACGGACCAGGC
CATGAGACCCCGATGACCATCCCTGAATTTTTTCGAGAGTCAGTCAACCGATTGGAACCTATCCAGCC
CTCGCATCCAAGAATGGCAAAAAGTGGGAAATTCTGAATTTCAACCAGTACTATGAGGCTGTGCGAAG
GCTGCAAAATCCTTGATCAAGCTGGGTTTGGAGCGTTTCCACGGAGTTGGTATCCTGGGGTTAACTCT
GCAGAGTGGTTTATCACTGCTGTTGGTGCCATCCTAGCCGGGGTCTTTGTGTTGGTATTTATGCCACC
AACTCTGCCGAGGTTTGTCAATATGTCATCACTCATGCCAAAGTGAACATCTTGCTGTTGAGAATGAT
CAACAGTTACAGAAAATCCTTTGATTCCACAGAGCAGCCTAGAGCCCCTAAAAGCGATCATCCAGTAC
AGACTGCCAATGAAGAAGAACAACAACCTGTACTCTTGGGATGATTTATGGAACCTGGCAGAAGTATC
CCTGACACCCAACTGGAGCAGGTATCGAGAGCCAGAAGGCGAATCAATGCGCAGTGCTCATCTACACT
TCAGGGACCACAGGCATACCCAAGGGAGTGATGCTCAGTCATGACAACATCACGTGGATTGCAGGAGCA
GTGACAAAGGACTTTAACTGACAGACAAGCATGAGACGGTGGTAGCTACCTCCCACTCAGCCATATT
GCAGCACAGATGATGGACATCTGGGTACCCATAAAGATTGGGGCGCTCACATACTTTGCTCAAGCAGAT
GCTCTCAAGGGCACCTTGGTAAGTACTCTAAAGGAGGTAAACCTACTGTCTTATTGGAGTGCCTCAA
ATTTGGGAGAAGATACATGAGATGGTGAAGAAAAATAGTGCCAAGTCCATGGGCTTGAAGAAGAAGGCA
TTCGTGTGGGCAAGAAACATTGGCTTCAAGGTCAACTCAAAAAGATGTTGGGGAATATAATACTCCC
GTGAGCTACCGCATGGCTAAGACTCTCGTGTTCAGCAAAGTCAAGACATCCCTTGGCTTGGATCACTGT
CACTCTTTTATCAGTGGGACTGCGCCCTCAACCAAGAGACTGCCGAGTTCTTTCTAAGCTTGGACATA
CCTATAGGCGAGTTGTATGGGTTGAGTGAGAGCTCGGGACCCACACGATATCCAACCAGAATAACTAC
AGGCTTCTAAGCTGTGGCAAGATCTTGACTGGGTGAAGAATATGCTGTTCCAGCAGAACAAGGATGGC
ATTGGGGAGATCTGCCTCTGGGGTAGGCACATCTTCATGGGCTATCTGGAAAGTGAGACTGAAACTACA
GAGGCCATCGATGATGAAGGCTGGCTACACTCTGGGGATCTGGGCCAGCTGGACGGTCTGGGTTTCTC
TATGTCACCGGCCACATCAAAGAAATCCTTATCACTGCTGGTGGTGAAAAATGTGCCCCCATTCTGTT
GAGACCTTGGTTAAGAAGAAGATCCCCATCATCAGTAACGCCATGTTAGTAGGAGATAAACTGAAGTTT
CTGAGCATGTTGCTGACGCTGAAGTGTGAGATGAATCAGATGAGCGGAGAACCTCTGGACAAGCTGAAC
TTCGAGGCCATCAACTTCTGTCGGGGTCTGGGCAGCCAGGCATCCACCGTACTGAGATTGTGAAGCAG
CAAGACCCCTGGTCTACAAGGCCATCCAGCAAGGCATCAATGCTGTGAACCAGGAAGCCATGAACAAT
GCACAGAGGATTGAAAAGTGGGTCTCTTGGAGAAGGACTTTTCCATCTATGGTGGAGAGCTAGGTCCA
ATGATGAAACTTAAGAGACATTTTGTAGCCAGAAATACAAAAACAAATTGATCACATGTACCACTGA
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001289177

Insert Size: 2001 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001289177.1</u>
RefSeq Size:	2798 bp
RefSeq ORF:	2001 bp
Locus ID:	81616
UniProt ID:	<u>Q5FVE4</u>
Cytogenetics:	19p13.3
MW:	74.4 kDa
Gene Summary:	Mediates activation of long-chain fatty acids for both synthesis of cellular lipids, and degradation via beta-oxidation. Able to activate long-chain fatty acids. Also able to activate very long-chain fatty acids; however, the relevance of such activity is unclear in vivo. Has increased ability to activate oleic and linoleic acid. May play a role in spermatogenesis. [UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a). Variants 1, 2, and 3 encode the same isoform.