

## Product datasheet for **SC337168**

### **RAD17 (NM\_001278622) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	RAD17 (NM_001278622) Human Untagged Clone
Tag:	Tag Free
Symbol:	RAD17
Synonyms:	CCYC; HRAD17; R24L; RAD17SP; RAD24
Vector:	pCMV6-Entry (PS100001)




[View online »](#)

Fully Sequenced ORF: >SC337168 representing NM\_001278622.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGAATCAGGTAACAGACTGGGTTGACCCATCATTTGATGATTTTCTAGAGTGTAGTGGCGTCTCTACT
ATTACTGCCACATCATTAGGTGTAATAACTCAAGTCATAGAAGAAAAATGGGCCTTCTACATTAGAA
AGCAGCAGATTTCCAGCGAGAAAAAGAGGAAATCTATCTTCCCTTAGAACAGATTTATGGTTTAGAAAA
TCAAAAGAAATATCTGTCTGAAAATGAACCATGGGTGGATAAATAAAACCAGAACTCAGCATGAACCT
GCTGTGCATAAAAAGAAAATTGAAGAAGTCGAAACCTGGTTAAAAGCTCAAGTTTTAGAAAAGCAACCA
AAACAGGGTGGATCTATTTTATAATAACAGGTCCTCCTGGATGTGGAAGACAACGACCTTAAAAATA
CTATCAAAGGAGCATGGTATTCAAGTACAAGAGTGGATTAATCCAGTTTTACCAGACTTCCAAAAAGAT
GATTTCAAGGGGATGTTAATACTGAATCAAGCTTCCATATGTTTCCCTATCAGTCTCAGATAGCAGTT
TTCAAAGAGTTTCTACTAAGAGCGACAAAGTATAACAAGTTACAATGCTTGGAGATGATCTGAGAACT
GATAAGAAGATAATTCTGGTTGAAGATTACCTAACCGTTTTATCGGGATTCTCATACTTTACATGAA
GTTCTAAGGAAGTATGTGAGGATTGGTCGATGTCCTCTTATATTTATAATCTCGGACAGTCTCAGTGG
GATAATAATCAAAGGTTATTGTTTCCCAAAGAAATTCAGGAAGAGTGTTCATCTCAAATATTAGTTTC
AACCTGTGGCACCAACAATTATGATGAAATTTCTTAATCGAATAGTGACTATAGAAGCTAACAAAGAT
GGAGGAAAAATTACTGTCCCTGACAAACTTCTCTAGAGTTGCTCTGTCAGGGATGTTCTGGTGATATC
AGAAGTGCAATAAACAGCCTCCAGTTTTCTTCTCAAAGGAGAAAAACACTTACGGCCAAGGAAAAAA
GGAATGTCTTTAAAAATCAGATGCTGTGCTGTCAAAATCAAACGAAGAAAAAACCTGATAGGGTTTTT
GAAAATCAAGAGGTCCAAGCTATTGGTGGCAAAGATGTTTCTCTGTTTCTCTTCAGAGCTTTGGGAAA
ATTCTATATTGTAAGAGCATCTTTAACAGAATTAGACTCACCTCGGTTGCCCTCTCATTTATCAGAA
TATGAACGGGATACATTACTTGTGAACTGAGGAGGTAGTAGAAATGTCACACATGCCTGGAGACTTA
TTAATTTATATCTTACCAAAACTACATAGATTTCTTCATGGAATTGATGATATTGTGAGAGCCAGT
GAATTTCTGAGTTTTGCAGATATCCTCAGTGGTACTGGAATACACGCTCTTTACTCAGGGAATATAGC
ACATCTATAGCTACGAGAGGTGTGATGCATTCCAACAAAGCCGAGGATATGCTCATTGCCAAGGAGGA
GGATCAAGTTTTCGACCCTTGCACAACTCAGTGGTTTCTAATAAATAAAAAGTATCGGGAAAAATTGC
CTGGCAGCAAAGCACTTTTCTGACTTCTGCCTACCAGCTTTATGCCTCCAACTCAGCTATTGCCA
TACCTTGCTCTACTAACCATTCCAATGAGAAATCAAGCTCAGATTTCTTTTATCCAAGATATTGGAAG
CTCCCTCTGAAGCGACACTTTGGAAGATTGAAAATGGAAGCCCTGACTGACAGGGAACATGGAATGATA
GACCCTGACAGCGGAGATGAAGCCAGCTTAATGGAGGACATTCTGCAGAGGAATCTCTGGGTGAACCC
ACTCAAGCCACTGTGCCGAACTGGTCTCTTCTTTGAGTCAAGATAGTGCCAGTGAAGTGCCTGCT
AGCCAGCCCCAGCCCTTTTCAGCCCAAGGAGACATGGAAGAAAACATAATAATAGAAGACTACGAGAGT
GATGGGACATAG
```

Restriction Sites: Sgfl-Mlul

**Plasmid Map:**


**ACCN:** NM\_001278622

**Insert Size:** 2013 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001278622.1](#)

**RefSeq Size:** 3067 bp

**RefSeq ORF:** 2013 bp

**Locus ID:** 5884

**UniProt ID:** [O75943](#)

**Cytogenetics:** 5q13.2

**Protein Families:** Druggable Genome

**MW:** 75.8 kDa

**Gene Summary:**

The protein encoded by this gene is highly similar to the gene product of *Schizosaccharomyces pombe rad17*, a cell cycle checkpoint gene required for cell cycle arrest and DNA damage repair in response to DNA damage. This protein shares strong similarity with DNA replication factor C (RFC), and can form a complex with RFCs. This protein binds to chromatin prior to DNA damage and is phosphorylated by the checkpoint kinase ATR following damage. This protein recruits the RAD1-RAD9-HUS1 checkpoint protein complex onto chromatin after DNA damage, which may be required for its phosphorylation. The phosphorylation of this protein is required for the DNA-damage-induced cell cycle G2 arrest, and is thought to be a critical early event during checkpoint signaling in DNA-damaged cells. Multiple alternatively spliced transcript variants of this gene, which encode four distinct protein isoforms, have been reported. Two pseudogenes, located on chromosomes 7 and 13, have been identified. [provided by RefSeq, Jul 2013]

Transcript Variant: This variant (9) has an alternate 5' terminal exon and lacks an internal exon in the 5' UTR, when compared to variant 1. The encoded protein (isoform 1) is identical to that encoded by variants 1, 5, 6, 7 and 8.