

Product datasheet for **SC337169**

GOLGA6L22 (NM_001271664) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GOLGA6L22 (NM_001271664) Human Untagged Clone
Tag:	Tag Free
Symbol:	GOLGA6L22
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC337169 representing NM_001271664.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGC
ATGCTGATGTGGCCCCAACCCACCTCCCTACCCATCCCCACCTCCCTACCCATCCCCACCTCCCTACC
CATCCCCACCTCCCTACCCATCCCCACCTCCCTACCCATCTATGATGTCCAAGAAACAGACAGAGC
AAATTGGCCGAGGCCAAGGAACAGTTGACAGACCATCATCCCCAGACCAACCCCTAGTGTTGGTACAGCA
GCAAGCGACACCAAAAAGAAGAAATAAATAATGGCACTAACCCCTGAGACAACCACTTCTGGTGGTTGC
CACTCGCCTGAGGATGAACAGAAGGCAAGCCACCAACATCAGGAAGCCCTAAGGAGGGAGCTAGAGGCC
CAGGTTTCATACCATACGAATCCTTACATGTCAGAAAAGTCTGAGTTTCTGAGTGGCACTCTACTACAGCCAG
CATGCTGTCAAGCAGTTGGAAGGAGAGGCCAGGGATCTGATCAGCCGCCTGCATGATTTCATGAAGTTT
GCAGGAGAGTTAGAGCAGGCTCTCTGCTGTCTGCTACACAGAAGAAGAAGGCGGATAGGTACATTGAG
GAGTTAACAAAGGAGAGGGACGCCCTGAGTCTGGAAGTGTACAGGAACACCATAACTGATGAGGAGCTG
AAGGAGAAAAATGCCAACTACAAGAAAACTTCAACTTGTAGAATCTGAAAAGTCTGAGATCCAGCTC
AACGTAAGGAGCTAAAAAGGAACTGGAGAGGGCCAAGCTCCTGCTGCCACAGCAGCAGCTGCAGGCG
GAGGCTGACCACTGGGTAAAGGAGCTGCAGAGTGTGTAGCAAGCTCCAGCCAGGTGGAAGAGAAC
GAGTTGTGGAACCGCCTGAACAGCAACAGGAGGAGAAGATGTGGAGGCAGGAGGAGAAGATACAGGAG
TGGGAGGAGAAGATACAGGAGCAGGAGGAGAAGATACGGGAGCAGGAGGAGAAGATACGGGAGCAGGAG
GAGAAGATGCGGAGGCAGGAGGAGATGATGTGGGAGAAGGAGGAGAAGATGCGGAGGCAGGAGGAGATG
ATGTGGGAGAAGGAGGAGAAGATGCGGAGGCAGGAGGAGATGATGTGGGAGAAGGAGGAGAAGATGCGG
AGGCTGGAGGAGATGATGTGGGAGAAGGAGGAGAAGATACGGGAGCTGGAAGAGAAGATGCACGAGCAG
GAGAAGATACGGGAGCAGGAAGAGAAGAGGCAGGAGGAGGAGAAGATACGCGAGCAGGAGAAGAGGCAG
GAGCAGGAGGCGAAGATGTGGAGGCAGGAGGAGAAGATACGGGAGCAGGAAGAGAAGATACGGGAGCAG
GAGAAAAAGATGTGGAGGCAGGAGGAGAAGATTACGAGCAGGAGAAGATACGGGAGGAGGAGAAGAGG
CAGGAGCAGGAGGAGATGTGGAGGCAGGAGGAGAAGATAAGGGAGCAGGAGGAGATATGGAGGCAAAAG
GAGAAGATGCACGAGCAGGAGAAGATACGGAAGCAGGAGGAGAAGGTGTGGAGGCAGGAGGAGAAGATG
CACGACCAGGAGGAGAAGATACGGGAGCAGGAGGAGAAGATGTGGAGGCAGGAGGAGAAGATAAGGGAG
CAGGAGGAGAAGATACGGGAGCAGGAGGAGAAGATACGAGAGCAGGAGGAGATGATGCAGGAACAGGAA
GAGAAGATGGGGAGCAGGAAGAGAAGATGCAAGAACAGGAGAAGATGCGGAGGCAGGAGGAGAAGATA
AGGGAGCAGGAGGAGAAGATACGGGAGCAGGAAGAGAAGATACGAGAGCAGGAGGAGAAGATATGGGAG
CAGGAGGAGAAGATACGAGAGCAGGAGGAGATGATGCAGGAACAGGAAGAGAAGATGTGGGAGCAGGAG
GAGAAGATGTGTGAGCAGGAAGAGAAGATGCAAGAACAGGAGGAGAAGATGCGGAGGCAGGAGGAGAAG
ATGTGGGAGCAGGAAGTGAAGTGTGCGGCAGCAGGAGGAGAAGATGCAGGAACACCAGGAGCACCTGGAA
GCTGCCATCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001271664

Insert Size: 2013 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001271664.2](#)

RefSeq Size: 3806 bp

RefSeq ORF: 2013 bp

Locus ID: 440243

Cytogenetics: 15q11.2

MW: 83.6 kDa