

Product datasheet for **SC337176**

Uromuroid (UMOD) (NM_001278614) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Uromuroid (UMOD) (NM_001278614) Human Untagged Clone
Tag:	Tag Free
Symbol:	Uromuroid
Synonyms:	ADMCKD2; ADTKD1; FJHN; HNFJ; HNFJ1; MCKD2; THGP; THP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC337176 representing NM_001278614.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGGCAGCCATCTCTGACTTGGATGCTGATGGTGGTGGCCTCTTGTTTCATCACAACGTCAGCC
ACTGACACCTCAGAAGCAAGAACAAGTACAACAGCCAGCCAGGTGGAGTTGGAGGACTCCGAGAGA
GGAGGCGACACTGAGCAGGGTCTGATGAAGATTTACCAGCCAGGGAAGATGGTCTCTGAATGTCAC
AGCAATGCCACCTGCACGGAGGATGAGGCCGTTACGACGTGCACCTGTCAGGAGGGCTTACCAGCGAT
GGCCTGACCTGCGTGGACCTGGATGAGTGCGCCATTCTGGAGCTCACAACGTCCTCCGCAACAGCAGC
TGCGTAAACACGCCAGGCTCCTTCTCTGCGTCTGCCCGAAGGCTTCCGCTGTGCGCCGGTCTCGGC
TGCACAGACGTGGATGAGTGCCTGAGCCTGGGCTTAGCCACTGCCACGCCCTGGCCACATGTGTCAAT
GTGGTGGCAGCTACTTGTGCGTATGCCCGCGGGCTACCGGGGGATGGATGGCACTGTGAGTGTCC
CCGGGCTCTGCGGGCCGGGTTGGACTGCGTGCCCGAGGGCGACGCGCTCGTGTGCGGGATCCGTGC
CAGGCGCACCGCACCTGGACGAGTACTGGCGCAGCACCAGTACGGGGAGGGCTACGCCGCGACACG
GACCTGCCGGCTGGTACCGCTTCGTGGGCCAGGGCGGTGCGCGCATGGCCGAGACTGCGTGCCAGTC
CTGCGCTGCAACACGGCCGCCCCATGTGGCTCAATGGCACGCATCCGTCCAGCGACGAGGGCATCGTG
AGCCGCAAGGCTGCGCGCACTGGAGCGGCCACTGCTGCCTGTGGGATGCGTCCGTCCAGGTGAAGGCC
TGTGCCGGCGGCTACTACGTCTACAACCTGACAGCGCCCCGAGTGTCACTGGCGTACTGCACAGAC
CCCAGTCCGTGGAGGGGACGTGTGAGGAGTGCAGTATAGACGAGGACTGCAAAATCGAATAATGGCAGA
TGGCACTGCCAGTGAACACGAGCTTCAACATCACTGATATCTCCCTCTGGAGCACAGGCTGGAATGT
GGGGCAATGACATGAAGGTGTCGCTGGGCAAGTGCCAGCTGAAGAGTCTGGGCTTCGACAAGGTCTTC
ATGTACTGAGTGACAGCCGGTGTCTGGGCTTCAATGACAGAGACAACCGGGACTGGGTGTCTGTAGTG
ACCCAGCCCGGGATGGCCCTGTGGGACAGTGTGACGAGGAATGAAACCCATGCCACTTACAGCAAC
ACCCTCTACCTGGCAGATGAGATCATATCCGTGACCTCAACATCAAAATCAACTTTGCATGCTCCTAC
CCCCTGGACATGAAAGTCAGCCTGAAGACCGCCCTACAGCCAATGGTCAGTGTCTAAACATCAGAGTG
GGCGGGACCGCATGTTACCGTGCAGTGGCGCTCTCCAGACCCCTTCTACACGCAGCCCTACCAA
GGCTCCTCCGTGACACTGTCCACTGAGGCTTTTCTCTACGTGGGACCATGTTGGATGGGGGCGACCTG
TCCCGATTTGACTGCTCATGACCAACTGCTATGCCACACCCAGTAGCAATGCCACGACCCCTGAAG
TACTTCATCATCCAGGACAGATGCCACACACTAGAGACTCAACTATCCAAGTGGTGGAGAATGGGGAG
TCCTCCAGGGCCGATTTCCGTCCAGATGTTCCGGTTTGTGAAACTATGACCTAGTCTACCTGCAC
TGTGAAGTCTATCTGTGACACCATGAATGAAAAGTGAAGCCTACCTGCTCTGGGACAGATTCCGA
AGTGGGAGTGTCATAGATCAATCCCGTGTCTGAACTGGGTCCCATCACACGAAAGGTGTCCAGGCC
ACAGTCTCAAGGGCTTTTAGCAGCTTGGGGCTCCTGAAAGTCTGGCTGCCTCTGCTTCTCTCGGCCAC
TTGACCCTGACTTTTCAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001278614

Insert Size: 2022 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278614.1](#)

RefSeq Size: 2589 bp

RefSeq ORF: 2022 bp

Locus ID: 7369

UniProt ID: [P07911](#)

Cytogenetics: 16p12.3

MW: 73.6 kDa

Gene Summary: The protein encoded by this gene is the most abundant protein in mammalian urine under physiological conditions. Its excretion in urine follows proteolytic cleavage of the ectodomain of its glycosyl phosphatidylinositol-anchored counterpart that is situated on the luminal cell surface of the loop of Henle. This protein may act as a constitutive inhibitor of calcium crystallization in renal fluids. Excretion of this protein in urine may provide defense against urinary tract infections caused by uropathogenic bacteria. Defects in this gene are associated with the renal disorders medullary cystic kidney disease-2 (MCKD2), glomerulocystic kidney disease with hyperuricemia and isosthenuria (GCKDHI), and familial juvenile hyperuricemic nephropathy (FJHN). Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jul 2013]

Transcript Variant: This variant (3) includes an additional in-frame exon in the 5' coding region, compared to variant 1, resulting in an isoform (b) that is longer than isoform a.