

Product datasheet for **SC337179**

FIGNL1 (NM_001287494) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: FIGNL1 (NM_001287494) Human Untagged Clone
Tag: Tag Free
Symbol: FIGNL1
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC337179 representing NM_001287494.
Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
ATGCAGACCTCCAGCTCTAGATCTGTGCACCTGAGTGAATGGCAGAAGAATTACTTCGCAATTACATCT
GGCATATGTACCGGACCGAAGGCAGATGCATACCGTGCACAGATATTACGCATTCAAGTATGCATGGGCA
AACTCTGAGATTTCCAGGTCTGTGCTACCAAAGTGTTCAAAAAATATGCAGAGAAATATTCTGCAATT
ATTGATTCTGACAATGTTGAATCTGGGTGAATAATTATGCAGAAAACATTTTAACTTTGGCAGGATCT
CAACAAACAGATAGTGACAAGTGGCAGTCTGGATTGTCAATAAATATGTTTTCAAAATGAGTAGTGTA
CAGAAGATGATGCAAGCTGGCAAAAAATTCAAAGACTCTCTGTTGGAACCTGCTCTGCATCAGTGGTA
ATCCATAAGGAGGCCACTGTCTTTGATCTTCTAAATTTAGTGTTTGTGGTAGTTCTCAAGAGAGTGAC
TCATTACCTAACTCAGCTCATGATCGAGACCGGACCCAAGACTTCCCGGAGAGCAATCGTTTGAAACTC
CTTCAGAAATGCCAGCCACCTATGGTGACTAACACTGCTAGGACTTGTCTACATTCTCAGCACCTGTA
GGTGAGTCAGCTACTGCAAAATCCATGTCACACCATTGTTTGAAAATGTCAAAAAGGAAAATCACAGC
TCTGCAAAAGAAAACATAGGACTTAATGTGTTCTTATCTAACCAGTCTTGTTCCTGCTGCCTGTGAA
AATCCACAGAGGAAGTCTTTTTATGGTTCTGGCACCATTGATGCACCTTCCAATCCAATACTGAATAAG
GCTTGTAGTAAAACAGAAGATAATGGCCAAAGGAGGATAGCAGCCTGCCTACATTTAAACTGCAAAA
GAACAATTATGGGTAGATCAGCAAAAAAGTACCACCAACCTCAGCGTGCATCAGGGTCTTCATATGGT
GGTGTAATAAAGTCTCTAGGAGCTAGTAGATCCCAGGGATACTTGGAAGTTTGTTCCTCTATACCC
AAGCAAGATGGGGGAGAGCAGAATGGAGGAATGCAATGTAAGCCTTATGGGGCAGGACCTACAGAACCA
GCATCCAGTTGATGAGCGTCTGAAGAACTTGGAGCCAAAGATGATTGAACTTATTATGAATGAGATT
ATGGATCATGGACCTCCAGTAAATTGGGAAGATATTGCAGGAGTAGAATTTGCTAAAGCCACCATAAAG
GAAATAGTTGTGTGGCCATGTTGAGGCCAGACATCTTTACTGGTTTAAAGGGACCCCTAAAGGAATT
TTGCTCTTTGGTCTCTGGGACTGGTAAACTCTAATGGCAAGTGCATTGCTAGTCAGTCTGGGGCA
ACATTTCTTAGCATCTCTGCTTCATCCTTAACTTCTAAATGGGTAGGTGAGGGGAGAAAATGGTCCGT
GCATTTGTTGCTGTTGCAAGGTGTCAGCAACCAGCTGTGATATTTATTGACGAAATTTGATTCCTTGTTA
TCTCAACGGGGAGATGGTGAGCATGAATCTTCTAGAAGGATAAAAACAGAATTTTATGTTCAATTAGAT
GGAGCAACAACATCTTCTGAAGATCGTATCCTAGTGGTGGGAGCAACAAATCGGCCACAAGAAATGAT
GAGGCTGCCCGGAGAAGATTGGTGAAAAGGCTTTATATCCCCTCCCAGAAGCTTCAGCCAGGAAACAG
ATAGTAATTAATCTAATGTCCAAGAGCAGTGTTCCTCAGTGAAGAAGAAATGAACAGATTGTACAG
CAGTCTGATGCGTTTTAGGAGCAGACATGACACAGCTTTCAGGGAGGCTTCTCTTGGTCTATTTCGC
AGTTTACAACTGCTGACATTGCTACCATAACACCGGATCAAGTTCGACCCATAGCTTACATTGATTTT
GAAAATGCTTTAGAACTGTGCGACCTAGTGTTCCTCAAAAGATTTAGAGCTTTATGAAAAGTGAAC
AAAATTTTGGTTGTGGAAGTAA
```



[View online »](#)

Restriction Sites: Sgfl-MluI

Plasmid Map:



ACCN: NM_001287494

Insert Size: 2025 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001287494.1](#)

RefSeq Size: 3381 bp

RefSeq ORF: 2025 bp

Locus ID: 63979

UniProt ID: [Q6PIW4](#)

Cytogenetics: 7p12.2

MW: 74.1 kDa

Gene Summary:

This gene encodes a member of the AAA ATPase family of proteins. The encoded protein is recruited to sites of DNA damage where it plays a role in DNA double-strand break repair via homologous recombination. This protein has also been shown to localize to the centrosome and inhibit ciliogenesis, and may regulate the proliferation and differentiation of osteoblasts. [provided by RefSeq, Oct 2016]

Transcript Variant: This variant (5) differs in the 5' UTR compared to variant 3. Variants 1-6 and 10-15 encode the same isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.