

## Product datasheet for **SC337219**

### ITGB6 (NM\_001282354) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ITGB6 (NM_001282354) Human Untagged Clone
Tag:	Tag Free
Symbol:	ITGB6
Synonyms:	A11H
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_001282354, the custom clone sequence may differ by one or more nucleotides

```

ATGGGGATTGAACTGCTTTGCCTGTTCTTTCTATTTCTAGGAAGGAATGATCACGTACAAGGTGGTGGCG
AGACTCTGCAGGTGCATGTCCGCCAGACTGAGGACTACCCGGTGGATTTGTATTACCTCATGGACCTCTC
CGCTCCATGGATGACGACCTCAACACAATAAAGGAGCTGGGCTCCCGGCTTTCCAAAGAGATGTCTAAA
TTAACCAGCAACTTTAGACTGGGCTTCGGATCTTTTGTGGAAAAACCTGTATCCCCTTTTGTGAAAAACA
CACCAGAGAAATTGCCAACCTTGCAAGTATTCCATACTTCTGTTTACCTACATTTGGATTCAAGCA
CATTTTGCATTGACAAAATGATGCTGAAAGATTCAATGAAATTGTGAAGAATCAGAAAAATTTCTGCTAAT
ATTGACACACCCGAAGGTGGATTTGATGCAATTATGCAAGCTGCTGTGTGAAGGAAAAAATTTGGCTGGC
GGAATGACTCCCTCCACCTCCTGGTCTTTGTGAGTGATGCTGATTCTCATTTTGGAAATGGACAGCAAAC
AGCAGGCATCGTCATTCTAATGACGGCTCTGCACTGGACAGCAAGAATGAATACTCCATGTCAACT
GTCTTGGAAATATCCAACAATTGGACAACCTATTGATAAACTGGTACAAAACAACGTGTTATTGATCTTCG
CTGTAACCAAGAACAAGTTTATATGAGAATTACGCAAACTTATCCTGGAGCTACAGTAGGTCT
ACTTCAGAAGGACTCCGGAACATTCTCCAGCTGATCATCTCAGCTTATGAAGAAGTGGCGTCTGAGGTG
GAACTGGAAGTATTAGGAGACACTGAAGGACTCAACTTGTCAATTTACAGCCATCTGTAACAACGGTACCC
TCTTCCAACACAAAAGAAATGCTCTCACATGAAAGTGGGAGACACAGCTTCTTCAGCGTGACTGTGAA
TATCCCACTGCGAGAGAAGAAGCAGGCACATTATCATAAAGCCTGTGGGGCTGGGGATGCCCTGGAA
TTACTTGTGAGCCAGAAATGCAACTGCGACTGTCAGAAAGAAGTGAAGTGAACAGCTCCAAATGTCACC
ACGGGAACGGCTCTTCCAGTGTGGGTGTGTGCCTGCCACCTGGCCACATGGGGCTCGCTGTGAGTG
TGGCGAGGACATGCTGAGCACAGATTCTGCAAGGAGGCCAGATCATCCCTCTGCAGCGGAAGGGGT
GACTGCTACTGTGGCAGTGTATCTGCCACTTGTCTCCATGGAACATTTATGGGCCTATTGGCAGT
GTGACAATTTCTCCTGCGTGAGACACAAAAGGGCTGCTCTGCGGAGGTAACGGCGACTGTGACTGTGGTGA
ATGTGTGTGCAGGAGCGGCTGGACTGGCGAGTACTGCAACTGCACCACCAGCACGGACTCCTGCGTCTCT
GAAGATGGAGTGTCTGCAGCGGGCGGGGACTGTGTTTGTGGCAAGTGTGTTTGCACAAAACCTGGAG
CCTCAGGACCAACCTGTGAACGATGTCCTACCTGTGGTACCCTGTAACCTAAACGGAGCTGCATTGA
GTGCCACCTGTCAGCAGCTGGCCAAGCCCGAGAAGAATGTGTGGACAAGTCAAAGTGTGGTGGCAGC
ATCAGTGAAGAAGAAGATTTCTCAAAGGATGGTTCTGTTTCTGCTCTCTGCAAGGAGAAAAATGAATGTC
TTATTACATTCTAATAACTACAGATAATGAGGGGAAAACCATCATTACAGCATCAATGAAAAAGATTG
TCCGAAGCCTCCAAACATTCCCATGATCATGTTAGGGGTTTCCCTGGCTATTCTTCTCATCGGGTTGTC
CTACTGTGCATCTGGAAGCTACTGGTGTCAATTCATGATCGTAAAGAAGTTGCCAAATTTGAAGCAGAAC
GATCAAAAGCCAAGTGGCAAACGGGAACCAATCCACTCTACAGAGGATCCACAAGTACTTTAAAAATGT
AACTTATAAACACAGGAAAAACAAAAGGTAGACCTTTCCACAGATTGCTAG
    
```

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_001282354

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001282354.1, NP_001269283.1</u>
<b>RefSeq Size:</b>	4408 bp
<b>RefSeq ORF:</b>	2082 bp
<b>Locus ID:</b>	3694
<b>Cytogenetics:</b>	2q24.2
<b>Protein Families:</b>	Druggable Genome, Transmembrane
<b>Protein Pathways:</b>	Arrhythmogenic right ventricular cardiomyopathy (ARVC), Dilated cardiomyopathy, ECM-receptor interaction, Focal adhesion, Hypertrophic cardiomyopathy (HCM), Regulation of actin cytoskeleton
<b>Gene Summary:</b>	<p>This gene encodes a protein that is a member of the integrin superfamily. Members of this family are adhesion receptors that function in signaling from the extracellular matrix to the cell. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. The encoded protein forms a dimer with an alpha v chain and this heterodimer can bind to ligands like fibronectin and transforming growth factor beta 1. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]</p> <p>Transcript Variant: This variant (3) lacks two in-frame exons in the 5' coding region compared to variant 1. The encoded isoform (b) is shorter compared to isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>