

Product datasheet for **SC337255**

FSD2 (NM_001281805) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FSD2 (NM_001281805) Human Untagged Clone
Tag:	Tag Free
Symbol:	FSD2
Synonyms:	SPRYD1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001281805, the custom clone sequence may differ by one or more nucleotides

```

ATGGAGGAGGAATCGGGGAGGAAGTGGGCTGGACAGGTCTACTCCAAGGATTTCCACTTTTACCACA
TGGACCTGTATGACTCTGAAGACAGACTGCACCTCTTTCCAGAAGAGAACTAGGATGAGGAAAGTAGT
CCAAGCTGAAATGGCCAATGAGTCAAGAGGAGCAGGGGATGGTAAGGCTCAAAGAGACCTTCAAGAGGAA
GTGGATGAACCTTGCCACTTATATGGACTTGAAGATGATCATGAATTAGGGGATGAGTTTGTGATGAAA
ACATACCCAGAACAGGGGTTTCAGAATATCCTCCTTATATGATGAAGAGGAGAGACCCAGCCAGGGAGCA
GAGAGACTGGAGACTTAGTGGAGAGGCAGCGGAGGCCGAGGACCTGGGCTTCGGAGGGTGGGGCTCAGCA
GGCCAGTGCCAGGACTTGCGGGAAGCCTATAGGTACACACACGGCCGTGCCAGCGAGGAGTATGAATGCT
ATGTCATCCCCGAGGAAGAGGATGAAGAAGAAGCTGCTGATGTCTTCTGTGCTCACTTGTAAAGACTCCAAT
AAGAGCTTTTTCAGAAGGTTTTTGTGAACATAAGGAGCATGAAGTATCCCACTCAATGAAGCACTGGAA
AGTGCCAAGGATGAAATTCACAAAAACATGTACAAATTGAAAAAGCAGATAATTGAGATGGAAAACTTTG
CAAATCATTTGGAGGAGTTTTTCATCAGTGTGGAGGAGAAATTTTGAAAAACAAGAACAACTTTGAGTC
ACATTACAACGAGATCTTGGAAACACTTGTCTAAAAATATGAGGAAAAAATACAAGCTCTAGGGGAGAAA
AAGAAAGAGAAGCTGGAAGCCTTATATGGACAAGTGGTCAGCTGTGGAGAAAACTAGATACTTGC AAAAG
AACTGATGGAAACAATAGAGGAGATGTGTACGAGGAGAAGGTGGATTTCATAAAGGATGCAGTGCTAT
GGCTGACAGACTCGGGAAATTCCTGAAAACAAAAACAGACGTGGAATCTCTGCACAGCCTGAGTTTGAA
GACCAGACCTTGATTTCTCTGATGTGGAGCAGCTGATGGGCTCCATTAAACACATTCCAGCTCCTTCTG
CTCCAGTCATAAATCCACAGGTCCTAAGTACGCCACAGGTTCTCAGTCCGAGTGTGCTGGAGCTTATA
CTCCGACGACACTGTGGAGAGCTATCAGCTGTCTACCGCCAGTGCAGGACAGCTCACCTGGGACGGAC
CAAGCAGCACCTTCTCCCCCATTATTTAAACCAAGAGATAAGGAGCTGTGAAGAGGCTGTGCTGATTT
GCTGGGAGTCTGGGAACCTGAATCCTGTGGACTGTACACTGTGGAGCTGACCCAGGCTGAAAGTCCAGA
AGCCTCGGGTGTAACTGAGTCTGTTGTGGGCATCCCGACCTGCGAGTCCGTGGTGCAGCTGCAGCCGGGG
CGGAGCTACATTATCTATGTGCGAGCCCTCAATATGGGGGGCCCCAGCGTGAGGAGCGAGCCAGCTACAG
TCCACACCATAGGAAGCTACTTTCGCCTAAACAAGGACACTTGCCATCCCTGGCTGACCATTTCTGAAGA
CGGACTTACGGCTGTACGAAGTGAAGGAGAACCCCCGCCAGGGAGCTGTACCCAGCGACACTCACTTC
ACCAGGTGTGTTGCTGTCATGGGAAATCTAATTCCAGTCCGAGGGCACCATTACTGGGAGGTGGAGGTGG
ATGAGCATTTGGACTACAGAGTTGGTGTGGCCTTTGCAGATGTCCGTAACAGGAGGATCTGGGAGCAAA
TTGCCTCTCCTGGTGCATGAGGCACACATTTGCATCATCAAGGCATAAGTATGAATTTCTGCACAACAGA
ACAACTCCAGATATAAGAATAACAGTTCACCAAGAAGATTGGCATTCTATTAGACTATGAACATTCAA
AGTTGTCAATTTTCAATGTGGACCTTTCTCAGCATCTATATACATTTAGTTGTGCTGCTTACGAATTTGT
GCATCCCTGTTTTTCTTTGAAAAGCCTGGGTGTCTAAAGGTACATAATGGCATTTCATGCCAAAACAT
GTCACCTTCTATAG
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001281805

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001281805.1</u> , <u>NP_001268734.1</u>
RefSeq Size:	6418 bp
RefSeq ORF:	2115 bp
Locus ID:	123722
UniProt ID:	<u>A1L4K1</u>
Cytogenetics:	15q25.2
Protein Families:	Druggable Genome
Gene Summary:	This gene encodes a protein that belongs to the FN3/SPRY family of proteins. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013] Transcript Variant: This variant (2) contains an alternate exon in the 5' UTR and lacks an alternate in-frame exon compared to variant 1. The encoded isoform (b) is shorter than isoform a. Variants 2 and 3 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.