

Product datasheet for **SC337317**

DRP1 (DNM1L) (NM_001278463) Human Untagged Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | DRP1 (DNM1L) (NM_001278463) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | DNM1L |
| Synonyms: | DLP1; DRP1; DVLP; DYMPLE; EMPF; EMPF1; HDYNIV; OPA5 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001278463, the custom clone sequence may differ by one or more nucleotides

```

ATGGAGGCGCTAATTCCTGTCATAAACAAGCTCCAGGACGCTTCAACACGGTGGGCGCCGACATCATCC
AGCTGCCTCAAATCGTCGTAGTGGGAACGACAGCAGCGGAAAGAGCTCAGTGCTAGAAAGCCTGGTGGG
GAGGGACCTGCTCCAGAGGACTGGAATTGTCACCCGGAGACCTCATTCTGCAACTGGTCCATGTT
TCACAAGAAGATAAACGGAAAACAACAGGAGAAGAAAATGGGGTGAAGCAGAAGAAATGGGGTAAATTTT
TTCACACCAAAAATAAGCTTTACACGGATTTTGATGAAATTCGACAAGAAATGAAAATGAAACAGAAAG
AATTTTCAGGAAATAATAAGGGAGTAAGCCCTGAACCAATTCATCTTAAGATTTTTTCACCCAACGTTGTC
AATTTGACACTTGTGGATTTGCCAGGAATGACCAAGGTGCCTGTAGGTGATCAACCTAAGGATATTGAGC
TTCAAATCAGAGAGCTCATTCTTCGTTTACAGTAATCCTAATTCATTATCCTCGCTGCTACTGCTGC
TAATACAGATATGGCAACATCAGAGGCACTTAAAATTTCAAGAGAGGTAGATCCAGATGGTCGCAGAACC
CTAGCTGTAATCACTAACTTGATCTCATGGATGCGGGTACTGATGCCATGGATGATTGATGGGAAGGG
TTATTCCAGTCAAACCTTGAATAATTGGAGTAGTTAACAGGAGCCAGCTAGATATTAACAACAAGAAGAG
TGTAACCTGATTCAATCCGTGATGAGTATGCTTTTCTTCAAAGAAATATCCATCTCTGGCCAATAGAAAT
GGAACAAAGTATCTTGCTAGGACTCTAAACAGGTTACTGATGCATCACATCAGAGATTGTTTACCAGAGT
TGAAAACAAGAATAAATGTTCTAGCTGCTCAGTATCAGTCTTTCTAAATAGCTACGGTGAACCCGTGGA
TGATAAAAGTGCTACTTTACTCCAATTATTACCAATTTGCCACAGAATATTGTAACACTATTGAAGGA
ACTGCAAAATATATTGAACTTCGGAGCTATGCGGTGGTCTAGAATTTGTTATATTTTCCATGAGACTT
TTGGGCGAACCTTAGAATCTGTTGATCCACTTGGTGGCCTAACACTATTGACATTTTGACTGCCATTAG
AAATGCTACTGGTCCCTCGTCTGCTTTATTTGTGCCTGAGGTTTCATTTGAGTTACTGGTGAAGCGGCAA
ATCAAACGTCTAGAAGAGCCAGCCTCCGCTGTGTGGAAGTGGTTCATGAGGAAATGCAAAGGATCATT
AGCACTGTAGCAATTACAGTACACAGGAATTGTTACGATTTCTAACTTCATGATGCCATAGTTGAAGT
GGTGACTTGTCTTCTCGTAAAAGGTTGCCTGTTACAAATGAAATGGTCCATAACTTAGTGGCAATTGAA
CTGGCTTATATCAACACAAAACATCCAGACTTTGCTGATGCTTGTGGGCTAATGAACAATAATATAGAGG
AACAAAGGAGAAACAGGCTAGCCAGAGAATTACCTTCAGCTGTATCACGAGACAAGTCTTCTAAAGTTCC
AAGTGCTTTGGCACCTGCCTCCAGGAGCCCTCCCCGCTGCTTCTGCTGAGGCTGATGGCAAGGTTGCA
TCTGGAGTGGTGGGTTGGAGATGGTGTCAAGAACCAACCACAGGCAACTGGAGAGGAATGCTGAAAA
CTTCAAAGCTGAAGAGTTATTAGCAGAAGAAAAATCAAACCCATTCCAATTATGCCAGCCAGTCCACA
AAAAGGTCATGCCGTGAACCTGCTAGATGTGCCAGTTCCTGTTGCACGAAAATCTGCTCGGGAACAG
CGAGATTGTGAGGTTATTGAACGACTCATTAAATCATATTTTCTATTGTCAGAAAGAATATTCAAGACA
GTGTGCCAAAGGCAGTAATGCATTTTTTGGTAAATCATGTGAAAGACACTCTTCAGAGTGAGCTAGTAGG
CCAGCTGTATAAATCATCCTTATTGGATGATCTTCTGACAGAATCTGAGGACATGGCACAGCGCAGGAAA
GAAGCAGCTGATATGCTAAAGGCATTACAAGGAGCCAGTCAAATATTGCTGAAATCCGGGAGACTCATC
TTTGGTGA
    
```

- Restriction Sites:** SgfI-MluI
- ACCN:** NM_001278463
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278463.1](#), [NP_001265392.1](#)

RefSeq Size: 4582 bp

RefSeq ORF: 2178 bp

Locus ID: 10059

UniProt ID: [O00429](#)

Cytogenetics: 12p11.21

Protein Pathways: Endocytosis, Fc gamma R-mediated phagocytosis

Gene Summary: This gene encodes a member of the dynamin superfamily of GTPases. The encoded protein mediates mitochondrial and peroxisomal division, and is involved in developmentally regulated apoptosis and programmed necrosis. Dysfunction of this gene is implicated in several neurological disorders, including Alzheimer's disease. Mutations in this gene are associated with the autosomal dominant disorder, encephalopathy, lethal, due to defective mitochondrial and peroxisomal fission (EMPF). Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jun 2013]
Transcript Variant: This variant (4) lacks an in-frame exon in the 3' coding region, compared to variant 1. The encoded isoform (4) is shorter, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.