

Product datasheet for **SC337369**

PSMD2 (NM_001278709) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PSMD2 (NM_001278709) Human Untagged Clone
Tag:	Tag Free
Symbol:	PSMD2
Synonyms:	P97; RPN1; S2; TRAP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001278709, the custom clone sequence may differ by one or more nucleotides

```

ATGAGTATGT CAGGAGAAGTGGCTAAGGAGTGGCAGGAGCTGGATGACGCAGAGAAGGTCACGGGGAGC
CTCTGCTCACTCTGGTGAAGGAAATCGTCCCCTATAACATGGCCCACAATGCAGAGCATGAGGCTTGCGA
CCTGCTTATGGAATTGAGCAGGTGGACATGCTGGAGAAGGACATTGATGAAAATGCATATGCAAAGGTC
TGCCTTTATCTACCAGTTGTGTGAATTACGTGCCTGAGCCTGAGAACTCAGCCCTACTGCGTTGTGCC
TGGGTGTGTTCCGAAAGTTAGCCGCTTCCCTGAAGCTCTGAGATTGGCATTGATGCTCAATGACATGGA
GTTGGTAGAAGACATCTTACCTCTGCAAGGATGTGGTAGTACAGAAACAGATGGCATTGATGCTAGGC
CGGCATGGGGTGTTCCTGGAGCTGAGTGAAGATGTCGAGGAGTATGAGGACCTGACAGAGATCATGTCCA
ATGTACAGCTCAACAGCAACTTCTTGGCCTTAGCTCGGGAGCTGGACATCATGGAGCCCAAGGTGCCTGA
TGACATCTACAAAACCCACCTAGAGAACAACAGGTTTGGGGCAGTGGCTCTCAGGTGGACTCTGCCCGC
ATGAACCTGGCCTCCTCTTTGTGAATGGCTTTGTGAATGCAGCTTTTGCCAAGACAAGCTGCTAACAG
ATGATGGCAACAATGGCTTTACAAGAACAAGGACCACGGAATGTTGAGTGCAGCTGCATCTTTGGGAT
GATTCTGCTGTGGGATGTGGATGGTGGCCTCACCCAGATTGACAAGTACCTGTACTCCTCTGAGGACTAC
ATTAAGTCAGGAGCTCTTCTTGCCTGTGGCATAGTGAAGTCTGGGGTCCGGAATGAGTGTGACCCTGCTC
TGGCACTGCTCTCAGACTATGTTCTCCACAACAGCAACACCATGAGACTTGGTTCATCTTTGGGCTAGG
CTTGGCTTATGCTGGCTCAAATCGTGAAGATGTCTAACACTGCTGCTGCCTGTGATGGGAGATCAAAG
TCCAGCATGGAGGTGGCAGGTGTACAGCTTTAGCCTGTGGAATGATAGCAGTAGGGTCTGCAATGGAG
ATGTAACCTCCACTATCCTTCAGACCATCATGGAGAAGTCAGAGACTGAGCTCAAGGATACTTATGCTCG
TTGGCTTCCTCTTGGACTGGGTCTCAACCACCTGGGAAGGGTGAGGCCATCGAGGCAATCCTGGCTGCA
CTGGAGGTTGTGTCAGAGCCATCCGCAAGTTTGGCAACACACTGGTGGATGTGTGATATGCAGGCT
CTGGGAATGTGCTGAAGGTGCAGCAGCTCTCCACATTTGTAGCGAACACTTTGACTCCAAAGAGAAGGA
GGAAGACAAGACAAGAAGGAAAAGAAAGACAAGGACAAGAAGGAAGCCCTGCTGACATGGGAGCACAT
CAGGGAGTGGCTGTTCTGGGGATTGCCCTATTGCTATGGGGGAGGAGATTGGTGCAGAGATGGCATTAC
GAACCTTTGGCCACTTGTGAGATATGGGGAGCCTACACTCCGGAGGGCTGTACCTTTAGCACTGGCCCT
CATCTCTGTTTCAAATCCAGACTCAACATCCTGGATACCCTAAGCAAATCTCTCATGATGCTGATCCA
GAAGTTTCTATAACTCCATTTTGGCATGGGCATGGTGGCAGTGGTACCAATAATGCCGCTGGCTG
CAATGTGCGCCAGTTAGCTCAATATCATGCCAAGGACCCAAACAACCTTTCATGGTGCCTTGGCACA
GGCCTGACACATTTAGGGAAGGGCACCTTACCCTCTGCCCTACCACAGCGACCGGCAGCTTATGAGC
CAGGTGGCCGTGGCTGGACTGCTCACTGTGCTTGTCTTTCTGGATGTTGAAACATTATTCTAGGCA
AATCACACTATGTATTGTATGGGCTGGTGGCTGCCATGCAGCCCCGAATGCTGGTTACGTTTGTAGGGA
GCTGCGGCCATTGCCAGTGTCTGTCCGTGTGGCCAGGCAGTGGATGTGGTGGGCCAGGCTGGCAAGCCG
AAGACTATCACAGGGTCCAGACGCATACAACCCAGTGTGTTGGCCACGGGGAACGGGCAGAATTGG
CCACTGAGGAGTTTCTTCTGTTACCCCATCTGGAAGGTTTGTATCCTTCGGAAGAACCCCAATTA
TGATCTCTAA
    
```

Restriction Sites: SgfI-MluI

ACCN: NM_001278709

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278709.1](#), [NP_001265638.1](#)

RefSeq Size: 2644 bp

RefSeq ORF: 2250 bp

Locus ID: 5708

UniProt ID: [Q13200](#)

Cytogenetics: 3q27.1

Protein Families: Druggable Genome

Protein Pathways: Proteasome

Gene Summary: The 26S proteasome is a multicatalytic proteinase complex with a highly ordered structure composed of 2 complexes, a 20S core and a 19S regulator. The 20S core is composed of 4 rings of 28 non-identical subunits; 2 rings are composed of 7 alpha subunits and 2 rings are composed of 7 beta subunits. The 19S regulator is composed of a base, which contains 6 ATPase subunits and 2 non-ATPase subunits, and a lid, which contains up to 10 non-ATPase subunits. Proteasomes are distributed throughout eukaryotic cells at a high concentration and cleave peptides in an ATP/ubiquitin-dependent process in a non-lysosomal pathway. An essential function of a modified proteasome, the immunoproteasome, is the processing of class I MHC peptides. This gene encodes one of the non-ATPase subunits of the 19S regulator lid. In addition to participation in proteasome function, this subunit may also participate in the TNF signalling pathway since it interacts with the tumor necrosis factor type 1 receptor. A pseudogene has been identified on chromosome 1. Alternative splicing results in multiple transcript variants of this gene. [provided by RefSeq, Jul 2013]

Transcript Variant: This variant (3) differs in the 5' UTR and the 5' coding region and initiates translation at an alternate start codon, compared to variant 1. It encodes isoform 3, which has a shorter and distinct N-terminus compared to isoform 1.