

Product datasheet for **SC337377**

LRR37A3 (NM_001303255) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRR37A3 (NM_001303255) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRR37A3
Synonyms:	LRR37; LRR37A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001303255, the custom clone sequence may differ by one or more nucleotides

```

ATGGAGTTTCACCGTGTGCCCAGGCTGGTGTGCGAACTCCTGAGCTCAGGCGATCCGCCCGCTCAGCTT
CCCAAAGTGTGGGATTATAGGCGTGAGCCACTGCACCTGGCACAGAAATTTCCAAGGAACTATATTTT
TTACATTGATGGAAATGTATGGAAGCATACAGTTGGACCGAGAACTAATTCTCAGAGAAATAACTTG
ACTGAATTACACAAGGATTCATTTGAAGGCTGCTATCCCTCCAGTATTTAATTCTCAATACAATCCTC
TGACAACCTGTTGAAGATCCATATCTCTTTAAATTGCCAGCATTAAATATCTAGACATGGGAACAACGCT
AGTCCCCTTACAACACTTAAGAACATTCTCATGATGACTGTTGAACTGGAAAACTGATCGTACCTAGC
CATATGGCCTGCTGCCTCTGCCAATTTAAAAACAGCATTGAGGCTGTCTGCAAGACAGTCAAGCTGCATT
GCAACAGTGCATGTCTGACAAACACCACACATTGTCTGAAGAAGCATCTGTAGGGAATCCAGAAGGAGC
GTTTCATGAAGGTGTTACAAGCCCGGAAGAATTACACAAGCACTGAGCTGATTATTGAGCCGGAGGAGCCC
TCAGACAGCAGTGGCATCAACTTGTGAGGCTTTGGGAGTGAGCAGCTAGACACCAATGACGAGAGTGATG
TTACCAGTACACTAAGTTACATCTTACCTTATTTCTCAGCGTTAACCTAGATGTGAAATCACTGTTACT
ACCGTTTCAATAAAGTCCCAACCAGGAAACAGGCTGGCAAGATTCAAAGTGTAGGCAAAAACCGGCGAG
AGACTGAATAGAGTCTCATGGGCCAAGGAGCATCCAGAAAAGGCACTTCAAAGAGGTAGGAAGGCAGA
GCATCAGGAGGGAACAGGGTGGCCAGGCATCTGTGGAGAACACTGCCGAAGAAAAAGGCTCGGAAGTCC
AGCCCCAAGGGAGCTGAAACAGCCTCACACACAGCAGGGGCTGAGAAGTTAGCGGGAAACGCCGCTAC
ACCAAGCCTTCGTTACCCAAGAGCATAAGGCAGCAGTCTCTGTGCTGAAACCTTCTCCAAGGGCGCGC
CTTCTACCTCCAGCCCTGCAAAAGCCCTACCACAGGTGAGAGACAGATGGAAAGACTTAACCCACGCTAT
TTCCATTTTAGAAAGTGCAAAAGGCTAGAGTTACAATATGAAGACATCTAAACCCATCGTACATTCAGA
AAAAATACCGCTTTCACAAAACCTCGCTCCCGCATGACCCAGAACACCCAAGTCAAAAAGAGTCCAA
AGGTGAGAAAGAAAAGTTATCTGAGTAGACTGATGCTCTCAAACAGGCTTCCGTTCTCTGACGGAAGAG
CCTCATAAATTCCTTCAACAAGGGCTTTTTTCATCCTTAAGAGACCTGAGTCTCAAGAAAATCCTTTT
CTGGAAGTATCTGCTCCTCAGAACATTTTATAGAAAACAATAATACAAAAGACACAACCTGCAAGAAATG
CCTTTGAGAAAAATGTTTTATGGAACACTAACATGCCAGAAGGAACCATCTCTGAAAACACAACTA
CAATCATCCTCCTGAGGCAGATTCGCTGGGACTGCATTCAACTAGGGCCAAGTGTAAACAACTGAG
ACAAAATGGGAATACAACAACGTGGGCACTGACCTGTCCCGGAGCCAAAAGCTTCAATTACCCATTGC
TCTCGTCCCAGGTGATCAGTTTGAATTCAGCTAACCCAGCAGTGCAGTCCGTATCCCCAACAAACA
TGTGAGAAGGCTCATTGCTCATGTTATCCGGACCTTGAAGATGGACTGCTCTGGGGCCATGTGCAAGT
ACCTGTGCCAAGCTCGTCTCCAGGACAGGCCACCTGATGAAGCTTCTCAGTGGGACGAGGAAGTAAAGG
CATCCAAGATAGAATGGGATACGGACCAATGGAAGACTGAGAACTACATTAATGAGAGCAGAGAAGCCCA
GAGTGAACAGAAAAGAGAAGTCGCTTGAAGTTCACAAAAGAAGTCCAGGATATGGCTATACAAAAAACTC
ATCTTGGCGTTAATTGTGACTGGAATAACGATTTTGAATTACTTCTCTGCCTCATTGAGATCTGTT
GTCACCGAAGGTCAATACAAGAAGATGAAGAAGGATTCTCAAGGGACAGCGAAGCCCAACGGAGGAGGA
GAGTGAAGCCCTGCCATAG
    
```

Restriction Sites: SgfI-MluI

ACCN: NM_001303255

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001303255.1](#), [NP_001290184.1](#)

RefSeq Size: 2601 bp

RefSeq ORF: 2259 bp

Locus ID: 374819

Cytogenetics: 17q24.1

Protein Families: Transmembrane