

Product datasheet for **SC337380**

ZNF167 (ZKSCAN7) (NM_001288590) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF167 (ZKSCAN7) (NM_001288590) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF167
Synonyms:	ZFP; ZNF64; ZNF167; ZNF448; ZSCAN39
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001288590, the custom clone sequence may differ by one or more nucleotides

```

ATGACCACTGCAGGCAGGGGAAATTTAGGCCTCATCCCCAGGAGCACTGCTTTCCAGAAGCAAGAGGGGC
GCCTGACTGTGAAGCAGGAGCCAGCAAACAGACCTGGGGCAGGGCAGCAGTCTCCAGAAGAATATCC
TCCTGTCTGCGAAATCTTCCGGCTACACTTCAGGCAATTGTGTTACCACGAGATGCTGGGCCGAGGAA
GCATTGAGCCGGCTTCGGGAGCTCTGCCGCTGGTGGCTCATGCCAGAGGTGCACACCAAGGAGCAGATCC
TGGAGCTGCTGGTCTTGAGCAGTTCTGAGCATCCTCCCTGGGAGCTCCGGACCTGGGTGCAGCTGCA
TCACCCTGAGAGTGGTGAGGAGGCTGTGGCTGTGGTGGAGGATTTCAGAGACACCTCAGTGGATCAGAG
GAGGTTTCAGCCCCAGCAGAGAAACAGGAAATGCATTTTGAGGAGACAACAGCTCTGGGTACAACAAAGG
AATCTCCTCTACCTACCCCTCAGTGGGGCTCAGCCCTGGAGCCACCTGGAGCCTCCTTATGACCC
AGGGACACACCACCTCCCAGTGGGGACTTCGCTCAATGTAATCTCCAGTTCCTACCCTCCTCAAGTG
GGGAACTCAGGAGACCAAGCAGGGGCACTGTACTTCGGATGGTCAGGCCCCAGGATACTGTGCATATG
AGGACCTATCTGTAGACTACACTCAGAAGAAATGAAAAGTCTCACACTCAGTCAGAGAGCCCTGCAGTG
GAACATGATGCCAGAAAATCACCATAGCATGGCTCCTTGGCAGGTGAGAACATGATGAAGGGTTCCAGAG
TTGACTCCAAAGCAGGAATTTTTAAAGGATCAGAGTCATCTAACAGGACATCAGGGGGACTCTTTGGGG
TGGTTCCTGGGGCAGCAGAGACTGGAGATGTTTGTGAAGATACTTTCAAGGAGTTAGAAGGACAAACCTC
AGATGAAGAAGGGAGCAGACTAGAAAATGATTTCTTGAAATAACAGATGAAGATAAGAAAAATCCACA
AAAGACAGATATGACAAATATAAGGAAGTTGGGAACATCCACCTCTGTCTTCCAGTCTGTTGAACATG
AAGGAGTTTTAAAGGGACAGAAATCCTATCGATGTGATGAATGTGGCAAAGCTTTCAATCGGAGTTCTCA
CCTTATTGGCCATCAGAGAATCCACACTGGAGAGAAACCTATGAGTGAATGAGTGTGGGAAGACCTTC
AGGCAACCTCCCAGCTCATTGTTTCATCTCAGAACCCACACAGGGGAAAAACCTATGAATGCAGTGAGT
GTGAAAAGGCCTATAGGCACAGCTCCCATCTCATTCAACACCAGAGACTCCATAATGGGGAGAAACCTTA
TAAATGTAATGAATGTGCAAAAGCCTTTACTCAGAGTTCCCGACTCACTGACCACCAGAGAACCCATACT
GGGGAGAAACCTTATGAATGCAATGAGTGTGGAGAGGCATTTCATTGAAGCAAAAGTCTTGCTCGACATC
AGGTCTGCACACTGGTAAGAAACCTTACAAATGCAATGAGTGTGGGAGAGCATTCTGTTCCAATAGAAA
TCTCATTGACCATCAGAGAATCCACACTGGGGAGAAAGCCTTATGAGTGTAGTGAATGTGGCAAAGCCTTC
AGTCGGAGTAAATGTCTTATTCGACATCAGAGCCTCCATACTGGGAAAAAGCCATACAAATGTAGTGAAT
GTGGGAAAGCCTTCAATCAGAACTCTCAACTCATTGAGCATGAGCGAATTCATACTGGAGAAAAACCTTT
TGAATGTAGCGAGTGTGGTAAGGCATTTGGTCTGAGTAAATGTCTTATTCGGCACCAGAGACTTCACACG
GGTGAAAAGCCCTATAAATGCAATGAGTGTGGAAAATCCTTCAATCAAACTCACACCTTATTATACACC
AGAGAATTCACACTGGTGAGAAACCTATGAATGTAATGAGTGTGGGAAGGTATTAGTTATAGTCCAG
CCTTATGGTACATCAGAGAACCCATACTGGGGAAAAACCTATAAATGCAATGATTGTGGGAAAGCTTTT
AGTGACAGCTCACAGCTTATTGTACACCAGAGAGTCCACACCGGAGAGAAACCTTATGAATGTAGTGAAT
GTGGGAAAGCCTTTAGTCAGCGTTCCACTTTAATCACCACCAGCGAACTCACACTGGAGAGAAGTCCTC
AGGTCTGGCTTGGTCAGTTTCTTAA

```

Restriction Sites: SgfI-MluI

ACCN: NM_001288590

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001288590.1, NP_001275519.1</u>
RefSeq Size:	3492 bp
RefSeq ORF:	2265 bp
Locus ID:	55888
UniProt ID:	<u>Q9P0L1</u>
Cytogenetics:	3p21.31
Protein Families:	Transcription Factors
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. Variants 1 and 3 both encode the same isoform (1).