

Product datasheet for **SC337385**

TSC22D2 (NM_001303264) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TSC22D2 (NM_001303264) Human Untagged Clone
Tag:	Tag Free
Symbol:	TSC22D2
Synonyms:	TILZ4a; TILZ4b; TILZ4c
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001303264, the custom clone sequence may differ by one or more nucleotides

```

ATGTCCAAGATGCCGGCCAAGAAGAAGAGCTGCTTCCAGATCACCAGTGTACCACGGCCCAGGTGGCCA
CTAGCATCACCGAGGACACCGAGAGCTTGGACGACCCGGACGAGTCACGCACAGAGGACGTCTCCTCCGA
GATTTTCGACGTCTCTCGGGCCACGGATTATGGCCCTGAGGAGGTCTGCGAGCGCAGCTCTCCGAAGAG
ACGCTTAACAATGTTGGGGATGCGGAGACTCCCGGGACCGTCTCCCAAACCTCCTCCTAGATGGGCAGC
TGGCAGCGGGCGGTGCTGCTCCCGCCAACGGAGGAGGAGTCGTTTCGGCCCGGAGCGTGTCTGGGGCGCT
CGCCAGTACCCTGGCGGGGCTGCCACTTCGGCCCCCGCCCCGGAGCACCCGGCGGGCCCCAGCTCGCG
GGCTCATCCGCCGGGCCAGTGACTGCAGCCCCATCTCAGCCTCCACCACATGTAGTTCCCGTTTTCGCG
TGATCAAGCTGGACCACGGGAGCGGAGAGCCCTATAGACGCGGCCGATGGACGTGTATGGAATACTATGA
GAGGGATTAGACAGCAGCGTCTGACTAGATCCGGGGATTGCATTAGACACAGCAGTACTTTTGACCAG
ACTGCGGAGCGGGACAGCGGCCCTGGGCGCCACCGGAGGGTCGGTGGTGGTAGTAGTGGCCTCCATGCAGG
GGGCGCACGGGCCGAGTCGGGAAGTACAGCTCCTTGACTGCTGTGTACAGCTACCCCGTCGGAGAA
AATGAGCCAGCCACTCCGGCCAGCCGAGAGTTTTAGCGTTGGGCAGCCACAGCCGCCGCCGCCACCC
GTAGGTGGGGCTGTGGCTCAAAGCTCGGCTCCGCTGCCGCGGTTCCCGGGAGCCGCGACCCGGGCCGAGC
CAATGATGGCAGCCGCGCAGCCAGCCAGCCCCAGGGAGCGGGGCCCGGGGGACAGACTCTGCCCGCGAC
GAATGTAACCCTGGCGCAGCCGGCTATGTCCCTGCCTCCGAGCCGGGCCCTGCAGTGGGCGCCCCCGCG
GCGCAGCAGCCCCAGCAGTTTCGGTATCCTCAGCCTCAGATACCGCCCGGACATTTGCTGCCCGTCCAGC
CCTCCGGCCAGAGTGAGTACCTGCAGCAGCAGTGGCCGGCTGCAGCCGCCAAGCCCGCGCAGCCCTC
TGCCACCGCGCCGCGCAGCGAGCCCGCCACGGCGGCCACCCTTCCCGTGGGCACCGGCCAGAATGCTTCC
TCGGTGGGCGCGCAGCTCATGGGCGCGTCTTCCAGCCAGCGAAGCCATGGCCCCCGGACGGGACCAG
CGCAAGGCGGGCAGGTGCGGCCTTGTAGCCGACTGGAGTGCCCCCGGCTACTGTGGGAGGCGTGGTGCA
GCCGTGCCTCGGTCTGCCGGGGTGGGACGCCCAGTCCGTGCCTCCGCGCAGATGGGTGGCAGTGGT
CCGCTGTCAGCCGTACCTGGTGGCCCTCAGCCGTGGTGCCTGGAGTTCCAAACGTGCCTGCAGCCGTGC
CCGCTCCAAGCGTGCCTAGTGTGTCTACCATTCTGTTACTATGCCAAATGTACCCGCGCCTCTGGCCCA
GTCGCAACAGCTGAGCAGCCATACGCCAGTGCAGGAGCAGCAGCATAATCCAGCATGTTGGGCTGCC
TTAGCGCCAGGCACACAGCGCACCAACAAGTCTACCACAGTCTGACCTAAGCCAGTTTCAAACCTCAGA
CCCAGCCTTTAGTCGGGCAAGTCGACGATACTAGAAGAAAATCAGAACCCCTACCTCAACCACCCTTTC
TCTCATTGCTGAAAATAAGCCTGTTGTGAAGCCGCTGTTGCAGATTCCTGGCAAACCCCTTCAGTTA
ACACCTATGAACAGTCTGGCCACCTCTGTATTCAGCATAGCTATTCCTGTTGATGGTGTGAAGACAGTG
CATCTGGGGGAGGTGTTGTAGCCATTGACAACAAAATAGAACAAGCAATGGATCTGGTGAAGGCAATTT
GATGTATGCAGTAAGAGAAGAAGTGAAGTTTTAAAGGAACAAATAAAAGAATTAGTTGAAAGAACTCT
TTACTTGAACGAGAAAATGCACTGTTAAAATCTTTTCAAGCAATGATCAATTATCCCAACTCCCAACCC
AACAGGCCAATCCTGGTAGCACTTCTCAACAGCAAGCAGTGATAGCACAGCCTCCGAGCCAACGCAACC
TCCACAGCAGCCGAATGTCTCCTCAGCATAA
    
```

Restriction Sites: SgfI-MluI

ACCN: NM_001303264

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001303264.1](#), [NP_001290193.1](#)

RefSeq Size: 4785 bp

RefSeq ORF: 2271 bp

Locus ID: 9819

UniProt ID: [O75157](#)

Cytogenetics: 3q25.1

Protein Families: Transcription Factors