

Product datasheet for **SC337394**

BARD1 (NM_001282543) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BARD1 (NM_001282543) Human Untagged Clone
Tag:	Tag Free
Symbol:	BARD1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001282543, the custom clone sequence may differ by one or more nucleotides

```

ATGCCGGATAATCGGCAGCCGAGGAACCCGCGAGCCGAGGATCCGCTCCGGGAACGAGCCTCGTTCGCGGC
CCGCCATGGAACCGGATGGTCGCGGTGCCTGGGCCACAGTCGCGCCGCGCTCGACCGCCTGGAGAAGCT
GCTGCGTGCTCGCGTTGTAATTGTGTAAGTGACTGCATTGGAAGTGGATGCCAGTGTGTACACCCCG
GCCTGGATACAAGACTTGAAGATAAATAGACAACCTGGACAGCATGATTCAACTTTGTAGTAAGCTTCGAA
ATTTGCTACATGACAATGAGCTGTCAGATTTGAAAGAAGATAAACCTAGGAAAAGTTTGTTAATGATGC
AGGAAACAAGAAGAAATCAATAAAAATGTGGTTTAGCCCTCGAAGTAAGAAAAGTCAGATATGTTGTGAGT
AAAGCTTCAGTGCAAAACCCAGCCTGCAATAAAAAAGATGCAAGTGCTCAGCAAGACTCATATGAATTTG
TTTCCCAAGTCTCCTGCAGATGTTTCTGAGAGGGCTAAAAAGGCTTCTGCAAGATCTGAAAAAAGCA
AAAAAGAAAACTTAGCTGAAATCAACCAAAAATGGAATTTAGAGGCAGAAAAAGAAGTGGTGAATTT
GACTCAAAGAGGAATCTAAGCAAAAGCTGGTATCCTTCTGTAGCCAACCATCTGTTATCTCCAGTCCCT
AGATAAATGGTGAATAGACTTACTAGCAAGTGGCTCCTTGACAGAATCTGAATGTTTTGGAAGTTAAC
TGAAGTCTCTTTACCATTGGCTGAGCAATAGAGTCTCCAGACACTAAGAGCAGGAATGAAGTAGTGACT
CCTGAGAAGTCTGCAAAAATTATCTTACATCTAAGAAATCTTTGCCATTAGAAAATAATGAAAAACGTG
GCCATCACAATAGACTTTCCAGTCCCATTTCTAAGAGATGTAGAACCAGCATTCTGAGCACCAGTGGAGA
TTTTGTTAAGCAAACGGTGCCCTCAGAAAATAACCATTGCGCTGAATGTTCTTACCACCTTCATGCAAA
CGTAAAGTTGGTGGTACATCAGGGAGGAAAAACAGTAACATGTCGGATGAATTCATTAGTCTTCCACCAG
GTACACCACCTTCTACATTAAGTAGTTCAAGTTACAGGCGAGTGATGTCTAGTCCCTCAGCAATGAAGCT
GTTGCCAATATGGCTGTGAAAAGAAATCATAGAGGAGAGACTTTGCTCCATATTGCTTCTATTAAGGGC
GACATACCTTCTGTTGAATACCTTTTACAAAATGGAAGTATCCAAATGTTAAAGACCATGCTGGATGGA
CACCATTGCATGAAGCTTGAATCATGGGCACCTGAAGGTAGTGAATTATTGCTCCAGCATAAGGCATT
GGTGAACACCACCGGTATCAAAATGACTCACCCTTACAGTGCAGCCAAGAATGGGCATGTGGATATA
GTCAAGCTGTTACTTTCTATGGAGCCTCCAGAAATGCTGTTAATATATTTGGTCTGCGGCCTGTCGATT
ATACAGATGATGAAAGTATGAAATCGCTATTGCTGCTACCAGAGAAGAATGAATCATCTCAGCTAGCCA
CTGCTCAGTAATGAACACTGGGCAGCGTAGGGATGGACCTCTTGTACTTATAGGCAGTGGGCTGCTTCA
GAACAACAGAAAATGCTCAGTGAGCTTGCAGTAATTCTTAAGGCTAAAAAATACTGAGTTTGACAGTA
CAGTAACCTCATGTTGTTGTTCTGGTGTGCAGTTCAAAGTACCTTGAAGTGTATGCTTGGGATTCTCAA
TGGATGCTGGATTCTAAAATTTGAATGGGTAAGCATGTCTACGAAGAAAAGTATGTGAACAGGAAGAA
AAGTATGAAATTCCTGAAGGTCCACGCAGAAGCAGGCTCAACAGAGAACAGCTGTTGCCAAAGCTGTTTG
ATGGATGCTACTTCTATTTGTGGGAACCTTCAAACACCATCAAAGGACAACCTTATTAAGCTCGTCAC
TGCAGGTGGGGCCAGATCCTCAGTAGAAAAGCCCAAGCCAGACAGTGACGTGACTCAGACCATCAATACA
GTCGCATACCATGCGAGACCCGATTCTGATCAGCGCTTCTGCACACAGTATATCATATGAAGATTTGT
GTAATTATCACCAGAGAGGGTTCGGCAGGGCAAAGTCTGGAAGGCTCCTTCGAGCTGGTTTATAGACTG
TGTGATGTCCTTTGAGTTGCTTCTCTTGACAGCTGA
    
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001282543

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282543.1](#), [NP_001269472.1](#)

RefSeq Size: 5466 bp

RefSeq ORF: 2277 bp

Locus ID: 580

UniProt ID: [Q99728](#)

Cytogenetics: 2q35

Protein Families: Druggable Genome

Gene Summary: This gene encodes a protein which interacts with the N-terminal region of BRCA1. In addition to its ability to bind BRCA1 in vivo and in vitro, it shares homology with the 2 most conserved regions of BRCA1: the N-terminal RING motif and the C-terminal BRCT domain. The RING motif is a cysteine-rich sequence found in a variety of proteins that regulate cell growth, including the products of tumor suppressor genes and dominant protooncogenes. This protein also contains 3 tandem ankyrin repeats. The BARD1/BRCA1 interaction is disrupted by tumorigenic amino acid substitutions in BRCA1, implying that the formation of a stable complex between these proteins may be an essential aspect of BRCA1 tumor suppression. This protein may be the target of oncogenic mutations in breast or ovarian cancer. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2013]

Transcript Variant: This variant (2) lacks an in-frame exon in the 5' coding region, compared to variant 1. The resulting isoform (2, also known as alpha) lacks an internal segment, compared to isoform 1. **Sequence Note:** This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.