

Product datasheet for **SC337399**

KIAA1586 (NM_001286274) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KIAA1586 (NM_001286274) Human Untagged Clone
Tag:	Tag Free
Symbol:	KIAA1586
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001286274, the custom clone sequence may differ by one or more nucleotides

```

ATGGGAGACCCGGGGTCGGAAATAATAGAATCTGTCCCTCCAGCTGGCCCTGAGGCATCTGAGTCAACAA
CGGATGAAAAATGAAGACGACATTCAGTTTGTGAGTATTCTGTTTCTAAAATGCCAAACGACAGGGTGA
TTTTTTCATTTTTTAAATGTGAAGAAGGTGAAAAACAGACAGAAAAATATGAAGTGAGCAAAAAATCAC
TGCAGATTGTCTAAGGCAAGGAACACATTTTCGAGTATATTGAACAACCAATCATTGAAGAAAAGCCAT
CACTTTTCATCAAAGAAAAGAAATAGATAATCTTGTGCTTCCAGATTGTTGAATGAAAAACAGCATTAT
GTTTACAGAACAATACAAATGGCTTGAATAAAAGAAGGTAAATTAGGATGTAAGGATTGTTTCAGCAGTT
CGGCATTTGGGATCGAAAGCAGAAAAGCATGTCCATGTGTCCAAGGAATGGATTGCATATTTAGTAACCC
CTAATGGCAGTAATAAACTACTAGGCAAGCTTCTCTACGAAAAAAAATAGGGAACATGATGTTTCTAA
AGCCCATGGTAAAATTCAGGATTGTTAAAGGAATCACTAATGATTCAATTTGTAATTTAGTGCATAAA
CAAAATAATAAAAAATTTGATGCTACTGTAAAAGTTTTCAATACTGTTTACAGTTTGTAAACATAACA
GACCTTTATCTGATATTGAGGGGCAAGAGAATTACAGGAAAAAATGGAGAGGTAAATTGTTTAAATAC
ACGTTACAGTGCAACAAGAAATAGCAGAACATATTGCAAAAAGAAATGAAGATGAAGATATTTAAGAAATATT
ATAGAAGAGAATGCCAAATCTGTATCATAATTGATGAGGCATCTACAGTTTCAAAGAAAACCCCTAG
TGATTTATCTCCAGTGCACAATTCAGTCAGCTCCTGCACCTGTTATGTTATTTGTGGCTTTAAAAGAATT
GGTGTCAACTATAGCAGAGTGTATTGTCAATACATTATTGACTACTTTAAATGATTGTGGTTTTACAAAT
GAATATTTGAAAGCAAATTTAATTGCATTTTGTCTGATGGTGCTAATACAATCCTGGGAAGAAAGTCTG
GAGTAGCTACAAAATTGTTAGAAAATTTTCTGAAATCATCATTGGAAGTGTAAATCATCGATTACA
ATTGTCACCTTGATGATTCTATATCCGAAATAAAACAAATTAATCATTTAAAAATATTTATTGATAAAAT
TATTCTATTTATCATCAACCTAATAAAAAATCAAACCAAGCTTCTAGGAAGTGTAGCTAAAGAACTTGAAA
CTGAAATTATTAAAAATTGGTCGAGTAATGGGACCAAGATGGGCGGCATGTAGTTTACAAGCTGCTACTGC
TGTATGGCATGCATATCCTATATTATATATGATTTTTCTCATTCTTACTCTGGTTTGGCGAAGAGATTA
GCTAACATTAATTTCTTACAAGACCTTGCTTTAATGATTGACATTCTTGAAGAATTTTCAGTACTTTCAA
CTGCATTACAGTCAAGATCACTAATATTAAGAAAGCAGAAAAATTTGATCAAACGTACCATAAGAGCTTT
GGAAAAATTTAAAAATGGTACTGGAAAGTATGAATCTCAAATGAAGATTTGATCAAGTCAGATAAGTTT
AAAGATATTCATTTAATAAAAAACAATAAATTTAATGCTCTTCTAGGAGTATATTACTAGACAATATAA
TTCAGCACATGAACCTACGCCTTTATCTGACAGAAACCATGAAGATATTTTAAATTTGATTTGCT
GGAACCTTCCACATGGCCTTATGAAGAAATAACTTCACCATGGATAGCTGGTGAAGAAACATTATTTTCA
TTGTGTAAAAATTTAAAAATATGAAGTTGATTTGAATGATTTTCGGGAATTTGTAATAATATATAAAAT
CAACAATGTTTCAATTCCTACAATATATACAAAGCTAAAAAGATAGTTAGCACCATTGCAATCAATAG
TGCTGAAGCTGAAAGGGGTTTCAATTTAATGAACATAATTTGTACAAGGGTGAGAAATAGTTTAAACAATA
GATCATGTATCAGATTTAATGACAATAAATTTACTGGGAAAGAAATTAGCAGATTGGGATGCAACACCGT
TTGTAAAAATCTTGGTCAAATTGCAACCACAGGTTGGCTACAGATACAAGAGTTCCGCAAAAGTCAACAAA
AGTCTTCCATGAGAATCAATTGGCTATATGGAACCTAAAAATAG
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001286274

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001286274.1](#), [NP_001273203.1](#)

RefSeq Size: 2883 bp

RefSeq ORF: 2283 bp

Locus ID: 57691

UniProt ID: [Q9HCI6](#)

Cytogenetics: 6p12.1

Gene Summary: E3 SUMO-protein ligase; facilitates UBE2I/UBC9-mediated SUMO2 modification of target proteins (PubMed:26524493).[UniProtKB/Swiss-Prot Function]
 Transcript Variant: This variant (2) lacks an alternate in-frame exon in the 5' coding region, compared to variant 1, resulting in an isoform (b) that is shorter than isoform a.