

Product datasheet for **SC337425**

ZNF585A (NM_001288800) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF585A (NM_001288800) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF585A
Synonyms:	Zfp27
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001288800, the custom clone sequence may differ by one or more nucleotides

```

ATGCCAGCTAATTGGACCTCACCTCAGAAATCCTCAGCCCTGGCTCCAGAGGATCATGGCAGCTCCTATG
AGGGATCAGTGTCTTCAGGGATGTGGCTATCGATTTTCAGCAGAGAGGAATGGCGGCACCTGGACCCTTC
TCAGAGAAACCTGTACCGGGATGTGATGCTGGAGACCTACAGCCACCTGCTCTCAGTAGGATATCAAGTT
CCTGAAGCAGAGGTGGTCATGTTGGAGCAAGGAAAGGAACCATGGGCACTGCAGGGTGAGAGGCCACGTC
AGAGCTGCCCAGGAGAGAAATTATGGGACCATAATCAATGTAGAAAAATCCTCAGTTATAACAAGTATC
CTCTCAACCTCAAAAAATGTATCCTGGGAGAAAGCTTATGAATGCGCCAAATTTGAAAAGATATTCACC
CAGAAGTCACAGCTCAAAGTACACCTGAAAGTTCTTGCAGGAGAAAAGCTCTATGTATGCATTGAATGTG
GGAAGGCTTTTGTACAGAAGCCAGAATTTATTATACACCAGAAAACCCATATGAGAGAGAAAACCTTTAA
ATGCAATGAATGTGGAATCCTTTTTCAAGTGTGTCCTCTTCAGGCATCAGAGAATTCATACCGGA
GAGAACTCTATGAATGCAGCCAGTGTGGAAAGGCTTCTTTATACTCAGATCTCAGTATACATGAGA
AAATTCATACTGGAGAGAGACCATGAATGCACTGACTGTGGCAAAGCATTACACAAAAGTCCACACT
CAAGATGCATCAGAAAATCCATACAGGCGAGAGATCCTACATCTGTATTGAATGCGGACAGGCCTTCATC
CAGAAGACCCATTTGATTGCACACCGAAGAATTCATACTGGAGAAAAACCATATGAGTGCAGTAAGTGTG
GCAATCCTTCATTTCCAAGTCACAACCTTCAGGTACATCAACGTGTTACACAAAGAGTGAAGCCCTATAT
ATGTACCGAATATGGGAAGGTCTTCAGCAATAATTCCAACCTCGTTACACATAAGAAAGTTCAAAGTAGA
GAGAAATCTTCCATATGTACTGAGTGTGGGAAGGCCTTACCTACAGGTGAGAGTTGATTATTCATCAGA
GAATTCACACTGGAGAGAAAACCTTATGAATGCAGTGTGGGAAAGCCTTCACTCAGAAGTCAGCACT
CACAGTGCATCAGAGAATTCATACAGGAGAAAAATCGTATATATGCATGAAATGTGGACTGGCCTTCATT
CAGAAGGCACACTTGATTGCACATCAAATAATTCATACTGGAGAGAAAACCTCATAAATGTGGTCACTGTG
GGAAATTGTTTACCTCCAAGTCGCAACTCCATGTTTCATAAACGAATTCACACAGGAGAAAAGCCCTATAT
GTGCAATAAATGTGGGAAGGCATTACCAACCGGTCAAATCTCATTACACATCAGAAAACTCATACAGGA
GAGAAATCTTATATATGTTCCAATGTGGAAAGGCCTTACCCAGAGGTGAGACTTGATTACACATCAGA
GAATCCATACTGGGAGAAAGCCTTATGAATGCAATACTGTGGAAAGCCTTCACTCAGAAGTCACACCT
CAATATACATCAGAAAATTCACACTGGAGAGAGACAGTATGAATGCCACGAATGTGGGAAGCCTTCAAC
CAGAAATCAATACTCATTGTTTCATCAGAAAATTCATACAGGAGAGAAAACCTATGTATGCACTGAGTGTG
GAAGAGCTTTATCCGCAAGTCAAACCTTATTACTCATCAAAGAATTCATACTGGAGAGAAAGCCTTATGA
ATGCACTGACTGTGGGAAGTCCTTACCTCCAAGTCTCAGCTCCTGGTGCATCAGCCAGTTCACACAGGA
GAGAAACCTATGTGTGTGCCGAGTGCGGGAAGGCCTTATGAGGAGGTCAAATCTCAGTAAGCACCAGA
AAACTCATACCGGAGAAAAGCCCTACATTTGTTCTGAATGTGGGAAGACCTTTCGACAGAAGTCAGAGTT
GATTACACATCAGAGAATTCATACTGGAGAGAAAACCTTATGAGTGCAGTGTGAGTGTGGGAAGTCTTCACT
AAAAATCAGAGTCCAAGTGCATCAGCGAATTCACACTGGAGAGAAAGCCTTACGTGTGTGCTGAGTGTG
GGAAGGCCTTACTGACAGGTCCAATTTGAATAAACATCAGACAACACACTGGAGACAAAACCTACAA
GTGTGGCATCTGTGGGAAGGCTTCGTTTCAGAAATCAGTGTTTCAGCGTCCATCAGAGCAGCCACGCTGA

```

Restriction Sites: SgfI-MluI

ACCN: NM_001288800

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001288800.1 , NP_001275729.1
RefSeq Size:	6720 bp
RefSeq ORF:	2310 bp
Locus ID:	199704
UniProt ID:	Q6P3V2
Cytogenetics:	19q13.12
Protein Families:	Transcription Factors
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (3) encodes the longer isoform (b). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.