

Product datasheet for **SC337497**

SMARCA6 (HELLS) (NM_001289069) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SMARCA6 (HELLS) (NM_001289069) Human Untagged Clone
Tag:	Tag Free
Symbol:	HELLS
Synonyms:	ICF4; LSH; Nbla10143; PASG; SMARCA6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001289069, the custom clone sequence may differ by one or more nucleotides

```

ATGCCAGCGGAACGGCCCGGGCAGCGCGGCTCGGAGGCTCCAGCAATGGTTGAACAACCTGGACTG
CTGTGATTACCCCGGCCATGCTAGAAGAGGAAGAACAGCTTGAAGCTGCTGGACTAGAGAGAGCGGAA
GATGCTGGAAAAGGCTCGCATGTCTTGGGATAGAGAGTGCACAGAAATTCGGTACCGTAGACTTCAACAT
TTGCTTGAAAAAGCAATATATACTCCAATTTTTATTGACGAAAATGGAACAGCAACAATTAGAGGAAC
AGAAGAAGAAAAGAAAATTGGAGAGAAAAAGGAGTCTTTAAAAGTTAAAAGGGTAAAAATTCAATTGA
TGCAAGTGAAGAGAAGCCAGTTATGAGGAAAAAAGAGGAAGAGAAGATGAATCATAACAATATTTTCAGAG
GTCATGTCAAAGAGGAAATTTGTCTGTGGCTAAAAAAAATAAAAAGGAGAATGAGGATGAAAACCTCT
CCTCTACTAATCTCTGTGTGGAAGATCTTCAGAAAAATAAAGATTGCAATAGTATAATTAAGATAGATT
GTCTGAAACGGTTAGGCAGAATACTAAATCTTTTTGACCCAGTCCGGAAGTGAATGGTCAGCCAGTA
CCTTTTCAACAACCAAAGCACTTACTGGAGGAGTGATGCGATGGTACCAAGTAGAAGGCATGGAATGGC
TTAGGATGCTTTGGGAAAATGGAATTAATGGCATTITAGCAGATGAAATGGGATTGGGTAAAGACAGTTCA
GTGCATTGCTACTATTGCATTGATGATTAGAGAGGAGTACCAGGACCTTTTCTGTCTGTGGCCCTTTG
TCTACACTTCTAACTGGATGGCTGAATTCAAAAGATTTACACCAGATATCCCTACAATGTTATATCATG
GAACCCAGGAGGAACGTCAAAAATTGCATTGCTATTGGAATACTTAATAGTAGATGAAGGACACAGGAT
TAAGAATATGAAGTGCCGTCTAATCAGGGAGTTAAAACGATTCAATGCTGATAACAAACTCTTTTGACT
GGTACTCCCTTGCAAAAACAATTTATCAGAACTTTGGTCATTGCTAAACTTTTTGTTGCCAGATGATTTG
ATGACTTGAAAAGCTTTGAGTCTTGGTTGACATCACTAGTCTTTCTGAAACTGCTGAAGATATTATTGC
TAAAGAAAAGAGAACAAGATGATTGCATATGCTGCACCAGATTTAACACCTTTCTATTGAGAAGACTG
AAGTCTGATGTTGCTCTTGAAGTTCCTCCTAAACGAGAAGTAGTCGTTTATGCTCCACTTCAAAGAAGC
AGGAGATCTTTTATACAGCCATTGTGAACCGTACAATTGCAAACATGTTTGGATCCAGTGAGAAAAGAAC
AATTGAGTTAAGTCTACTGGTCGACAAAACGACGAACTAGAAAATCAATAAATTACAGCAAAATAGAT
GATTTCCCTAATGAATTGAAAAACTGATCAGTCAAATACAGCCAGAGGTGGACCGAGAAAAGAGCTGTTG
TGGAAGTGAATATCCCTGTAGAATCTGAAGTTAATCTGAAGCTGCAGAATAAATGATGCTACTTCGTAA
ATGTTGTAATCATCCATATTTGATTGAATATCCTATAGACCCTGTTACACAAGAATTTAAGATCGATGAA
GAATTGGTAACAAATTCTGGGAAGTTCTTGATTTGGATCGAATGCTGCCAGAACTAAAAAAGAGGTC
ACAAGGTGCTGCTTTTTTCAAAATGACAAGCATGTTGGACATTTTATGATGATTACTGCCATCTCAGAGA
TTTCAACTTCAGCAGGCTTGATGGTCCATGTCTTACTCAGAGAGAGAAAAAACATGCACAGCTTCAAC
ACGGATCCAGAGGTGTTTATCTTCTAGTGAGTACACGAGCTGGTGGCCTGGGCATTAATCTGACTGCAG
CAGATACAGTTATCATTATGATAGTGATTGGAACCCCACTCGGATCTTCAGGCCCAGGATAGATGTCA
TAGAATTGGTCAGACAAAGCCAGTTGTTGTTTATCGCCTTGTACAGCAAATACTATCGATCAGAAAATT
GTGGAAAGAGCAGCTGCTAAAAGGAACTGGAAAAGTTGATCATCCATAAAAATCATTCAAAGGTGGTC
AGTCTGGATTAATCTGTCTAAGAATTTCTTAGATCCTAAGGAATTAATGGAATTAATAAATCTAGAGA
TTATGAAAGGGAATAAAAAGGATCAAGAGAGAAGGTCATTAGTGATAAAGATCTAGAGTTGTTGTTAGAT
CGAAGTGATCTTATTGATCAAATGAATGCTTCAGGACCAATTAAGAGAAGATGGGGATATTCAAGATAT
TAGAAAATCTGAAGATTCCAGTCTGAATGTTTGTAA
    
```

Restriction Sites: SgfI-MluI

ACCN: NM_001289069

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001289069.1](#), [NP_001275998.1](#)

RefSeq Size: 3067 bp

RefSeq ORF: 2421 bp

Locus ID: 3070

UniProt ID: [Q9NRZ9](#)

Cytogenetics: 10q23.33

Gene Summary: This gene encodes a lymphoid-specific helicase. Other helicases function in processes involving DNA strand separation, including replication, repair, recombination, and transcription. This protein is thought to be involved with cellular proliferation and may play a role in leukemogenesis. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jan 2014]
Transcript Variant: This variant (4) uses an alternate in-frame splice site compared to variant 1. The encoded protein (isoform 4) is shorter, compared to isoform 1.