

Product datasheet for **SC337519**

ZNF43 (NM_001256653) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF43 (NM_001256653) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF43
Synonyms:	HTF6; KOX27; ZNF39L1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001256653, the custom clone sequence may differ by one or more nucleotides

```

ATGCCCGACATCCTGGAAGCTGGGAAATGGGACCATTGACATTTATGGATGTGGCCATAGAATTCTGTC
TGGAGGAGTGGCAATGCCTGGCATTGCACAGCAGAATTTATATAGGAATGTGATGTTAGAGAACTACAG
AAACCTGGTCTTCTGGGTATTGCTGTCTCTAAGCCAGACCTGATCACCTGTCTGGAGCAAGAAAAAGAG
CCTTGGGAGCCTATGAGGAGACATGAAATGGTAGCCAAAACCCCAAGTTATGTGTTCTATTTTACCCAAG
ACTTTTGGCCAGAGCAGCATATAAAAGATCCTTTCCAAAAGCGACACTGAGAAGATATAAAAACTGTGA
ACATAAAAAATGTACATTTAAAAAAGACCATAAAAGTGTGGATGAGTGAAGGTGCACAGAGGAGGTTAT
AATGGATTTAACCAATGTTTGCAGCTACCCAGAGCAAAATATTTCTATTTGATAAATGTGTGAAAGCCT
TTCATAAATTTTCAAATTCAAACAGACATAAGATAAGCCATACTGAAAAAAAACCTTTCAAATGCAAAGA
ATGTGGCAAAATCATTTTGCATGCTTCCACATCTAGCTCAACATAAAATAATTCATACCAGAGTGAATTC
TGCAAATGTGAAAAATGTGGAAAAGCTTTTAACTGCCCTCAATCATCACTAAACATAAGAGAATTAATA
CTGGAGAGAAACCCTACACATGTGAAGAATGTGGCAAAGCTTTTAAATGGTCCTCAGCCTTACTACACA
TAAAAAAAATTATACTAGATACAAACTCTACAAATGTGAAGAATGTGGCAAAGCTTTTAAACAAGTCCTCA
ATCCTTACTACCCATAAGATAATTCGCACTGGAGAGAAATTTCTACAAATGTAAAGAATGTGCCAAAGCTT
TTAACCAATCCTCAAACCTTACTGAACATAAGAAAATTCATCCTGGAGAGAAACCTTACAAATGTGAAGA
ATGTGGCAAAGCCTTTAACTGGCCCTCAACTCTTACTAAACATAAGAGAAATTCATACTGGAGAGAAACCC
TACACATGTGAAGAATGTGGCAAAGCCTTTAAACAGTTCTCAAACCTTACTACACATAAGAGAATCCATA
CTGCAGAGAAATTTCTATAAATGTACAGAATGTGGTGAAGCTTTTAGCCGGTCCTCAAACCTTACTAAACA
TAAGAAAATTCATACTGAAAAGAAACCCCTACAAATGTGAAGAATGTGGCAAAGCTTTTAAAGTGGTCTCA
AAGCTTACTGAACATAAGTTAACTCATACTGGAGAGAAACCCCTACAAATGTGAAGAATGTGGCAAAGCCT
TTAACTGGCCCTCAACCCTTACTAAACATAACAGAATTCATACTGGAGAGAAACCCCTACAAATGTGAAGT
ATGTGGCAAAGCCTTTAAACAGTTCTCAAACCTTACTACACATAAGAGAATTCATACTGCAGAAAAACCG
TACAAATGTGAAGAATGTGGCAAAGCTTTTAGCCGGTCCTCAAACCTTACTAAACATAAGAAAATTCACA
TTGAAAAGAAACCCCTACAAATGTGAAGAATGTGGCAAAGCTTTTAAAGTGGTCTCAAAGCTTACTGAACA
TAAGATAACTCATACTGGAGAGAAACCCCTACAAATGTGAAGAATGTGGCAAAGCTTTTAAACATTCTCA
ATCCTTACCAAACATAAGAGGATTCATACTGGAGAGAAACCCCTACAAGTGTGAAGAATGTGGCAAAGCTT
TTACCAATCCTCAAACCTTACTACACATAAGAAAATTCATACTGGAGAGAAATTCATAAATGTGAAGA
ATGTGGCAAAGCTTTTACCAATCTTCAAACCTTACTACACATAAAAAAATTCATACTGGAGGAAAACCC
TACAAATGTGAAGAATGTGGCAAAGCTTTTAAACAGTTCTCAACTTACTAAACATAAGATAATTCACA
CTGAGGAGAAACCCCTACAAATGTGAAGAATGTGGCAAAGCCTTTAAGTGGTCTCAAACCTTACTAAACA
TAAGATAATTCATACTGGAGAGAAACCCCTACAAATGTGAAGAATGTGGCAAAGCTTTTAAACTGTCTCA
ACCCTTCTACACATAAGATTATTCATACTGGAGAGAAACCCCTACAAATGTGAAAAATGTGGCAAAGCTT
TTAACCGATCCTCAAACCTTATTGAACATAAGAAAATTCATACTGGAGAGCAACCCTACAAATGTGAAGA
ATGTGGCAAAGCATTTAACTATTCCTCACACCTTAATACACATAAGAGAATTCATACTAAAGAGCAACCC
TACAAATGTAAAGAATGTGGCAAAGCTTTCAACCAATATTCAAACCTTACTACACATAACAAAATTCATA
CTGGAGAGAAACTCTACAAACCTGAAGATGTGACAGTGATTTTGACAACACCTCAAACCTTTTCAAACAT
AAAATAA
    
```

Restriction Sites: Sgfl-Mlul

ACCN: NM_001256653

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001256653.1</u> , <u>NP_001243582.1</u>
RefSeq Size:	5249 bp
RefSeq ORF:	2457 bp
Locus ID:	7594
Cytogenetics:	19p12
Protein Families:	Transcription Factors
Gene Summary:	<p>This gene belongs to the C2H2-type zinc finger gene family. The zinc finger proteins are involved in gene regulation and development, and are quite conserved throughout evolution. Like this gene product, a third of the zinc finger proteins containing C2H2 fingers also contain the KRAB domain, which has been found to be involved in protein-protein interactions. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (5) differs in the 5' UTR and initiates translation at an upstream start codon, compared to variant 1. The encoded isoform (3) is longer and has a distinct N-terminus, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>