

Product datasheet for **SC337520**

RABEP1 (NM_001291581) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RABEP1 (NM_001291581) Human Untagged Clone
Tag:	Tag Free
Symbol:	RABEP1
Synonyms:	RAB5EP; RABPT5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001291581, the custom clone sequence may differ by one or more nucleotides

```

ATGGCGCAGCCGGGCCCGCTTCCAGCCTGACGAGGATCTGAAGAGGCAAAATGCAGTATTACAAGCTG
CACAAGATGATTTGGGACACCTTCGAACCCAGCTGTGGGAAGCTCAAGCAGAGATGGAGAATATTAAGGC
GATTGCCACAGTCTCTGAGAACCAAGCAAGAAGCTATAGATGAAGTGAAGAACAGTGGAGAGAAGAA
GTTGCTTCACTTCAGGCTGTTATGAAAGAAACAGTTCGTGACTATGAGCACCAGTCCACCTTAGGCTGG
AGCAGGAGCGAACACAGTGGGCACAGTATAGAGAATCCGCAGAGAGGAAATAGCTGATTTAAGAAGAAG
GCTGTCTGAAGGTCAAGAGGAGGAAAATTTAGAAAATGAAATGAAAAAGGCCAAGAGGATGCTGAGAAA
CTTCGGTCCGTTGTGATGCCAATGGAAGGAAATTCAGCTTTGAAGGATAAACTGACAGAGGCTGAAG
ACAAAATTAAGAGCTGGAGGCCCAAAGGTTAAAGAAGTGAATCATTATCTGGAAGCTGAGAAATCTTG
TAGGACTGATCTAGAGATGTATGTAGCTGTTTTGAATACTCAGAAATCTGTTCTACAGGAAGATGCTGAG
AAACTGCGAAAGAATTGCATGAAGTTTCCATCTCTGGAGCAAGAGCGACAACAACACAACAGTTAA
AACATACGTGGCAGAAGGCCAATGACCAGTTTCTGGAATCTCAGCGTTTACTGATGAGAGACATGCAGCG
AATGGAGATTGTGCTAATTCAGAACAGCTCCGACAAGTTGAAGAAGTGAAGAAGAAAGATCAGGAGGAT
GATGAACAACAAAGACTCAATAAGAGAAAGGATCACAAAAAGCAGATGTTGAGGAAGAAAATAAAAATAC
CAGTAGTGTGTGCTTTAACTCAAGAAGAATCTTCAGCCCAGTTATCAAATGAAGAGGAGCATTAGACAG
CACCCGTGGCTCAGTTCATTCTTAGATGCAGGCTTGCTGTTGCCATCTGGAGATCCTTTAGTAAATCG
GACAATGACATGTTTAAAGATGGACTCAGGAGAGCACAGTCTACAGACAGCTTGGAACTCGGGCTCAT
TGCAATCCAAAGCTTTAGGCTATAACTACAAGCAAAAATCTGCTGGAACTGGACGAGTCAAGTTTGG
ACCAGTGGTAGGAGCAGATTCAGTGTCTGAGAATTTGATACTGCATCCCTTGGGTCCTCCAGATGCCA
AGTGGGTTTATGTTAACCAGAGATCAGGAAAGAGCAATCAAGGCGATGACACCAGAACAAGAAGAGACAG
CGTCCCTCTCCAGCGTTACCCAGGGCATGGAGAGTGCCTATGTGTCCCTAGTGGTTATCGTTTAGT
TAGTGAACAGAAATGGAATCTCTTGCAAGAGGTACATAATGCTGGAATAAACTTGGTAGCGTTGT
GATATGTGTTCCAATTACGAAAAACAGTTACAAGGAATTCAGATTCAGGAGGCTGAAACGAGAGACCAGG
TGAAAAAACTACAGCTGATGCTAAGGCAAGCTAATGACCAGTTAGAGAAGACAATGAAAGATAAGCAGGA
GCTGGAAGACTTCATAAAGCAAAGCAGCGAAGATTCGAGTACCAGATCTCTGCACTCGTCTAAGAGCC
CAGGCCTCCGAGATCTTACTGAAGAGTTACAGCAGGGGCTTTCCAGGCAAAGAGGGATGTTCCAGAAC
AGATGGCGGTGCTGATGCAGTACGGAACAGGTTTCAGAAGAGCTGGTGAAGTTACAGAAAGATAATGA
CAGTCTCCAGGAAAGCACAGCCTGCATGTGTCATTACAGCAAGCAGAAGACTTCATCTCCAGACACT
ACAGAGGCACTGCGGGAGTTGGTATTAATAACCGTGAGGACATCATAATGTGCGGACAGCAGCAGACC
ACGTAGAAGAAAAGCTGAAGGCTGAGATACTTTTCTAAAAGAGCAGATCCAAGCAGAACAGTGTAAAA
AGAAAACTTTGAAGAACTCTGCAACTAGAAAATAGAAAATGCAAGGAGGAAATAGCTTCTATTTCTAGC
CTAAAAGCTGAATTAGAAAGAATAAAGTGGAAAAAGGACAGTTGGAGTCCACATTAAGAGAGAAGTCTC
AACAGCTTGAGAGTCTTCAGGAAATAAAGATCAGTTTGGAAAGAGCAGTTAAAGAAAGAGACTGCTGCTAA
GGCTACCGTTGAACAATAATGTTTGAAGAGAAGAATAAAGCTCAGAGATTACAGACAGAATTAGATGTC
AGTGAGCAAGTCCAGAGGATTTTGTAAAGCTTTCACAGACCCTCAGGTGCAGTTAGAGCGGATCCGGC
AAGCTGACTCCTTGAGAGAATCCGGGCAATTCGAATGATACTAACTGACAGACATTAACCAGCTTCC
TGAGACATGA
    
```

Restriction Sites: SgfI-MluI

ACCN: NM_001291581

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001291581.1</u> , <u>NP_001278510.1</u>
RefSeq Size:	5816 bp
RefSeq ORF:	2460 bp
Locus ID:	9135
UniProt ID:	<u>Q15276</u>
Cytogenetics:	17p13.2
Protein Pathways:	Endocytosis
Gene Summary:	<p>Rab effector protein acting as linker between gamma-adaptin, RAB4A and RAB5A. Involved in endocytic membrane fusion and membrane trafficking of recycling endosomes. Involved in KCNH1 channels trafficking to and from the cell membrane (PubMed:22841712). Stimulates RABGEF1 mediated nucleotide exchange on RAB5A. Mediates the traffic of PKD1:PKD2 complex from the endoplasmic reticulum through the Golgi to the cilium (By similarity). [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) lacks an alternate in-frame exon in the 5' coding region, compared to variant 1, resulting in an isoform (3) that is shorter than isoform 1.</p>