

Product datasheet for **SC337556**

SP140 (NM_001278451) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SP140 (NM_001278451) Human Untagged Clone
Tag:	Tag Free
Symbol:	SP140
Synonyms:	LYSP100; LYSP100-A; LYSP100-B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001278451, the custom clone sequence may differ by one or more nucleotides

```

ATGCCCCAGCAGGGCCAGCAGGGGCAGATGGCAAGTGGAGACAGCAATCTCAACTTCAGGATGGTCCGAC
AGATCCAGAACGTAGAGGGTCAGAACCTGCAGGAGCAGGTTTGCCTGAGCCCATTTTCAGTTCTTCAG
AGAAAACAAGGTGGAGATTGCAAGTGAATAACAAGGCCATTTCCCTTTCCTTATGGCCCTCCGAGACCGC
TCCTTCATCTCCGAGCAGATGTATGAACATTTTCAAGAAGCTTTTAGAAACCTGGTCCCAGTGACAAGAG
TGATGTATTGTGTAAGTCAAGTGGAGAAGACATTTGGCTGGTCACATCTGGAAGCATTGTTTCAGCAG
GATTAACCTGATGGCCTATCCTGATTTAAACGAGATTTACAGAAGCTTCCAGAATGTATGCTATGAACAC
TCACCTCTCCAAATGAATAATGTAAACGATTTAGAAGATAGACCCAGATTACTACCATATGGTAAACAAG
AGAACAGCAATGCCTGTCATGAAATGGATGATATAGCAGTGCCTCAGGAAGCCTTGAGCTCCTCGCCAAG
GTGTGAGCCAGGTTTCTTTCAGAGTCTTGTGAGCAGTTAGCTCTCCAAAGGCTGGTGGAGGAGATGCT
GAAGATGCACCCAGCCTACTACCAGGTGGGGGAGTGTCTGTAACCTTGCTATACAAATAGATGAAGGAG
AATCAGAAGAAAATGCCAAAGTTACTGCCTTATGATACAGAAGTTCTAGAAAAGCAACGGGATGATAGATGC
GGCAAGGACATACAGCACAGCACCAGGGGAGAAAACAGGGAGAGGAGGAAGGCAGGAACAGTCCCAGAAAA
AGAAAACCAAGACAAGGAGAAGTACCAAGAGAGTCCAGAGGGAAGAGACAAAGAGACCTTTGATCTAAAAA
CTCCCAAGTCACTAATGAAGGAGAACCAGAGAAGGGGCTCTGTCTACTACCAGGTGAAGGAGAAGAGGG
CAGTGATGACTGTTTCAGAAAATGTGTGATGGAGAAGAGCCCAAGGAGCCTTAGCTCCCTAGCAAGATGT
GGGTGAGTTTCTGTTTTATCTGCAGAGACCTTTGATCTAAAGACTCCCAAGTCACTAATGAAGGAGAAC
CAGAGAAGGAGCTCAGTCTACTACCAGGTGAAGGAGAAGTGTCTAGTGAAGTAAAAATCACCAATGAA
TGAAGAAGGAGAATCAGAAGAGCTTGTCTTACGCTGTATATGATAATGTACCAGGAGCGGAGCAATCA
GCATATGAAAATGAGAAGTGTCTCTGTGTCATGTGTTTTCTCAGAAGAGGTGCCAGGAAGCCAGAAAACA
GGACGGAAAAGTATCAAGCGTGTGGCACAATGGATACTGTGGATATTGCAAAACAACCTCCACTTTGGGAAA
ACCCAAGAGGAAAAGAAGAAAAAGAGGGGGCATGGCTGGAGCAGAATGAGAATGAGAAGGCAGGAAAAC
AGCCAACAAAATGATAATAGCAAAGCCGACGGCCAGGTGGTCTCCAGTGAAGAAAGGCGAACGTGAATC
TGAAAGACCTTTCCAAGATTAGGGGGAGAAAGAGAGGCAACCTGGAACCCGCTTCACTCAGAGTGACAG
AGCTGCACAGAAAAGAGTCCGATCAAGAGCTTCAAGAAAGCACAAAGATGAAACTGTGGATTTTAAGGCT
CCTTTGCTTCCAGTGACCTGTGGTGGGGTGAAGGGAATTTACATAAGAAGAAAATGCAGCAAGGAATCT
TGGTGAAGTGTATACAGACTGAGGATGGAAAATGGTTCACCCCAAGGAAATTTGAAATCAAAGGAGGCCA
TGCAAGATCAAAGAACTGGAGGCTGAGTGTGCGCTGTGGCGGGTGGCCCTACGATGGCTGATGGAGAAT
GGATTTCTGCCTGATCCTCCAAGAATACGTTACAGGAAAAAAAAGAGAATACTGAAGTCTCAAAACAATA
GCTCAGTTGACCTTGTATGAGAAAACCTGGATGAGTGTGAGGTGTGCCGGGACGGAGGGGAGCTGTTCTG
TTGCGACACTTGTTCAGAGTCTTCCATGAGGACTGTCACATCCCGCTGTGGAAGCTGAGAGGACCCCG
TGGAAATTCATCTTCTGCAGGATGAAGGAGTCTCCGGGAAGCCAACAGTGTGTCAGGAATCTGAGGTCC
TGGAGAGGCAGATGTGTCCTGAGGAACAGTTGAAATGTGAGTTCCTCCTTGAAGTCTATTGCTGTTT
TGAGAGCTCCTTTTTGCCAAGATTCCATACTATTATTATTAGAGAGGCGTGTCAAGGCTGAAGGAG
CCCATGTGGTTGGATAAAAATCAAGAAAAGGCTGAATGAGCACGGTTACCCCAAGTGGAGGGGTTGTAC
AAGACATGCGCCTCATCTTCCAGAACCACAGGGCCTTTACAAGTACAAGGATTTTGGCCAAATGGGATT
TAGACTGGAGGCTGAGTTGAGAAGAATTTCAAGGAAGTGTGCTATTACAGAAAACAAATGGGAACAAT
TGA
    
```

Restriction Sites: Sgfl-Mlul

ACCN: NM_001278451

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001278451.1</u> , <u>NP_001265380.1</u>
RefSeq Size:	3189 bp
RefSeq ORF:	2523 bp
Locus ID:	11262
UniProt ID:	<u>Q13342</u>
Cytogenetics:	2q37.1
Protein Families:	Druggable Genome, Transcription Factors
Gene Summary:	<p>This gene encodes a member of the SP100 family of proteins, which are share common domains including an N-terminal homogeneously staining region domain followed by a SP100/autoimmune regulator/NucP41/P75/deformed epidermal autoregulatory factor domain, a plant homeobox zinc finger, and a bromodomain. The encoded protein is interferon-inducible and is expressed at high levels in the nuclei of leukocytes. Variants of this gene have been associated with multiple sclerosis, Crohn's disease, and chronic lymphocytic leukemia. Alternative splicing results in multiple variants. [provided by RefSeq, Aug 2016]</p> <p>Transcript Variant: This variant (3) lacks an alternate in-frame exon compared to variant 1. The resulting isoform (3) has the same N- and C-termini but is shorter compared to isoform 1.</p>