

Product datasheet for **SC337578**

alpha Glucosidase II (GANAB) (NM_001278192) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	alpha Glucosidase II (GANAB) (NM_001278192) Human Untagged Clone
Tag:	Tag Free
Symbol:	GANAB
Synonyms:	G2AN; GIIA; GLUII; PKD3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001278192, the custom clone sequence may differ by one or more nucleotides

```

ATGGCGGCGGTAGCGGCAGTGGCGGCGCGTAGGAGGCGGCTTCTGTCTCTGGTCGTGATGAGAACAGTG
TGGAGTTAACCATGGCTGAGGGACCCTACAAGATCATCTTGACAGCACGGCCATTCCGCCTTGACCTACT
AGAGGACCGAAGTCTTTTGTAGTGTCAATGCCCGAGGACTCTGGAGTTTGAGCATCAGAGGGCCCT
AGGGTCTCTTTCTCGGATAAGGTTAATCTCACGCTTGGTAGCATATGGGATAAGATCAAGAACCTTTTCT
CTAGGCAAGGATCAAAAGACCCAGCTGAGGGCGATGGGGCCAGCCTGAGGAAACACCCAGGGATGGCGA
CAAGCCAGAGGAGACTCAGGGGAAGGCAGAGAAAGATGAGCCAGGAGCCTGGGAGGAGACATTCAAACCT
CACTCTGACAGCAAGCCGTATGGCCCCATGTCTGTGGGTTTGGACTTCTCTCTGCCAGGCATGGAGCATG
TCTATGGGATCCCTGAGCATGCAGACAACCTGAGGCTGAAGGTCAGTGGGTGGGGAGCCATATCGCCT
CTACAATTTGGATGTGTTCCAGTATGAGCTGTACAACCAATGGCCTTGTATGGGCTGTGCCTGTGCTC
CTGGCACACAACCCTCATCGGACTTGGGCATCTTCTGGCTCAATGCTGCAGAGACCTGGGTTGATATAT
CTTCCAACACTGCCGGGAAGACCCTGTTTGGGAAGATGATGGACTACCTGCAGGGCTCTGGGGAGACCC
ACAGACAGATGTTGCTGGATGTGAGAGACTGGCATATTGACGCTTCTCTGCTGCTGGGCCCTCCATC
TCTGATGTTTTCCGGCAATATGCTAGTCTCACAGGAACCCAGGCGTTGCCCCACTCTTCTCCCTCGGCT
ACCACCAGAGCCGTTGGAACACCGGGACGAGGCTGATGTGCTGGAAGTGGATCAGGGCTTTGATGATCA
CAACCTGCCCTGTGATGTGATCTGGCTAGACATTGAACATGCTGATGGCAAGCGGATTTTACCTGGGAC
CCCAGTCGCTTCCCTCAGCCCCGACCATGCTTGAGCGCTTGGCTTCTAAGAGGGCGGAAGCTGGTGGCCA
TCGTAGACCCCCACATCAAGGTGACTCCGGCTACCGAGTTCACGAGGAGCTGCGGAACCTGGGGCTGTA
TGTTAAACCCGGGATGGCTCTGACTATGAGGGCTGGTGTGCTGGCCAGGCTCAGCTGGTACCCTGACTTC
ACTAATCCACGATGAGGGCCTGGTGGCTAACATGTTCACTATGACAATTATGAGGGCTCAGTCCCA
ACCTCTTTGTCTGGAATGACATGAACGAACCATCTGTGTTCAATGGTCCCTGAGGTACCATGCTCAAGGA
TGCCACGATTATGGGGCTGGGAGCACCGGGATGTGCATAACATCTATGGCCTTTATGTGCACATGGCG
ACTGCTGATGGGCTGAGACAGCGCTCTGGGGCATGGAACGCCCTTTGTCTGGCCAGGGCCTTCTTCG
CTGGCTCCAGCGCTTTGGAGCCGTGTGGACAGGGGACAACACTGCCAGTGGGACCATTGGAAGATCTC
TATTCCTATGTGCTCAGCTTGGGGCTGGTGGGACTTTCCTTCTGTGGGGCGGATGTGGGTGGCTTCTC
AAAAACCCAGAGCCAGAGCTGCTTGTGCGCTGGTACCAGATGGGTGCTTACCAGCCATTCTCCGGGCAC
ATGCCACTTGGACTGGGCGACGAGACCCATGGCTGTTACCATCTCAGCACAATGATATAATCCGAGA
TGCCCTGGGCCAGCGATATTCTTTGCTGCCCTTCTGGTACACCCTTTATATCAGGCCCATCGGGAAGGC
ATTCTGTGATGAGGCCCTGTGGGTGAGTACCCTCAGGATGTGACTACCTTCAATATAGATGATCAGT
ACTTGCTTGGGGATGCGTTGCTGGTTCACCCTGTATCAGACTCTGGAGCCATGGTGTCCAGGTCTATCT
GCCTGGCAAGGGGAGGTGTGGTATGACATTCAAAGCTACCAGAAGCATCATGGTCCCAGACCCTGTAC
CTGCCTGTAACCTAAGCAGTATCCCTGTGTTCCAGCGTGGAGGGACAATCGTGCCTCGATGGATGCGAG
TGCGGCGGTCTTCAGAATGTATGAAGGATGACCCATCACTCTCTTTGTTGCACTTAGCCCTCAGGGTAC
AGCTCAAGGAGAGCTCTTCTGGATGATGGGCACAGTCAACTATCAGACTCGCCAAGAGTTCTCTGCTG
CGTCGATTCTCATTCTCTGGCAACACCTTGTCTCCAGCTCAGCAGACCCCTGAAGGACACTTTGAGACAC
CAATCTGGATTGAGCGGGTGGTATAATAGGGGCTGAAAGCCAGCAGCTGTGGTACTCCAGACAAAAGG
ATCTCCAGAAAGCCGCTGTCTTCCAGCATGACCCCTGAGACCTCTGTGTTGGTCCCTGCGCAAGCCTGGC
ATCAATGTGCATCTGATTGGAGTATTACCTGCGATAA
    
```

Restriction Sites: Sgfl-Mlul

ACCN: NM_001278192

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001278192.1, NP_001265121.1</u>
RefSeq Size:	3684 bp
RefSeq ORF:	2559 bp
Locus ID:	23193
UniProt ID:	<u>Q14697</u>
Cytogenetics:	11q12.3
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Metabolic pathways, N-Glycan biosynthesis
Gene Summary:	<p>This gene encodes the alpha subunit of glucosidase II and a member of the glycosyl hydrolase 31 family of proteins. The heterodimeric enzyme glucosidase II plays a role in protein folding and quality control by cleaving glucose residues from immature glycoproteins in the endoplasmic reticulum. Expression of the encoded protein is elevated in lung tumor tissue and in response to UV irradiation. Mutations in this gene cause autosomal-dominant polycystic kidney and liver disease. [provided by RefSeq, Jul 2016]</p> <p>Transcript Variant: This variant (4) lacks three consecutive exons in the 5' coding region but maintains the reading frame, compared to variant 3. The resulting isoform (4) lacks an internal segment, compared to isoform 3.</p>