

Product datasheet for **SC337586**

Nucleolar protein 14 (NOP14) (NM_001291978) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nucleolar protein 14 (NOP14) (NM_001291978) Human Untagged Clone
Tag:	Tag Free
Symbol:	Nucleolar protein 14
Synonyms:	C4orf9; NOL14; RES4-25; RES425; UTP2
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001291978, the custom clone sequence may differ by one or more nucleotides

```

ATGGCGAAGGCGAAGAAGGTGCGGGGCGCAAGGAAGGCCTCCGGGGCGCCGGCGGGAGCGCGAGGGGGCC
CGGCGAAGGCCAACTCCAATCCGTTTCGAGGTGAAAGTTAACAGGCAGAAGTTCAGATCCTGGGCCGGAA
GACGCGCCACGACGTGGGACTGCCCGGGGTGTCTCGCGCACGGGCCCTCAGGAAGCGTACACAGACTTTA
CTAAAAGAGTACAAAGAAAGGGATAAATCCAATGTATTACAGAGATAAACGCTTCGGAGAATACAACAGCA
ACATGAGCCCCGAGGAGAAGATGATGAAGAGTTTGTCTGGAACAGCAGCGACATCATGAGAAAAAAG
CATCTACAATCTAAATGAAGATGAAGAATTGACTCATTATGGCCAGTCTTTGGCAGACATCGAGAAGCAT
AATGACATTGTGGACAGTGACAGCGATGCTGAGGATCGAGGAACGTTGTCTGCTGAGCTGACTGCTGCC
ACTTTGGAGGAGCGGTGGGCTCCTTCACAAGAAGACTCAACAGGAAGGCGAGGAGCGGGAGAAACCGAA
GTCCCGGAAAGAGCTGATTGAAGAGCTCATTGCCAAGTCAAAACAAGAGAAGAGGGAGAGACAAGCTCAA
CGAGAAGATGCCCTCGAGCTCACGGAAGCTAGACCAAGACTGGAAAGAAATTCAGACTCTCCTGTCCC
ACAAAACCTCCCAAGTCAGAGAACAGAGACAAAAGGAAAAACCAAGCCCGATGCATATGACATGATGGT
TCGCGAGCTTGGCTTTGAAATGAAGGCGCAGCCCTCTAACAGGATGAAGACGGAGGCAGAATTGGCAAAAG
GAAGAGCAGGAGCACCTCAGGAAGCTGGAGGCTGAGAGACTTCGAAGAATGCTTGGAAAGGATGAGGATG
AAAATGTTAAGAAACCAAAACATATGTCAGCAGATGATCTGAATGATGGCTTCGTCTAGATAAAGATGA
CAGGCGTTTGTCTTCTACAAAGATGGAAGATGAATGTCGAGGAAGATGTCCAGGAAGAGCAAGCAAG
GAAGCCAGTGACCCTGAGAGCAACGAGGAAGAAGGTGACAGTTTCAGGCGGGGAGGACACAGAGGAGAGCG
ACAGCCCAGATAGCCACTTGGACCTGGAATCCAACGTGGAGAGTGAGGAAGAAAACGAGAAGCCAGCAAA
AGAGCAGAGGCAGACTCCTGGGAAAGGGTTGATAAGCGGCAAGGAAAGAGCTGGAAAAGCTACCAGAGAC
GAGCTGCCCTACAGTTTCGACGCCCTGAATCCTATGAGGAAGTGAATCTCTGTTGTTAGGAAGATCGA
TGAAGAGCAGCTTTTGGTGGTGGAGAGAATTCAGAAAGTGAACACCCAGTCTCGCAGAAAGGAAACAA
AGCAAAATTAGAAAACTGTTTGGCTTTCTTTGGAATACGTTGGCGATTGGCTACAGATGACCCACCA
GACCTCACAGTCATTGATAAGTTGGTTGTGCACTTATATCATCTTGGCAGATGTTTCTGAATCTGCAA
GTGACGCTATCAAATTTGTTCTCCGAGATGCGATGCATGAGATGGAAGAAATGATTGAGACCAAGGCCG
GGCGGCATTGCCAGGTTGGATGTGCTCATTTATTTGAAATCACTGGGCTGCTATTTCAAATCTCCGAC
TTCTGGCACCCAGTGGTGACCCCTGCCCTCGTGTGCCTCAGTCAGCTGCTCACCAAGTGCCCATCCTGT
CCCTCCAGGACGTGGTGAAGGGCTGTTCTGTGCTGCCTGTTCTGGAGTATGTGGCTTTGTCCCAGAG
GTTTATACCTGAGCTTATTAATTTTCTTCTGGGATTCTTTACATAGCAACTCAAACAAAGCAAGCCAA
GGTTCCACTCTGGTGCACCCCTTCAGAGCGCTTGGGAAGAACTCGGAAGTCTCGTGGTGTCTGCTAGAG
AGGATGTGGCCACGTGGCAGCAGAGCAGCCTCTCCCTCCGCTGGGCGAGTAGACTGAGGGCCCCAACTTC
GACAGAGGCCAATCACATCCGACTGTCTGCCTGGCTGTGGGCCTGGCCCTGCTGAAGCGCTGCGTGCTC
ATGTACGGGTCCCTGCCATCCTTCCACGCCATCATGGGCCTCTCCAAGCCCTCCTACGGATCACCTGG
CGGACTGCAGCCACCCGAGGAGCTCCAGGAGCTGTGTGAGAGCACACTGACCGAAATGGAAGCCAGAA
GCAGCTCTGCCGGCCGCTGACCTGTGAGAAGAGCAAGCCTGTCCCACTGAAGCTTTTACACCCCGGCTG
GTCAAAGTCTCGAGTTTGAAGAAAAACAAGGCAGTAGTAAGGAGGAACAGGAAAGGAAGAGGCTGATCC
ACAAACACAAGCGTGAATTTAAAGGGGCCGTTTCGAGAAATCCGCAAGGACAATCAGTTCTTGGCGAGGAT
GCAACTCTCAGAAATCATGGAACGGGATGCGGAAAGAAAGCGGAAAGTAAAGCAGCTTTTAAACAGCCTG
GCTACACAGGAAGGCGAATGGAAGGCTCTGAAGAGGAAAAAGTTCAAAAAATAA

```

Restriction Sites: SgfI-MluI

ACCN: NM_001291978

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001291978.1, NP_001278907.1</u>
RefSeq Size:	3672 bp
RefSeq ORF:	2574 bp
Locus ID:	8602
UniProt ID:	<u>P78316</u>
Cytogenetics:	4p16.3
Gene Summary:	This gene encodes a protein that plays a role in pre-18s rRNA processing and small ribosomal subunit assembly. The encoded protein may be involved in the regulation of pancreatic cancer cell proliferation and migration. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2014] Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (1). Variants 1 and 2 encode the same isoform (1).