

Product datasheet for **SC337593**

Hsp105 (HSPH1) (NM_001286504) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hsp105 (HSPH1) (NM_001286504) Human Untagged Clone
Tag:	Tag Free
Symbol:	HSPH1
Synonyms:	HSP105; HSP105A; HSP105B; NY-CO-25
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001286504, the custom clone sequence may differ by one or more nucleotides

```

ATGGCAACGGCGGCCGTTCTCCGGGGACCGGCTGCCCATTTGGGTAGAATCTTTCCAGAAGGCTCGAGAAG
AAGGAAGCGGAAGTGGCAGCTGGAGGGGCGGTTGAGGCGCCGGTCAGTCATATCATTTGGATCAAAAAA
TAGAACAAATCGGAGTTGCAGCCAAAAATCAGCAAATCACTCATGCAAAACAATACGGTGTCTAACTTCAA
AGATTTTATGGCCGAGCATTCAATGACCCCTTCATTCAAAGGAGAAGGAAAACTTGAGTTACGATTTGG
TTCCATTGAAAAATGGTGGAGTTGGAATAAAGGTAATGTACATGGGTGAAGAACATCTATTTAGTGTGGA
GCAGATAACAGCCATGTTGTTGACTAAGCTGAAGGAAACTGCTGAAAACAGCCTCAAGAAACAGTAACA
GATTGTGTTATTTTCACTCCCTCCTTTTACAGATGCTGAGAGGCGATCTGTGTTAGATGCTGCACAGA
TTGTTGGCCTAACTGTTAAGACTTATGAATGACATGACAGCTGTTGCTTTGAATTACGGAATTTATAA
GCAGGATCTCCCAAGCCTGGATGAGAAACCTCGGATAGTGGTTTTTGTGATATGGGACATTCAGCTTTT
CAAGTGTCTGCTTGTGCTTTTAAACAAGGAAAAATTGAAGTACTGGGAACAGCTTTTGATCTTTCTTAG
GAGGAAAAAACTTCGATGAAAAGTTAGTGAACATTTTTGTGCAGAAATTTAAACTAAGTACAAGTTGGA
TGCAAAATCCAAAATACGAGCACTCCTACGTCTGTATCAGGAATGTGAAAAACTGAAAAAGCTAATGAGC
TCTAACAGCACAGACCTTCCACTGAATATCGAATGCTTTATGAATGATAAAGATGTTTCCGAAAAGATGA
ACAGGTCACAATTTGAAGAACTCTGTGCTGAACCTCTGCAAAAGATAGAAGTACCCCTTTTACTGTT
GGAACAAACTCATCTCAAAGTAGAAGATGTGAGTGCAGTTGAGATTGTTGGAGGCGCTACACGAATTTCCA
GCTGTGAAGGAAAGAATTGCCAAATCTTTGAAAAAGATATTAGCACAACTCAATGCAGATGAAGCAG
TAGCCAGAGGATGTGCATTACAGTGTGCAATACTTTCCCGGCATTTAAAGTTAGAGAATTTCCGTCAC
AGATGCAAGTTCCTTTTCCAATATCTCTGATCTGGAACCATGATTGAGAAGATACTGAAGGTGTTTGA
GTCTTTTAGTCGAAACCATGCTGCTCCTTTCTCAAAGTTCTCACCTTTCTGAGAAGGGGGCCTTTTGAGC
TAGAAGCTTTCTATTCTGATCCCCAAGGAGTTCCATATCCAGAAGCAAAAATAGGCCGCTTTGTAGTTCA
GAATGTTTCTGCACAGAAAGATGGAGAAAAATCTAGAGTAAAAGTCAAAGTGCAGTCAACACCCATGGC
ATTTTACCATCTCTACGGCATCTATGGTGGAGAAAGTCCCAACTGAGGAGAATGAAATGTCTTCTGAAG
CTGACATGGAGTGTCTGAATCAGAGACCACCAGAAAACCCAGACACTGATAAAAAATGCCAGCAAGACAA
CAGTGAAGCTGGAACACAGCCCCAGGTACAAACTGATGCTCAACAAACCTCACAGTCTCCCCCTCACCT
GAACCTACCTCAGAAGAAAAAAAATCCCAGATGCTGACAAAGCAAAATGAAAAAAAGTTGACCAGCCTC
CAGAAGCTAAAAAGCCCAAAAATAAAGGTGGTGAATGTTGAGCTGCCTATTGAAGCCAATTTGGTCTGGCA
GTTAGGGAAGACCTTTAACATGTATATTGAGACAGAGGGTAAGATGATAATGCAAGATAAATTTGGAA
AAAGAAAAGGAATGATGCTAAAAATGCAGTTGAGGAATATGTGTATGAGTTGAGAGACAAGCTGTGTGGAC
CATATGAAAAATTTATATGTGAGCAGGATCATCAAAATTTTTTGGAGACTCCTCACAGAAACTGAAGACTG
GCTGTATGAAGAAGGAGAGGACCAAGCTAAACAAGCATATGTTGACAAGTTGGAAGAATTAATGAAAAT
GGCACTCCAGTTAAAGTTCGGTTTCAGGAAGCTGAAGAACGGCCAAAAATGTTTGAAGAACTAGGACAGA
GGCTGCAGCATTATGCCAAGATAGCAGCTGACTTCAGAAATAAGGATGAGAAATACAACCATATTGATGA
GTCTGAAATGAAAAAGTGGAGAAGTCTGTTAATGAAGTATGGAATGGATGAATAATGTCATGAATGCT
CAGGCTAAAAAGAGTCTTGATCAGGATCCAGTTGTACGTGCTCAGGAAATTAACAAAAAATCAAGGAAT
TGAACAACACATGTGAACCCGTTGTAACACAACCGAAACCAAAAATTAATCACCCTCAACTGGAAGAAG
TCCAAATGGCCCAAAATTTGATAAAAAAGGAAGAAGATTTAGAAGACAAAAACAATTTTGGTGTGAACT
CCACATCAGAATGGTGAATGTTACCCTAATGAGAAAAATCTGTTAATATGGACTTGGACTAG
    
```

Restriction Sites: SgfI-MluI

ACCN: NM_001286504

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001286504.1</u> , <u>NP_001273433.1</u>
RefSeq Size:	5221 bp
RefSeq ORF:	2583 bp
Locus ID:	10808
UniProt ID:	<u>Q92598</u>
Cytogenetics:	13q12.3
Protein Families:	Stem cell - Pluripotency
Gene Summary:	<p>This gene encodes a member of the heat shock protein 70 family of proteins. The encoded protein functions as a nucleotide exchange factor for the molecular chaperone heat shock cognate 71 kDa protein (Hsc70). In addition, this protein plays a distinct but related role as a holdase that inhibits the aggregation of misfolded proteins, including the cystic fibrosis transmembrane conductance regulator (CFTR) protein. Elevated expression of this protein has been observed in numerous human cancers. [provided by RefSeq, Mar 2017]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR, lacks a portion of the 5' coding region, and initiates translation at an alternate start codon, compared to variant 1. The encoded isoform (3) has a distinct N-terminus and is longer than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>