

Product datasheet for **SC337612**

MCK10 (DDR1) (NM_001297652) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MCK10 (DDR1) (NM_001297652) Human Untagged Clone
Tag:	Tag Free
Symbol:	DDR1
Synonyms:	CAK; CD167; DDR; EDDR1; HGK2; MCK10; NEP; NTRK4; PTK3; PTK3A; RTK6; TRKE
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001297652, the custom clone sequence may differ by one or more nucleotides

```

ATGGGACCAGAGGCCCTGTCATCTTTACTGCTGCTGCTCTTGGTGCCAAGTGGAGATGCTGACATGAAGG
GACATTTTTGATCCTGCCAAGTGCCGCTATGCCCTGGGCATGCAGGACCGACCATCCCAGACAGTGACAT
CTCTGTTCCAGTCTCTGGTCAGATTCCACTGCCGCCGCCACAGCAGGTTGGAGAGCAGTGACGGGGAT
GGGGCCTGGTGCCCGCAGGGTCGGTGTTCCTCAAGGAGGAGGAGTACTTGCAGGTGGATCTACAACGAC
TGCACCTGGTGGCTCTGGTGGGCACCCAGGGACGGCATGCCGGGGCCTGGGCAAGGAGTCTCCCGGAG
CTACCGGCTGCGTTACTCCCGGGATGGTCGCGCTGGATGGGCTGGAAGGACCGCTGGGGTCAGGAGGTG
ATCTCAGGCAATGAGGACCCTGAGGGAGTGGTCTGAAGGACCTTGGGCCCCCATGGTTGCCGACTGG
TTCGTTCTACCCCGGGCTGACCGGGTCATGAGCGTCTGTCTGCGGGTAGAGCTCTATGGCTGCCTCTG
GAGGGATGGACTCTGTCTTACACCGCCCTGTGGGCAGACAATGTATTTATCTGAGGCCGTGTACCTC
AACGACTCCACCTATGACGGACATACCGTGGCGGACTGCAGTATGGGGTCTGGCCAGCTGCCAGATG
GTGTGGTGGGGCTGGATGACTTTAGGAAGAGTCAGGAGCTGCGGGTCTGGCCAGGCTATGACTATGTGG
ATGGAGCAACCACAGTCTCCAGTGGCTATGTGGAGATGGAGTTTGGAGTTTGGCCGGCTGAGGGCCTC
CAGGCTATGCAGGTCCACTGTAACAACATGCACACGCTGGGAGCCCGTCTGCCTGGCGGGGTGGAATGTC
GCTTCCGGCGTGGCCCTGCCATGGCCTGGGAGGGGAGCCATGCCCCACAACCTAGGGGGCAACCTGGG
GGACCCAGAGCCCGGGCTGTCTCAGTGCCCTTGGCGGGCGTGGCTCGCTTCTGCAGTGCCGCTTC
CTCTTGGCGGGCCCTGGTACTCTTACAGGAAATCTCCTTATCTCTGATGTGGTGAACAATTCCTCTC
CGGCACTGGGAGGCACCTTCCGCCAGCCCTGGTGGCCGCTGGCCACCTCCCACTTCCAGCAG
CTTGGAGCTGGAGCCAGAGGCCAGCAGCCGTGGCAAGGCCGAGGGGAGCCCGACCCCATCTCATC
GGCTGGTGGCCATCATCTGCTCTGCTCATATTGCCCTCATGCTCTGGCCAGTCTGGCCAGTCTGGC
GCAGGCTCTCAGCAAGGCTGAACGGAGGGTGTGGAAAGAGGAGCTGACGGTTCACTCTCTGTCCCTGG
GGACACTATCCTCATCAACAACCGCCAGGTCCTAGAGAGCCACCCCGTACCAGGAGCCCGGCCCTCGT
GGGAATCCGCCCCACTCCGCTCCCTGTGTCCCAATGGCTCTGCCTACAGTGGGACTATATGGAGCCTG
AGAAGCCAGGCGCCCGCTTCTGCCCCACCTCCCAGAACAGCGTCCCCATTATGCCGAGGCTGACAT
TGTTACCTGCAGGGCGTACCGGGGGCAACACCTATGCTGTGCCTGCACTGCCCCAGGGGCAGTCGGG
GATGGGCCCCAGAGTGGATTTCCCTCGATCTCGACTCCGCTTCAAGGAGAAGCTTGGCGAGGGCCAGT
TTGGGGAGGTGCACCTGTGTGAGGTCGACAGCCCTCAAGATCTGGTTAGTCTTGATTTCCCTTAATGT
GCGTAAGGGACACCTTTGCTGGTAGTGTCAAGATCTTACGGCCAGATGCCACCAAGAATGCCAGGAAT
GATTTCTGAAAGAGGTGAAGATCATGTGAGGCTCAAGGACCCAAACATCATTCCGCTGCTGGCGTGT
GTGTGCAGGACGACCCCTCTGCATGATTACTGACTACATGGAGAACGGCGACCTCAACCAGTTCCTCAG
TGCCACACAGCTGGAGGACAAGGCAGCCGAGGGGGCCCTGGGGACGGGCAGGCTGCGCAGGGGGCCACC
ATCAGTACCCAAATGCTGCTGCATGTGGCAGCCAGATCGCCTCCGGCATGCGCTATCTGGCCACTCA
ACTTTGTACATCGGGACCTGGCCACGGGAAGTGCCTAGTTGGGGAAAATTTACCATCAAAATCGCAGA
CTTTGGCATGAGCCGGAACCTCTATGCTGGGACTATTACCGTGTGCAGGGCCGGCAGTGTGCCATC
CGCTGGATGGCCTGGGAGTGCATCTCATGGGGAAGTTCAGACTGCGAGTGACGTGTGGCCTTTGGTG
TGACCTGTGGGAGGTGCTGATGCTCTGTAGGGCCAGCCCTTTGGGCAGCTCACCGACGAGCAGGTGATC
CGAGAACGCGGGGAGTTCTTCCGGGACCGGGCCGAGGTGTACCTGTCCCGGCCGCTGCCTGCCCG
CAGGGCCTATATGAGCTGATGCTTCGGTCTGGAGCCGGGAGTCTGAGCAGCGACCACCTTTTCCAGC
TGCATCGGTTCTGGCAGAGGATGCACTCAACACGGGTGA
    
```

Restriction Sites: Sgfl-Mlul

ACCN: NM_001297652

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001297652.1</u> , <u>NP_001284581.1</u>
RefSeq Size:	3920 bp
RefSeq ORF:	2631 bp
Locus ID:	780
UniProt ID:	<u>Q08345</u>
Cytogenetics:	6p21.33
Protein Families:	Druggable Genome, Protein Kinase, Transmembrane
Gene Summary:	<p>Receptor tyrosine kinases play a key role in the communication of cells with their microenvironment. These kinases are involved in the regulation of cell growth, differentiation and metabolism. The protein encoded by this gene belongs to a subfamily of tyrosine kinase receptors with homology to Dictyostelium discoideum protein discoidin I in their extracellular domain, and that are activated by various types of collagen. Expression of this protein is restricted to epithelial cells, particularly in the kidney, lung, gastrointestinal tract, and brain. In addition, it has been shown to be significantly overexpressed in several human tumors. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Feb 2011]</p> <p>Transcript Variant: This variant (7) differs in the 5' UTR compared to variant 1. It encodes isoform 1 (also known as DDR1a). Variants 1, 7, and 8 encode the same isoform (1).</p>