

## Product datasheet for **SC337634**

### delta 2 Catenin (CTNND2) (NM\_001288716) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	delta 2 Catenin (CTNND2) (NM_001288716) Human Untagged Clone
Tag:	Tag Free
Symbol:	CTNND2
Synonyms:	GT24; NPRAP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_001288716, the custom clone sequence may differ by one or more nucleotides

```

ATGAGCAGCATGAGTTCAGCAGAAGAGCAGTTCAGTGGCAGTCACAAGATGGTCAAAAAGATATCGAAG
ATGAGCTTACAACAGGTCTCGAGCTGGTGGACTCCTGTATTAGGTCACACTACAGGAATCAGGAATACTTGA
CCCACAGGATTATTCTACAGGTGAAACTGGTCCCGAGCCTCATACAGCAGCCAGCATGGGCACCTGGGC
CCAGAGTTGCGGGCCCTGCAGTCCCAGAACACCACATAGATCCCATCTATGAAGACCGCGTCTATCAGA
AGCCCCCTATGAGGAGTCTCAGCCAGAGCCAGGGGACCCTCTGCCGCCAGCACACACCGGCACCTACCG
CACGAGCACAGCCCCATCTTCCCCTGGTGTGACTCCGTCCCCTTGACAGCGCACAGGCAGCCAGCACGGC
CCACAGAATGCCGCCGCGCCACCTTCCAGAGGGCCAGCTATGCCGCCGGCCAGCCTCAAATTACGCGG
ACCCCTACCGACAGCTGCAGTATTGTCCCTCTGTTGAGTCTCCATACAGCAAATCCGGCCCTGCTCTCC
GCCTGAAGGCACCTTGGCCAGGTCCCCGTCCATTGATAGCATTAGAAAAGATCCCAGAGAATTTGGATGG
AGAGACCCGGAAGTCCCGAAGTGATTGATGTTGCAGCACCAGTTCCCTCGGTCCAGTCAACGCGG
CAGCCTACTTGAACACCTCTGTTTTGGAGACAACAAAATTAAGCCGAGATAAGGAGACAAGGAGGCAT
CCAGCTCCTGGTGGACCTGTTGGATCATCGGATGACCGAAGTCCACCGTAGTGCCTGTGGAGCTCTGAGA
AACCTGGTGTATGGGAAGGCCAACGATGATAACAAAATTGCCCTGAAAAACTGTGGTGGCATCCCAGCAC
TGGTGAGGTTACTCCGCAAGACGACTGACCTGGAGATCCGGGAGCTGGTACAGGAGTCTTTTGGAACT
CTCCTCATGCGATGCACTCAAAATGCCAATCATCCAGGATGCCCTAGCAGTACTGACCAACGCGGTGATT
ATCCCCCACTCAGGCTGGGAAAATTCGCCTCTTCAGGATGATCGGAAAATACAGCTGCATTATCACAGG
TGCTGCGTAACGCCACCGGGTGCCTAAGGAATGTTAGTTCCGCCGGAGAGGAGGCCCGCAGAAGGATGAG
AGAGTGTGATGGGCTTACGGATGCCTTGCTGTACGTGATCCAGTCTGCGTGGGGAGCAGTGAGATCGAT
AGCAAGACCGTTGAAAACCTGTGTGTCATTTTAAAGAACCTCTCGTACCGGCTGGCGGCAGAAACGTCCTC
AGGGACAGCACATGGGCACGGACGAGCTGGACGGGCTACTCTGTGGCGAGGCCAATGGCAAGGATGCTGA
GAGCTCTGGGTGCTGGGGCAAGAAGAAGAAGAAAAGAAATCCCAAGATCAGTGGGATGGAGTAGGACCT
CTTCCAGACTGTGCTGAACCACCAAAAGGGATCCAGATGCTGTGGCACCCATCAATAGTCAAACCTACC
TCACACTGCTCTCTGAGTGTCAAATCCAGACACGCTGGAAGGGGCGGCAGGCGCCTGCAGAACTTGGC
TGCAGGGAGCTGGAAGTGGTCAATATATCCGAGCCGCTGTCCGAAAAGAGAAAGGCTGCCCATCCTC
GTGGAGCTGCTCCGAATAGACAATGACCGTGTGGTGTGCGCGGTGGCCACTGCGTGCAGAACATGGCCT
TGGACGTCAGAAAATAGGAGCTCATCGGCAATACGCCATGCGAGACCTAGTCCACAGGCTTCCAGGAGG
GAACAACAGCAACAACACTGCAAGCAAGGCCATGTCGGATGACACAGTGACAGCTGTCTGTGCACACTG
CAGCAAGTGATTACCAAGAACATGGAGAACGCCAAGGCCTTACGGGATGCCGGTGGCATCGAGAAGTTGG
TCGGCATCTCAAAGCAAGGAGATAAACTCTCAAAGTGGTCAAGGCTGCATCTCAGGTCCTCAA
CAGCATGTGGCAGTACCGAGATCTGAGGAGTCTCTACAAAAGGATGGATGGTCAACAATACCACTTTGTA
GCCTCGTCTTCAACCATCGAGAGGGACCGCAAAGGCCCTACTCCTCCTCCCGCAGCCCTCCATCTCCC
CTGTGCGCGTGTCTCCCAACAACCGCTCAGCAAGTGGCCAGCTTACCTCGGGAATGATCAGCCTCAA
AGAAAGGAAAACAGACTACGAGTGCACCGGCAGCAACGCCACCTACCACGGAGCTAAAGGCGAACACACT
TCCAGGAAAGATGCCATGACAGCTCAAAACACTGGAATTTCACTTTGTATAGGAATTTCTATGGTGGC
CCGCTGAAGACATCAAACACAACCAGTTTCAGCACAGCCAGTCCCACAGGAGCCAGCAGAAAAGATTA
CGAGACCTACCAGCCATTTCAAGATTCCACAAGAAATTACGATGAGTCTTCTTCGAGGACAGGTCCAC
CATCGCCCTCCCGCCAGCGAGTACACCATGCACCTGGGTCTCAAGTCCACCGGAACTACGTTGACTTCT
ACTCAGCTGCCCGTCCCTACAGTGAAGTGAAGTATGAAACGAGCCACTACCCGGCCTCCCGGACTCCTG
GGTGTGA
    
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001288716

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001288716.1</a> , <a href="#">NP_001275645.1</a>
<b>RefSeq Size:</b>	4551 bp
<b>RefSeq ORF:</b>	2667 bp
<b>Locus ID:</b>	1501
<b>UniProt ID:</b>	<a href="#">Q9UQB3</a>
<b>Cytogenetics:</b>	5p15.2
<b>Protein Families:</b>	Druggable Genome
<b>Gene Summary:</b>	<p>This gene encodes an adhesive junction associated protein of the armadillo/beta-catenin superfamily and is implicated in brain and eye development and cancer formation. The protein encoded by this gene promotes the disruption of E-cadherin based adherens junction to favor cell spreading upon stimulation by hepatocyte growth factor. This gene is overexpressed in prostate adenocarcinomas and is associated with decreased expression of tumor suppressor E-cadherin in this tissue. This gene resides in a region of the short arm of chromosome 5 that is deleted in Cri du Chat syndrome. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Dec 2013]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR and coding region, uses a downstream start codon, and lacks two alternate exons in the 5' coding region, compared to variant 1. The resulting protein (isoform 3) has a shorter and distinct N-terminus compared to isoform 1.</p> <p>Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>