

Product datasheet for **SC337638**

PER3 (NM_001289864) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PER3 (NM_001289864) Human Untagged Clone
Tag:	Tag Free
Symbol:	PER3
Synonyms:	FASPS3; GIG13
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001289864, the custom clone sequence may differ by one or more nucleotides

```

ATGGTTGCCATACACCAAAAAGTTTTGAAGTATGCAGGGCATCCTCCCTTGAACATTCTCCCATTGAT
TTTGTACTCAAACGGAGACTACATCATACTGGATTCCAGTTGGTCCAGCTTTGTGAATCCCTGGAGCCG
GAAGATTTCTTTCATCATTGGTCGGCATAAAGTTCGAACGAGCCCACTAAATGAGGATGTTTTGTACT
AAAATTAAGATGAACGATAATGACAAAAGACATAACAGAATTACAAGAACAAATTTACAACTTCTCT
TACAGCCAGTTACGCTGAGCGTGTCCAGCGGCTACGGGAGCCTGGGGAGCAGCGGGTGCAGGAGCAGCT
TGTGAGCATCGCCTCCTCCAGTGAAGCCAGTGGGCACCGTGTGGAGGAGACGAAGGCGGAGCAGATGACC
TTGCAGCAGGTCTATGCCAGTGTGAACAAAATTAATACTGGGTGAGCAGCTCTACATTGAGTCAATGA
CCAAATCATCATTCAAGCCAGTGAAGCCAGTGAAGCCAGTGAAGCCAGTGAAGCCAGTGAAGCCAGTGA
TGGTGGTGAATGTAAAGACCTTTACTTCTCCACCAAACTGAAAAACAATAGTGTGTACTGAGCCCC
TGTGAGGATTTGAGGAACGATGAGCACAGCCCATCTATCAACAGATCACTGTATGCAGAGTGCATCA
GATACCTGAAGAGCTACAACATTCAGCTTTGAAAAGAAAGTGTATCTCCTGTACAAATACAATCTTTC
CTCCTCAGAAGAAGACAAACAGAACCAAGGCAGATGATGTCCAAGCCTTACAAGCTGGTTTGAATC
CCAGCCATACCTAAATCAGAAATGCCAACAAATGGACGGTCCATAGACACAGGAGGAGGACTCCACAGA
TCCTGTCCACGGCGATGCTGAGCTTGGGGTCGGGCATAAGCCAATGCGGTTACAGCAGCACCATTGTCCA
TGTCCCACCCAGAGACAGCCAGGGATGCTACCTCTTCTGTGAGCCCTGGACCCTGAACATGCAGCCA
GCCCTTTGACCTCGGAAGAATTTAAACACGTGGGGCTCACAGCGGCTGTCTGTGAGCGCACACCCAGA
AGGAAGAGCAGAATTATGTTGATAAATCCGAGAAAAGATCCTGTATCACCCCTACAGCTCCTATCTTCA
GCAAGAAAGCAGGAGCAAAGCTAAATATTCATATTTTCAAGGAGATTCTACTTCAAGCAGACGCGGTCCG
GCCGCTGCAGGAAAGGGAAGCACAAGCGGAAGAAGCTGCCGGAGCCGCCAGACAGCAGCTCGAACA
CCGGCTCTGGTCCCCGAGGGGAGCGCATCAGAACGCACAGCCCTGCTGCCCTCCGCGGCTCCTCTCC
GCACACCTCGAGCCCGACCTTCCACCTGCCGCCATGGTGCCAGCCAGGCCCTTACCTCGTCCCAGCT
TTTCCCCTCCCAGCCGCGACCTCACCCGGAAGAGAATACGCAGCCCCGAACTGCACCGGAAGGCTGC
ATGGGCTGCCCTTGTCCGAGGGCTTGCAGCCTTACCCAGCTTTCCCTTTTCTTACTTGGATACTTTTAT
GACCGTTTTCTGCCTGACCCCTGTCTGTCTCTGTTGTCGCCATCGTTTTTGGCATGCCATTCTCTG
GGGGCAGCAGCTTCTGCGATATCACCTCAATGTCGTGAGCAATGAGTCCAATCTGGACCCACCCC
CTTCAGTACCAGCCAAAGGAGAGAGGAGAAAAGTGGGAGGCACAAAGCGAGGGGCACCCGTTCCATTAC
TTCGAGAAGCAGCTCACCTTGCAGTAAACTTACTTCAAGGAAGAGATGCCAGACCTCTGAATCTCCA
GATCAGATGAGAAGGAACACGTGCCACAACTGAGTATTGTGTTACAGGCAACAATGGCAGTGAGAGCA
GTCCTGTACTACCGGTGCACTGTCCACGGGTGACCTCCCAGGAGAATCCATCCCATCCTACTGCCAG
CGCTCTGTCCACAGGATCGCCTCCCATGAAGAATCCATCCCATCCTACTGCCAGCGCTCTGTCCACAGGA
TCGCCTCCCATGAAGAATCCATCCCATCCTACTGCCAGCACACTGTCATGGGATTGCTCCCAGCAGGA
CTCCATCCCATCCTACTGCCACTGTTCTGTCCACGGGGTCACTCCCAGCGAATCCCATCCAGAAGTGG
TTCAGCAGCATCAGGAAGCAGCGACAGCAGTATACCTTACTAGTAGTGTATTTCTTCTAAAATCTCC
CAAAATGGGCAGCAATCTCAGGACGTACAGAAAAAAGAAACATTTCCCTAATGTGCCGAAGAGCCCATCT
GGAGAATGATACGGCAGACACCTGAGCGCATTCTCATGACATACCAGGTACCTGAGAGGGTTAAAGAAAT
TGTAATAAAGAACCTGGAAAAGCTAGAAAAGTATGAGGCAGCAGCAGCCAGTTTCTCATGGGCAA
AAGGAGGAGCTGGCTAAGGTGTATAATTGATTCAAAGCCAGACTGCACTCAAGAAATCGACATTAAG
CCTGTGCACTTGTGAAAATGAAGATTCAGCTGATGGTGCGCCACATCCTGTGTCAGGTTCTGGTAGA
AGACAGCTGTGA
    
```

Restriction Sites: SgfI-MluI
ACCN: NM_001289864

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001289864.1 , NP_001276793.1
RefSeq Size:	6394 bp
RefSeq ORF:	2673 bp
Locus ID:	8863
UniProt ID:	P56645
Cytogenetics:	1p36.23
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Circadian rhythm - mammal
Gene Summary:	<p>This gene is a member of the Period family of genes and is expressed in a circadian pattern in the suprachiasmatic nucleus, the primary circadian pacemaker in the mammalian brain. Genes in this family encode components of the circadian rhythms of locomotor activity, metabolism, and behavior. This gene is upregulated by CLOCK/ARNTL heterodimers but then represses this upregulation in a feedback loop using PER/CRY heterodimers to interact with CLOCK/ARNTL. Polymorphisms in this gene have been linked to sleep disorders. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2014]</p> <p>Transcript Variant: This variant (5) differs in the 5' UTR and coding sequence and uses an alternate in-frame splice junction at the 5' end of an exon compared to variant 1. The resulting isoform (5) is shorter at the N-terminus and lacks a single internal aa compared to isoform 1.</p>