

## Product datasheet for **SC337646**

### **PACE4 (PCSK6) (NM\_001291309) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	PACE4 (PCSK6) (NM_001291309) Human Untagged Clone
Tag:	Tag Free
Symbol:	PCSK6
Synonyms:	PACE4; SPC4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_001291309, the custom clone sequence may differ by one or more nucleotides

```

ATGCCTCCGCGCGCGCCCTGCGCCCGGGCCCCGGCCGCGCCCGGGCCGCGCCGCCACCGACACCG
CCGCGGGCGCGGGGGCGCGGGGGCGCGGGGGCGCGGGGGCCCGGGCGGGCCCGGGTTCGGCCGCTCGCGCCGCG
TCCTGGCGCTGGCTGCTGCTGCTGGCGCTGCCTGCCCGCTGCTCCGCGCCCGCGCGCCCGTCTAC
ACCAACCACTGGGCGGTGCAAGTCTGGGCGGCCCGGCCGAGGCCGACCGCGTGGCGGCGGCGCACGGGT
ACCTCAACTTGGGCCAGATTGGAACCTGGAAGATTACTACATTTTTATCACAGCAAAACCTTTAAAAG
ATCAACCTTGAGTAGCAGAGGCCCTCACACCTTCTCAGAATGGACCCCAAGGTGAAATGGCTCCAGCAA
CAGGAAGTGAACGAAGGTGAAGAGACAGGTGCGAAGTGACCCGAGGCCCTTTACTTCAACGACCCCA
TTTGGTCAACATGTGGTACCTGCATTGTGGCGACAAGAACAGTCGCTGCCGGTCGGAATGAATGTCCA
GGCAGCGTGAAGAGGGGCTACACAGGAAAAACGTGGTGGTACCATCCTTGATGATGGCATAGAGAGA
AATCACCTGACCTGGCCCAAATATGATTCTACGCCAGCTACGACGTGAACGGCAATGATTATGACC
CATCTCCACGATATGATGCCAGCAATGAAAATAAACACGGCACTGTTGTGCGGGAAGATTGCTGCTTC
AGCAAACAATTCCTACTGCATCGTGGGCATAGCGTACAATGCCAAAATAGGAGGCATCCGCATGCTGGAC
GGCGATGTCACAGATGTGGTCGAGGCAAAGTCGCTGGGCATCAGACCCAACTACATCGACATTTACAGTG
CCAGCTGGGGGCCGACGACGACGGCAAGACGGTGGACGGGCCCGGCCGACTGGCTAAGCAGGCTTTCGA
GTATGGCATTAAAAAGGGCCGCGCAGGGCCTGGGCTCCATTTTCGTCTGGGCATCTGGGAATGGCGGGAGA
GAGGGGGACTACTGCTCGTGGATGGCTACACCAACAGCATCTACACCATCTCCGTCAGCAGCGCCACCG
AGAATGGTACAAGCCCTGGTACCTGGAAGAGTGTGCCTCCACCTGGCCACCACCTACAGCAGTGGGGC
CTTTATGAGCGAAAAATCGTCACCACGGATCTGCGTCAGCGCTGTACCGATGGCCACACTGGGACCTCA
GTCTCTGCCCCCATGGTGGCGGCATCATCGCCTGGCTCTAGAAGCAAAGAGCATCCCCTTAGTGCAGG
TGCTGCGGACTACGGCCCTGACCAGCGCCTGCGCGGAGCACTCGGACCAGCGGGTGTCTACTTGGAGCA
CGTGGTGGTTCGCACCTCCATCTCACACCCACGCCGAGGAGACCTCCAGATCTACCTGGTTTCTCCCTCG
GGAACCAAGTCTCAACTTCTGGCAAAGAGGTTGCTGGATCTTTCCAATGAAGGGTTTACAACTGGGAAT
TCATGACTGTCCACTGCTGGGAGAAAAGGCTGAAGGGCAGTGGACCTTGGAAATCCAAGATCTGCCATC
CCAGGTCCGCAACCCGGAGAAGCAAGGGAAGTTGAAAGAATGGAGCCTCATACTGTATGGCACAGCAGAG
CACCCGTACCACACCTTCAGTGCCATCAGTCCCCTCGCGGATGCTGGAGCTCTCAGCCCCAGAGCTGG
AGCCACCAAGGCTGCCCTGTCACCCTCCAGGTGGAAGTTCTGAAGATGAGGAAGATTACACAGCTCA
ATCCACCCAGGCTCTGTAATATTTACAGACCAGTGTGTGCCATCCGGAGTGTGGTGACAAGGCTGT
GATGGCCCAATGCAGACCAAGTCTTGAAGTGCCTCCACTTACGCTGGGAGTGTCAAGACCAGCAGGA
AGTGCGTGAGTGTGTGCCCTTGGGCTACTTTGGGGACACAGCAGCAAGACGCTGTGCCCGGTGCCACAA
GGGGTGTGAGACCTGCTCCAGCAGAGCTGCGACGCAAGTGCCTGTCTTGGCCCGCGGGTTCTATACCAC
CAGGAGATGAACACCTGTGTGACCCTCTGTCTGCAGGATTTTATGCTGATGAAAGTCAAAAAATTGCC
TTAAATGCCACCAAGCTGTAAAAAGTGCCTGGATGAACCTGAGAAATGTACTGTCTGTAAAGAAGGATT
CAGCCTTGACGGGGCAGCTGCATTCTGACTGTGAGCCAGGCACCTACTTTGACTCAGAGCTGATCAGA
TGTGGGAATGCCATCACACCTGCGGAACCTGCGTGGGGCCAGGCAGAGAAGAGTGCACTTCACTGTGCGA
AAAACCTCCACTTCCACGACTGGAAGTGTGTGCCAGCCTGTGGTGGAGGCTTCTACCCAGAAGAGATGCC
GGGCTTGCCCAAAAAGTGTGCGAAGGTGTGACGAGAAGTCTTGGAGCTGTGCAGGCTCCAGCAGGAAC
TGTAGCAGGTGTAAAGACGGCTTACACAGCTGGGGACCTCCTGCATCACCACACACGTGCAGCAACG
CTGACGAGACATTCTGCGAGATGGTGAAGTCCAACCGGCTGTGCGAACGGAAGCTCTTATTCAAGTCTG
CTGCCGACGTCCTCCTGGCCGGTAA
    
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001291309

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001291309.1</a> , <a href="#">NP_001278238.1</a>
<b>RefSeq Size:</b>	4349 bp
<b>RefSeq ORF:</b>	2688 bp
<b>Locus ID:</b>	5046
<b>UniProt ID:</b>	<a href="#">P29122</a>
<b>Cytogenetics:</b>	15q26.3
<b>Protein Families:</b>	Druggable Genome, Protease, Secreted Protein
<b>Gene Summary:</b>	<p>This gene encodes a member of the subtilisin-like proprotein convertase family, which includes proteases that process protein and peptide precursors trafficking through regulated or constitutive branches of the secretory pathway. The encoded protein undergoes an initial autocatalytic processing event in the ER to generate a heterodimer which exits the ER and sorts to the trans-Golgi network where a second autocatalytic event takes place and the catalytic activity is acquired. The encoded protease is constitutively secreted into the extracellular matrix and expressed in many tissues, including neuroendocrine, liver, gut, and brain. This gene encodes one of the seven basic amino acid-specific members which cleave their substrates at single or paired basic residues. Some of its substrates include transforming growth factor beta related proteins, proalbumin, and von Willebrand factor. This gene is thought to play a role in tumor progression and left-right patterning. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Feb 2014]</p> <p>Transcript Variant: This variant (9) lacks two consecutive exons in the central coding region but maintains reading frame, compared to variant 1. The encoded isoform (i) is shorter, compared to isoform PACE4-AI.</p>