

## Product datasheet for **SC337652**

### **SAP97 (DLG1) (NM\_001290983) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	SAP97 (DLG1) (NM_001290983) Human Untagged Clone
Tag:	Tag Free
Symbol:	DLG1
Synonyms:	DLGH1; hdlg; SAP-97; SAP97
Vector:	<u>pCMV6 series</u>



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_001290983, the custom clone sequence may differ by one or more nucleotides

```

ATGCCGGTCCGGAAGCAAGATACCCAGAGAGCATTGCACCTTTTGGAGGAATATCGTTCAAACCTAAGCC
AAACTGAAGACAGACAGCTCAGAAGTTCCATAGAACGGGTATTAACATATTTTCAGAGCAACCTCTTTCA
GGCTTTAATAGATATTCAGAATTTTGAAGTGACCTTACTGGATAATCCAAAATGTATAGATCGTTCA
AAGCCGCTCTGAACCAATTC AACCTGTGAATACTTGGGAGATTTCCAGCCTTCCAAGCTCTACTGTGACTT
CAGAGACTGCAAGCAGCCTTAGCCCTAGTGTAGAGAAATACAGGTATCAGGATGAAGATACACCTCC
TCAAGAGCATATTTCCCAAAAATCACAATGAAGTGATAGGTCCAGAATTGGTTCATGTCTCAGAGAAG
AACTTATCAGAGATTGAGAATGTCCATGGATTTGTTTCTCATTCTCATATTTACCAATAAAGCCAACAG
AAGCTGTTCTCCCTCTCCTCCACTGTCCCTGTGATCCCTGTCCCTGCCAGTCCCTGCTGAGAATACTGT
CATCTACCCACCATACCACAGGCAATCCTCCCCAGTACTGGTCAACACAGATAGCTTGGAAACACCA
ACTTACGTTAATGGCACAGATGCAGATTATGAATATGAAGAAATCACACTTGAAGGGGAAATTCAGGGC
TTGGTTTCAGCATTGCAGGAGGTACGGACAACCCACACATTGGAGATGACTCAAGTATTTTCATTACCAA
AATTATCACAGGGGGAGCAGCCGCCAAGATGGAAGATTGCGGGTCAATGACTGTATATTACGAGTAAAT
GAAGTAGATGTTCTGTGATGAACACATAGCAAAGCAGTTGAAGCGTTGAAAGAAGCAGGGTCTATTGTAC
GCTTGTATGTAAGAAAGAAAGGAAACAGTGTGAGAAAAAATATGGAATAAAGCTCATTAAAGTCTTAA
AGGTCTTGGGTTTAGCATTGCTGGAGGTGTTGGAATCAGCATATTCCTGGGGATAATAGCATCTATGTA
ACCAAAAATAATGAAGGAGGTGCAGCACATAAGGATGGCAAACTTCAGATTGGAGATAAACTTTTAGCAG
TGAATAACGTATGTTTGAAGAAGTTACTCATGAAGAAGCAGTAACTGCCTTAAAGAACACATCTGATTT
TGTTTATTTGAAAGTGGCAAAACCCAAAGTATGTATATGAATGATGGCTATGCACCACCTGATATCACC
AACTCTTCTTCAGCCTGTTGATAACCATGTTAGCCCATCTTCTTCTTGGGCCAGACACCAGCATCTC
CAGCCAGATACTCCCAAGTTTCTAAAGCAGTACTTGGAGATGATGAAATACAAGGGAACCTAGAAAAGT
TGTTCTTCATCGTGGCTCAACGGCCTTGGTTTCAACATTGTAGGAGGAGAAGATGGAGAAGGAATATTT
ATTTCTTTATCTTAGCCGGAGGACCTGCTGATCTAAGTGGAGAGCTCAGAAAAGGAGATCGTATTATAT
CGGTAACAGTGTGACCTCAGAGCTGCTAGTCATGAGCAGGCAGCAGCTGCATTGAAAAATGCTGGCCA
GGCTGTCACAATTGTTGCACAATATCGACCTGAAGAATACAGTCGTTTTGAAGCTAAAATACATGATTTA
CGGGAGCAGATGATGAATAGTAGTATTAGTTCAGGGTCAGGTTCTCTTGAAGTAAAGTAAAGTAAAGT
TCTATGTCAGAGCCCTTTTGTATGACAAGACTAAAGACAGTGGGCTTCCAGTCAGGGACTGAACTT
CAAATTTGGAGATATCCTCCATGTTATTAATGCTTCTGATGATGAATGGTGGCAAGCCAGGCAGGTTACA
CCAGATGGTGAAGCGATGAGGTGCGAGTGATTCCAGTAAACGCAGAGTTGAGAAGAAAGAACGAGCCC
GATTA AAAACAGTGAAATTC AATTTCAAACAGGAGATAAAGGGGAGATCCCTGACGACATGGGATCAAA
AGGCCTGAAGCATGTAACCTTCTAATGCCAGCGATAGTGAAGTAGTTACCGTGGTCAAGAAGAATACGTC
TTATCTTATGAACCAGTGAATCAACAAGAAGTTAATTACTCGACCAGTGCATATTGGGACCTATGA
AAGACAGGATAAATGATGACTTGATCTCAGAATTTCTGACAAAATTTGGATCCTGTGTTCTCATAACAAC
TAGACAAAACGAGATTATGAGGTAGATGGAAGAGATTATCATTTTGTGACTTCAAGAGAGCAGATGGAA
AAAGATATCCAGGAACATAAATTCATTGAAGCTGGCCAGTATAACAATCATCTATATGGAACAAGTGTTG
AGTCTGTACGAGAAGTAGCAGAAAAGGGCAACACTGTATCCTTGTGATGTCTGGAAATGCCATAAAGAG
ATTACAGATTGCACAGCTTTACCTATCTCCATTTTATTAAACCCAAATCCATGGAAAATATCATGGAA
ATGAATAAGCGTCTAACAGAAGAACAAGCCAGAAAAACATTTGAGAGAGCCATGAAACTGGAACAGGAGT
TACTGAACATTTACAGCTATTGTACAGGGGGATACGCTGGAAGACATTTACAACCAAGTGAACAGAT
CATAGAAGAACAATCTGGTTCTTACATCTGGGTTCCGGCAAAAAGAAAAGCTATGA
    
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001290983

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001290983.1</a> , <a href="#">NP_001277912.1</a>
<b>RefSeq Size:</b>	5046 bp
<b>RefSeq ORF:</b>	2715 bp
<b>Locus ID:</b>	1739
<b>UniProt ID:</b>	<a href="#">Q12959</a>
<b>Cytogenetics:</b>	3q29
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	T cell receptor signaling pathway
<b>Gene Summary:</b>	<p>This gene encodes a multi-domain scaffolding protein that is required for normal development. This protein may have a role in septate junction formation, signal transduction, cell proliferation, synaptogenesis and lymphocyte activation. A multitude of transcript variants deriving from alternative splicing and the use of multiple alternate promoter have been observed, including some splice variants that may be specific to brain and other tissues. An upstream uORF may regulate translation at some splice variants of this gene. [provided by RefSeq, Sep 2018]</p> <p>Transcript Variant: This variant (6) differs in the 5' UTR, compared to variant 1. Both variants 1 and 6 encode the same isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>