

## Product datasheet for **SC337660**

### **SALL2 (NM\_001291446) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	SALL2 (NM_001291446) Human Untagged Clone
Tag:	Tag Free
Symbol:	SALL2
Synonyms:	COLB; HSAL2; p150(Sal2); Sal-2; ZNF795
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_001291446, the custom clone sequence may differ by one or more nucleotides

```

ATGTCTCGGCGAAAGCAGCGGAAACCCCAACGTTAATCTCGGACTGCGAAGTCCCAGCGCTCTGAGA
ACGGTGATGCTAGCGAGGAGATCACCCCCAAGTCTGTGCCAAGTCTGCGCACAATTTACTGACCCAAC
TGAATTCCTCGCCACCAGAACGCATGTTCTACTGACCTCCTGTAATGGTGATAATTGGGGCCAGGAG
AACCCCAACAACCTTTCGGCCTCCTCTGAACCCGGCCTGAGGGTCACAATAATCCTCAGGTCATGGACA
CAGAGCATAGCAACCCCCAGATTCTGGGTCTCCGTGCCACCGATCCCACCTGGGGCCAGAGAGGAG
AGGAGAGGAGTCTTCAGGGCATTCTCTGGTCTGCTGCCACAGGGGCAGAAACGCCAAGCAGGCTTCTTC
CACCTTTACCACCCTGGGGTACAGCATCCTTTCTCTGCTGGAGGGTTGGGCGAAGCCACAAACCCA
CCCCTGCCCTTCCCAGCCTTGCAGGCAGCACAGATCAGCTGATTGCCTCGCCTCATCTGGCATTCCC
AAGCACCACGGGACTACTGGCAGCACAGTGTCTGGGGCAGCCGAGGCTTGGAGCCACTGCCTCCCCA
GGGCTCCTGAAGCCAAAGAATGGAAGTGGTGGAGCTGAGCTACGGAGAAGTATGGGTCCCTTGAGAAGC
CTGGTGAAGGCACAAATGCCGCTTCTGTGCCAAAGTATTTGGCAGTGACAGTGCCTGCAGATCCACCT
TCGTTCCACACGGGTGAGAGGCCCTATAAGTGCAATGTCTGTGAAACCGTTTTACCACCCGTGGCAAC
CTCAAAGTGCATTTCCACCGGCATCGTGAGAAGTACCCACATGTGCAGATGAACCCACACCCAGTACCAG
AGCACCTAGACTATGTCATTACCAGCAGTGGCTTGCCTTATGGTATGTCCGTGCCACCAGAGAAGGCCGA
GGAGGAGGCAGCCACTCCAGGTGGAGGGTTGAGCGCAAGCCTCTGGTGGCCTCCACAACAGCACTCAGT
GCCACAGAGAGCCTGACTCTGCTCTCCACCAGTGCAGGCACAGCCACGGCTCCAGGACTCCCTGCTTTCA
ATAAGTTTGTGCTCATGAAAGCAGTGAACCCAAGAATAAAGTGTATGAAAACACCCCCAGGGAGTGA
GGGCTCAGCCATCAGTGGAGTGGCAGAAAGTAGCACGGCAACTCGCATGCAACTAAGTAAAGTTGGTGACT
TCACTACCAAGCTGGGCACTGCTTACCAACCCTTCAAGTCCACTGGCAGCTTCCCCTTCCCCTTGTGC
TAGAGCCCTTGGGGCCCTCACCTCTGAGACATCAAAGCTGCAGCAACTGGTAGAAAAGATTGACGGCA
AGGAGCTGTGGCGGTGACCTCAGCTGCCTCAGGAGCCCCACCACCTCTGCCCTGCACCTTCATCTCA
GCCTCTTCTGGACCTAACCAAGTGTGTCATCTGTCTCCGAGTGCTTAGCTGTCTCGGGCCCTACGCCTTC
ATTATGGCCAACATGGAGGTGAGAGGCCCTTCAAATGCAAAGTGTGTGGCAGAGCCTTCTCCACCAGGGG
TAATCTGCGTGACATTTCTGGGCCACAAGGCCAGTCCAGCTGCCCGGGCACAGAATTCCTGCCCATC
TGCCAGAAGAAGTTCACCAATGTGTCACTCTGCAGCAGCATGTCGGATGCACCTGGGGGGCCAGATCC
CCAACGGTGGTACTGCACTCCCTGAAGGTGGAGGAGCTGCTCAGGAGAATGGCTCCGAGCAATCTACAGT
CTCCGGGGCAGGAGTTTCCCCAGCAGCAGTCCAGCAGCCATCACCGAAGAGGAGTTGTCTGAGGAG
GAGGAAGAGGAGGATGAGGAAGAAGAGGAAGATGTGACTGATGAAGATTCCTGGCAGGGAGAGGCTCAG
AGAGTGGAGGTGAGAAGGCAATATCAGTGAGAGGTGATTCAGAAGAGGCATCTGGGGCAGAGGAGGAGT
GGGGACAGTGGCGGCAGCAGCCACAGCTGGGAAGGAGATGGACAGTAATGAGAAAATACTCAACAGTCT
TCTTTGCCACCACCACCACCTGACAGCCTGGATCAGCCTCAGCCAATGGAGCAGGGAAGCAGTGGTG
TTTTAGGAGGCAAGGAAGAGGGGGCAAACCGGAGAGAAGCTCAAGTCCGGCATCAGCACTACCCACAGA
AGGGGAAGCCACCAGCCTGACCTTGGTAGAGGAGCTGAGCCTGCAGGAGGCAATGAGAAAGGAGCCAGGA
GAGAGCAGCAGCAGAAAGCCTGCGAAGTGTGTGGCCAGGCTTCCCCTCCAGGCAGCTCTGAGGAGC
ATCAGAAGACCCACCCCAAGGAGGGCCGCTCTCACTTGTGTTTTCTGCAGGCAGGCTTCTTGGAGCG
GGCTACCTCAAGAAGCATATGCTCCTGGCACACCACCAGAACCAGTATGTGGCATTCTGTCAAATGGC
CTGCCATGAAGCCCTGGAATTCAGCTCCACCTCCACTACCACTCCAAGCCTGGCCCCACCAGTGTGT
TTGGCCTAGGAAGTGTGGCTGGGAAGTGCCTCCAACAATGGGATCCAGGGAAGCCAAGGAGAAGACAGC
CCCCCTCTATTTAGCCTCCTGCACCAAGGCAGTGCCTGAGAAGCCCATCATAGACAAGAAGTAG
    
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001291446

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001291446.1</a> , <a href="#">NP_001278375.1</a>
<b>RefSeq Size:</b>	3968 bp
<b>RefSeq ORF:</b>	2727 bp
<b>Locus ID:</b>	6297
<b>UniProt ID:</b>	<a href="#">Q9Y467</a>
<b>Cytogenetics:</b>	14q11.2
<b>Protein Families:</b>	Druggable Genome
<b>Gene Summary:</b>	<p>This gene encodes a protein containing multiple zinc finger domains. The encoded protein functions in optical fissure closure during development of the eye in the embryo. Mutations in this gene are associated with ocular coloboma. [provided by RefSeq, Jul 2016]</p> <p>Transcript Variant: This variant (2) lacks two segments in the 3' coding region, which results in a frameshift, compared to variant 1. The encoded isoform (b) has a distinct C-terminus and is shorter than isoform a.</p>