

Product datasheet for **SC337709**

POLR1B (NM_001282774) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	POLR1B (NM_001282774) Human Untagged Clone
Tag:	Tag Free
Symbol:	POLR1B
Synonyms:	A135; RPA2; RPA135; Rpo1-2; TCS4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001282774, the custom clone sequence may differ by one or more nucleotides

```

ATGGATCCTGGCAGCCGGTGGCGGAACCTGCCAGCGGGCTAGCCTAAGCACTTGACTGACCCTCTT
ATGGAATCCC CGGGAACAGCAAAAAGCAGCGTTGCAGGAGCTGACGCGGGCGCAGCTGGAGTCCTCAA
CTACGCTGTGCACGAGGGTCTCGGCCTCGCGGTGCAGGCTATACCTCCCTTTGAATTTGCTTTCAAAGT
GAGCGTATCTCTTTACTATTCTGGATGCTGTTATCAGTCCACCTACAGTTCCAAAAGGGACCATCTGCA
AAGAGGCCAATGTTTATCCAGCAGAATGCCGGGCCGAAGGAGTACCTACCGTGGGAAGTTGACAGCTGA
TATCAACTGGGCAGTGAATGGAATCTCAAAAGGAATCATTAAAGCAGTTTCTTGGCTATGTTCCCATCATG
GTGAAATCCAAGCTTTGCAACTTACGTAACCTTCCCCACAAGCCCTCATTGAGCACCATGAGGAGGCAG
AGGAAATGGGGGGCTATTTATAATCAATGGCATTGAAAAAGTCATCCGAATGTTGATTATGCCTCGGAG
AAATTTTCCATTGCAATGATAAGACAAAATGGAAAACCAGAGGGCCTGGTTTACTCAGTATGGAGTT
TCAATGCACTGTGTGAGGGAAGAACATTCGCTGTCAATATGAACCTCCACTACTGGAAAAATGGCACTG
TTATGTTGAACCTTATTTACCGAAAAGAACTGTTCTTTCTTCTTTGGGATTTGCACTTAAAGCACTTGT
CAGCTTTTCTGATTATCAGATCTTTAGGAGCTCATCAAAGGAAAAGAGGATGATTCTTTCTTAGGAAC
TCTGTTTCTCAGATGTTAAGGATTGTAATGGAAGAGGGTGTTCGACACAAAAACAGGTCCTTAACTACC
TAGGTGAATGCTTCAGAGTAAAACCTCAATGTTCTGACTGGTACCCAAATGAGCAAGCTGCGGAGTTCTT
GTTTAAACCAGTGCATCTGTATCCACTTGAATCCAATACTGAAAAGTTTTATAGCTTTGTCTCATGACG
CGAAAGCTCTTTGCTTTAGCCAAAGGAGAGTGCATGGAGGACAATCCTGATAGTTTGGTGAACCAGGAAG
TCCTCACACCGGGTCAGCTTTCCTTATGTTCTGAAGGGTCACTCCCATTGATGGAGCTCCCCACCGA
TCATACAGTGTGAGAGAGAAAAGAATTCCTCCCTGGATGGAAGTGGTCTTATACCCATGACAGGAAAA
CCAAGTCTGTACCCAGGATTGTTCTTTTTACCCTCTTGTAGACTGGTACGGCCTGTGAGGAACTTAG
CATTGGGCAAAGAAGAGCTAATTGGAACATGGAACAGATCTTCATGAATGTCGCTATCTTTGAGGATGA
AGTTTTTCTGAGGATTACCACACACCAGGAACCTTTCCACACAGCCTGCTGAGTGTGATTGCCAATTC
ATCCCTTCTCTGATCACAACCAGAGTCCACGGAACATGTACCAATGCCAGATGGGTAAGCAAACCTATGG
GCTTTCCACTTCTCACTTATCAAGACCGATCGGATAACAACTGTATCGTCTCAGACTCCTCAGAGTCC
CTTGGTGAGACCCTCCATGTATGATTATTATGACATGGATAACTCCAATTGGGACCAATGCCATCGTT
GCTGTGATTTCTTACACTGGCTATGATATGGAAGATGCCATGATTGTGAATAAGGCCTCTGGGAACGAG
GCTTTGCCATGGAAGTGTCTACAAGTCTGAGTTCATAGACCTCTGAAAAAATTAACAAGGAGATAG
TAGCCTGGTGTGGCATCAAACCTGGTGACCCACGCTTCTGCAGAAGTTAGATGACGATGGATTGCCG
TTTATAGGAGCAAACCTGCAGTACGGAGATCCGTATTACAGCTACCTCAACCTCAACACCGGGAAAGTT
TTGTGATGTAATAAGAGTAAAGAAAATTGTGTTGTGGATAACATCAAAGTGTGCAAGTAAAGCACTGG
GAGTGGAAAATCAAGTGTGTTTGCATCACTATGAGAGTGCCTCGGAACCCAACTATCGGAGATAAATTT
GCCAGTCGCCATGGGCAGAAGGGCATTTAAGCAGATTGTGGCCGGCTGAGGACATGCCTTTTACTGAGA
GTGGGATGGTCCCAGACATCTGTTCAATCCCCATGGTTTTCCATCCCGCATGACCATTGGGATGTTAAT
TGAGAGTATGGCCGGGAAGTCTGCAGCTTTGCATGGTCTCTGCCATGATGCTACACCCTCATCTTCTCA
GAGGAGAACTCGGCCTTAGAATACTTTGGTGAGATGTTAAAGGCTGCTGGCTACAATTTCTATGGCACC
AGAGGTTATATAGTGGCATCAGTGGGCTAGAAGTGGAAAGCAGACATCTTCATAGGAGTGGTTTTATTATCA
GCGCTTACGCCATATGGTCTCAGACAAAATTTCAAGTAAGGACAACCTGGAGCCCGAGACAGAGTACCAAC
CAGCCTATTGGGGGAAGAAATGTCCAGGGTGAATCCGTTTTGGGAGATGGAACGGGATGCGCTTTTAG
CTCATGGTACATCTTTCTCTTCCATGACCGCCTCTTCAACTGCTCAGATCGGTCCGTAGCCCATGTGTG
TGTGAAGTGTGGCAGTTTACTCTCTCACTGTTGGAGAAGCCACCCCTTCTTGGTCTGCCATGCGCAAC
AGAAAATACAACGTACTCTGTGTAGTCGAGTGACACTATCGATACTGTTTCTGTGCCTTATGTTTTTC
GGTATTTGTAGCTGAAGTGGCAGCTATGAACATCAAAGTGAAGTGGATGTTGTTAA
    
```

Restriction Sites: Sgfl-RsrII
ACCN: NM_001282774

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001282774.1 , NP_001269703.1
RefSeq Size:	4878 bp
RefSeq ORF:	2859 bp
Locus ID:	84172
UniProt ID:	Q9H9Y6
Cytogenetics:	2q14.1
Protein Families:	Transcription Factors
Protein Pathways:	Metabolic pathways, Purine metabolism, Pyrimidine metabolism, RNA polymerase
Gene Summary:	<p>Eukaryotic RNA polymerase I (pol I) is responsible for the transcription of ribosomal RNA (rRNA) genes and production of rRNA, the primary component of ribosomes. Pol I is a multisubunit enzyme composed of 6 to 14 polypeptides, depending on the species. Most of the mass of the pol I complex derives from the 2 largest subunits, Rpa1 and Rpa2 in yeast. POLR1B is homologous to Rpa2 (Seither and Grummt, 1996 [PubMed 8921381]).[supplied by OMIM, Mar 2008]</p> <p>Transcript Variant: This variant (4) lacks two alternate exons that result in a short frameshift in the central coding region, and also uses an alternate splice site that restores the reading frame in the remaining coding region, compared to variant 1. The encoded isoform (4) is shorter than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>