

Product datasheet for **SC337719**

KCNH6 (NM_001278919) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KCNH6 (NM_001278919) Human Untagged Clone
Tag:	Tag Free
Symbol:	KCNH6
Synonyms:	ERG-2; ERG2; hERG-2; HERG2; Kv11.2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001278919, the custom clone sequence may differ by one or more nucleotides

```

ATGCCGGTCCGCAGGGGCCACGTCGCTCCCCAAAACACTTACCTGGACACCATCATCCGCAAGTTCGAGG
GCCAAAGTCGGAAGTTCCTGATTGCCAATGCTCAGATGGAGAAGTGCCCATCATTTACTGCAACGACGG
CTTCTGCGAAGTCTTCGGCTACTCCCGAGTGGAGGTGATGCAGCAACCCGCACCTGCGACTTCCTCACA
GGCCCCAACACACCAAGCAGCGCGTGTCCCGCTAGCGCAGGCCCTGCTGGGGGCTGAGGAGTGAAGG
TGGACATCCTCTACTACCGCAAGGATGCCTCCAGCTTCCGCTGCCTGGTAGATGTGGTGCCCGTGAAGAA
CGAGGACGGGGTGTATCATGTTCACTTCAACTTCGAGGACCTGGCCAGCTCCTGGCCAAGTGCAGC
AGCCGCAGCTTGTCCAGCGCTGTTGTCCAGAGCTTCTGGGCTCCGAGGGCTCATGGCAGGCCAG
GCGGACCAGGGCCAGGCACAGGCAGGGGCAAGTACAGGACCATCAGCCAGATCCCACAGTTCACGCTCAA
CTTCGTGGAGTTCAACTTGAGAAGCACCCTCCAGCTCCACCACGGAGATTGAGATCATCGCGCCCAT
AAGGTGGTGGAGCGGACACAGAAGTCACTGAGAAGTCAACCAGTCTGTCCCTGGGCGGGATGTGC
TGCCGGAGTACAAGCTGCAGGCGCGCGCATCCACCGCTGGACCATCCTGCACTACAGCCCTTCAAGGC
CGTGTGGGACTGGCTCATCTGTGCTGTGCTACACGGCTGTCTTACGCCCTACTACGCCCTTCTC
CTGCTCAGCGATCAGGACGAATCAGGGCTGGGGCTGCAGCTATACCTGCAGTCCCCTCACTGTGGTGG
ATCTCATCGTGGACATCATGTTGCTGTTGGACATCGTCACTTCCGCACCACCTATGTCAACACCAA
TGATGAGGTGGTCCAGCCACCCCGCCGATCGCCGTCCTACTTCAAGGGCTGGTTCCTCATTGACATG
GTGGCCGCATCCCTTTCGACCTCTGATCTTCCGCACTGGTCCGATGAGACCACAACCCGATGGGGC
TATTGAAGACAGCGCGGTGTGCGGCTGGTGCAGTGCAGGAAAGTGGACCGCTACTCTGAGTATGG
GGCGGTGTGCTTCTTGTCTATGTGCACCTTCGCGCTCATAGCGCACTGGCTGGCTGGACAGCTGGTAC
GGTCAAGCAATGTGGAGCGGCCCTACCTAGAACAAGATCGGCTGGCTGGACAGCTGGTGGTGCAGC
TTGGCAAGCGCTACAACGGCAGCGACCCAGCCTCGGGCCCTCGGTGCAGGACAAGTATGTCACAGCCCT
CTACTTACCTTTCAGCAGCTCACCAGCGTGGGCTTCGGCAATGTCTCGCCCAACACCAACTCCGAGAAG
GTCTTCTCCATCTGCGTATGCTCATCGGCTCCCTGATGTACGCCAGCATCTTCGGGAACGTGTCCGCGA
TCATCCAGCGCTGACTCGGGCACCAGCGCTACCACACGCAGATGCTGCGTGTCAAGGAGTTCATCCG
CTTCCACCAGATCCCAACCCACTGCGCCAGCGCTGGAGGAGTATTTCCAGCACGCTGGTCTACACC
AATGGCATTGACATGAACGCGGTGCTGAAGGGCTTCCCGAGTGCCTGCAGGCTGACATGCTCCTGCACC
TGCACCGCGCACTGCTGCAGCACTGCCAGCTTTCAGCGCGCCGCAAGGGCTGCCTGCGCGCGTACG
CGTCAAGTTCAAGACCACCCACGCGCCCTGGGGACACGCTGGTGCACCTCGGGACGCTGCTCCTCACC
CTCTACTTTCATCTCCCGAGGCTCCATCGAGATCCTGCGCGACGAGTGGTGGTGGCCATCCTAGGAAAGA
ATGACATCTTTGGGAACCCGTACGCTCCATGCCAGCCAGGCAAGTCCAGTGCAGACGTGCGGGCTCT
GACCTACTGCGACCTGCACAAGATCCAGCGGGCAGATCTGCTGGAGGTGCTGGACATGTACCCGGCCTTT
GCGGAGAGCTTCTGGAGTAAGCTGGAGGTACCTTCAACTGCGGGACGCAGCCGGGGTCTCCTCAT
CCCCCGACAGGCTCCTGGCAGCAAGACCACCAAGGTTTCTTCTCAGTGACAACCAAGTCAAGTGCAGC
CCCTCCCTGAGCATCTCAGATGCATCTGGCTTGGCTGAGCTACTGCAGGAAATGCCCAAGGCAC
AGCCCCAAGCCCTCAGGAAGACCCAGATTGCTGGCTTGAAGTGGGCTCCAGGCTAGAGCAGCTCC
AGGCCAGATGAACAGGCTGGAGTCCCGGCTGCTCAGACCTCAGCCGATCTTGCAGTCTCCTCAGAA
GCCATGCCCGAGGGCCACGCCAGTACATTTCTGGAAGCCCTGCCTCCAATGACCTGGCTTGGTCTCT
ATAGCCTCGGAGACGACGAGTCCAGGGCCAGGCTGCCCGAGGGCTTTCTGCCTCCTGCACAGACCCAA
GCTATGGAGACTTGGATGACTGTAGTCCAAAGCACAGGAACTCCTCCCCAGGATGCCTCACCTGGCTGT
GGCAACGGACAAAACCTTGGCACCATCCTCAGAACAGGAACAGCCTGAGGGGCTCTGGCCACCCCTAGCC
TCACCTTACATCCCCTGGAAGTACAAGGACTCATCTGTGGTCCCTGCTTCTCCTCCCTCCCTGAACACC
TTGGCTGTGTTCCCAAGCAGCTGGACTTCCAGAGACATGGCTCAGATCCTGGATTGCAGGGAGTTGGG
CCACTGA
    
```

Restriction Sites: SgfI-MluI
ACCN: NM_001278919

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001278919.1 , NP_001265848.1
RefSeq Size:	3783 bp
RefSeq ORF:	2877 bp
Locus ID:	81033
UniProt ID:	Q9H252
Cytogenetics:	17q23.3
Protein Families:	Druggable Genome, Ion Channels: Other, Transmembrane
Gene Summary:	<p>Voltage-gated potassium (Kv) channels represent the most complex class of voltage-gated ion channels from both functional and structural standpoints. Their diverse functions include regulating neurotransmitter release, heart rate, insulin secretion, neuronal excitability, epithelial electrolyte transport, smooth muscle contraction, and cell volume. This gene encodes a member of the potassium channel, voltage-gated, subfamily H. This member is a pore-forming (alpha) subunit. Alternative splicing results in multiple transcript variants that encode different isoforms. [provided by RefSeq, Jul 2013]</p> <p>Transcript Variant: This variant (3) lacks an alternate in-frame exon in the 3' coding region, compared to variant 1. It encodes isoform 3, which is shorter, compared to isoform 1.</p>