

Product datasheet for **SC337720**

VPS35L (NM_001300743) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	VPS35L (NM_001300743) Human Untagged Clone
Tag:	Tag Free
Symbol:	VPS35L
Synonyms:	C16orf62; EC97
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001300743, the custom clone sequence may differ by one or more nucleotides

```

ATGGCGCCGGCGGGAGTTGTAACCCGGGCTGTCCGGAGCGGGGAGCTGCCCTCACAAGCATGGCGTCAG
CTGAAAATGAAGCCTGTGCTGTGCGGAGCGTCGCTGCCCTCACAAGCATGGCGTCTGCAGAAGGTGCT
GTGCGGCCGCTGCGGGGCGGCCAGCTGCCCTTCACAAACATGGCGGCCGAGGGGTGCGGGGAGTGGCGGG
GTAAGGATGGGAAGCCGAGCAGACGGCCCCAGAACAAGCGGTATGTGACTGGGAAGATGGCCGTCTTTC
CTTGGAAGTCCAGGAATAGGAAGTACAAAGCTGAATTTGCATCATGCCGACTGGAGGCTGTACCATTGGA
GTTTGGGACTATCACCTCTGAAACCCATAACTGTACAGAGTCAAAGACAAAGAAAGTGAACCGGAAA
GGAAGCACTTCTCCACGTCTCTCTCTCTCCAGCTCCGTGGTGGACCCGCTGAGCAGCGTCTCGATG
GGACTGACCCCTCTCCATGTTTGCAGCCACTGCTGACCCCGCAGCCTTGGCAGCTGCCATGGACAGCTC
CAGAAGGAAACGTGATAGAGATGATAACTCCGTTGTAGGATCGGATTTTGAGCCTTGACCAACAAACGG
GGAGAAATCCTTGCCCGGTACCACTACCGAAAAGCTGTCTATTAATCTGTTTATGGGACTGAAAAAG
GCAAAGCTGGGACTGCCACATTGGCAATGTCAGAGAAGGTGCGGACCCGCTGGAGGAGCTGGATGACTT
TGAGGAGGGTTCCAAAAAGGAGCTGTTGAAGTCTGACTCAGCAGGATTACGTGAACCGCATAGAGGAGCTC
AACCAATCGCTGAAGGATGCCTGGGCTCAGACCAGAAAGTGAAGGCTCTAAAAATAGTCATCCAGTGTT
CAAAGCTTCTTTAGACACCAAGTGTATTCAGTTCTACCAAGCAAAATTTGTCTTATCACCGACATACT
TGATACATTTGGAAGCTCGTGTACGAGCGCATCTTTCCATGTGTGGATAGCCGACGCTCTTACCA
GATCACTTTTCTCCAGAGAATGCAATGACACGGCCAAGGAAACATGCCTAAATTTGGTTTTTCAAGATTG
CCTCCATCAGGGAAGTCAATCCAAGATTTACGTGGAGGCATCCATCCTGAAATGAACAAATTCCTCTC
CAAAACGGGAATTTAGAGTGCCTGCCCCGTTGACATGCATGATCAGAGGGATCGGAGACCCACTAGTG
TCGGTGTATGCCGTGCCTACCTGTGCGGGCTGCTGACCGAGATGATGGAAGGTGTAAGAACTAG
GAAACAATGCCTTGCTGTTGAATTCGTGATGTCTGCCTTCCGGGCTGAGTTCATCGCCACAAGGTCTAT
GGATTTTATTGGCATGATTAAAGAGTGTGATGAATCTGGTTTCCCAAGCATCTTCTTTTTCGATCACTG
GGATTAAGCTTGCCCTTGCTGATCCTCTGAGAGTGACCGACTTCAGATTCTCAACGAAGCTTGGAAG
TCATCACTAAGCTGAAGAAGCCACAGGACTACATTAATTTGTGCCGAAGTGTGGGTGGAATACACCTGCAA
GCATTTACGAAACGAGAGGTGAATACCGTTTGGCAGATGTCATCAAGCACATGACTCCAGATCGTGCA
TTTGAAGATTCTACCCAGCTTCAGTTAATAATTAAGAAAGTATTGCCCACTTCCATGACTTCTCAG
TTCTTTTCTCAGTGAAAAATTTCTGCCGTTTCTGGACATGTTCCAAAAAGAGAGTGTGCGGGTGGAGGT
TTGCAAAATGCATCATGGACGCCTTTATCAAGCATCAACAAGAGCCACCAAGACCCGGTCATCTTGAAT
GCCCTTTTGCATGTTTGCAAGACCATGCATGACTCTGTGAATGCACTCACTCTTGAGGATGAGAAAAGAA
TGCTGTCATATTTGATTAATGGATTTATAAAAAATGGTTTCCTTTGGCCGTGATTTGAACAACAGCTGAG
TTTTATGTTGAGTCCAGGTGATGTTTGAATCTGGAGCCTGTTCTTGTGAGTTGATTCATGCCTGT
GTTGCCTACTGCTTCATACCATCCCCTCCCTGGCGGGCATCTTCACACGTCTCAATCTCTACCTGCATT
CTGGTCAGGTGGCCTTGCCCAACAGTGCTCTCCCAAGCTGATGCTTTTTTCAAAGCCGCTATAAGCCT
TGTTCCGGAAGTCCAAAGATGATTAATATTGATGGGAAGATGCGGCCATCGGAATCGTTCTCTTCTGGAA
TTCCTCTGCAATTTCTTTTCTACTTTATTAATAGTTCGGATCATCCTGAACATGGGGTCTGTTTCTTG
TTCGAGAGCTTCTCAACGTGATCCAGGACTACCTGGGAGGACAACAGCGATGAGAAAAATCCGCATCTA
CACCTGCGTCTGCATCTCCTCTCCGCCATGAGCCAGGAGACGTACCTTTACCACATAGACAAAGTGGAC
TCCAACGACAGCCTCTACGGGGGAGACTCCAAGTTCCTGGCAGAAAAACAACAAGCTGTGTGAGACGGTGA
TGGCTCAGATCTAGAGCATCTGAAAACCTGGCCAAGGACGAGGCCCTGAAGCGCCAGAGCTCGTTGGG
CCTTTCTTTTAAACAGCATCTTGGCCATGGGGACCTACGCAACAACAAGCTCAACAGCTCTCCGTC
AACCTGTGGCACCTGGCAGAGGACGGCTGTGCAGACACCAGGACCATGGTGAACCGCTAGAATACA
TCAAGAAGCAAAGCAAACAACAGACATGACTCATCTGACGGAGCTGGCCCTCAGACTCCCTCTGCAAAAC
AAGGACCTGA
    
```

Restriction Sites: SgfI-MluI
ACCN: NM_001300743

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001300743.1, NP_001287672.1</u>
RefSeq Size:	3615 bp
RefSeq ORF:	2880 bp
Locus ID:	57020
Cytogenetics:	16p12.3
Gene Summary:	Acts as component of the retriever complex. The retriever complex is a heterotrimeric complex related to retromer cargo-selective complex (CSC) and essential for retromer-independent retrieval and recycling of numerous cargos such as integrin alpha-5/beta-1 (ITGA5:ITGB1) (PubMed:28892079). The recruitment of the retriever complex to the endosomal membrane involves CCC and WASH complexes (PubMed:28892079). In the endosomes, drives the retrieval and recycling of NxxY-motif-containing cargo proteins by coupling to SNX17, a cargo essential for the homeostatic maintenance of numerous cell surface proteins associated with processes that include cell migration, cell adhesion, nutrient supply and cell signaling (PubMed:28892079). Involved in copper-dependent ATP7A trafficking between the trans-Golgi network and vesicles in the cell periphery; the function is proposed to depend on its association with the CCC complex and cooperation with the WASH complex on early endosomes. Seems not to be required for CCC complex stability (PubMed:25355947). [UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) lacks alternate in-frame exons in the coding region, compared to variant 1. It encodes isoform 2 which is shorter than isoform 1.