

Product datasheet for **SC400001**

MIRLET7A1 Human MicroRNA Expression Plasmid (MI0000060)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Human
Vector: pCMV-MIR (PCMVMIR)
Precursor: UGGGAUGAGGUAGUAGGUUGUUAUAGUUUAGGGUCACACCCACCACUGGGAGAUAAUAUACAAUCUACUGUCUUUCCUA
Mature Sequence: UGAGGUAGUAGGUUGUUAUAGUU

Full insert: >MI0000060 GTGGTAAGAGGGTGATTTCTGAATCTTGTAATACATGGTTTTAGGAGCGGATTCAGATAACCAAGCATT
TAAAATACTATTAATGAAATACAGGAAATGAAACCACAGCATAGATTATGCATGTAGCCAAAATGTTTCAG
TAAACTTCATTTTCAACGTAAGTGAATGAAAATGGTCTAATACTATTTTTCTTATCACTCACACAGGAA
ACGAGGATTACCGAGGAGGAAAAAAGCCTTCTGTGGTGCTCAACTGTGATTCTTTTACCATTACC
CTGGATGTTCTTCACTGTGGGATGAGGTAGTAGGTTGTATAGTTTTAGGGTCACACCCACCACTGGGA
GATAACTATACAATCTACTGTCTTTCTAACGTGATAGAAAAGTCTGCATCCAGGCGGTCTGATAGAAAG
TCAGTTAACTAATTGTACAATATTTAAGATTAAGTCTGTTAAAGAGATGTAGTGCAGCATTTGTTTATG
GCCTGGAAATAAATTAATTTAGAGATAAAGTCTGTAGCAAGTACACTGGATGGGGGTGGGAAACCTTTT
GCTTCTTGTCTTATTTCTGTGTCAGAATAAATGATTTTTTTATTTTGATTATGCTGATAATTTTAT
GTTGAAATTTCTTTTCGAAAGAGATTGTAATTTCCATTCCAGAAG

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0000060](#)

Synonyms: LET7A1 let-7a-1 MIRNLET7A1

Locus ID: 406881



[View online »](#)