

## Product datasheet for **SC400003**

### MIRLET7A3 Human MicroRNA Expression Plasmid (MI0000062)

#### Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVIR)
Precursor:	GGGUGAGGUAGUAGGUUGUUAUAGUUUGGGGUCUGCCCGCUAUGGGUAUACUAUACAUCUACUGUCUUUCCU
Mature Sequence:	UGAGGUAGUAGGUUGUUAUAGUU
Full insert:	>MI0000062 CCTCCTTCTCTCTTCCCTCGTCCCTCTCTGGTTGGTGGTTTCCTGCAGGCTCCCGTCCCTGCTGGTGCTG GCCACAATGTCCCACTCCAGGGTTTCGGCGTCCCAGCCCCCTGCGCCACCGCGCTGCCCGCCAGAA TCCCTGTGCCCTTGGTGCCTGGCCCTGCCGAGCCTCGAGCCCCTGTTCTCCTCAGCCCTCTTCCCTCCC GCGTCCCAGGAGGTGCCTCTGGAAGCCACGGAGTCCCATCGGCACCAAGACCGACTGCCTTTGGGGTG AGGTAGTAGGTTGTATAGTTTGGGGCTCTGCCCTGCTATGGGATAACTATACAATCTACTGTCTTTCCTG AAGTGGCTGTAATATCTGCGGTGGACAGAGCGTCTGGAACCCTGGCTGGGAGCGGGCAGGGCCAGGTTTG GGGCAGCCTTGGCAGCAGTCCGGGGCAGGGGCCGCTACACTGAGAAGTCTGACAGGCCTAGGTGCCAC TTGCTGTGTGACCTTGGACAGGCCCTGATCTCTCTGGGTCTCAGTTTCTCCTCTGTAATAATGGAGGCA AATGAGGATGGAAGGAGATGCAGTGTGGAGCATCGAGGGCAGAGGAGAGCTGAGCCGACCC
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	<a href="#">MI0000062</a>
Synonyms:	LET7A3 let-7a-3 MIRNLET7A3
Locus ID:	406883



[View online »](#)