

Product datasheet for **SC400004**

MIRLET7B Human MicroRNA Expression Plasmid (MI000063)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	CGGGGUGAGGUAGUAGGUUGUGUGGUUUUCAGGGCAGUGAUGUUGCCCCUCGGAAGAUACUAUACAACCUACUGCCUUCCUG
Mature Sequence:	UGAGGUAGUAGGUUGUGUGUU
Full insert:	>MI000063 AAGGCACCTTCTCTGAGCAGGAAGTGAGAGGAACAGCTCTGCATACACTGGGTCCCACATGGCACAATCT GAAGGCAGACAGTGGCTCCTCTGTACCTGGGGAACTGAGGCCAGAGAGCCAGGGACTTCCCAAGACCA GCCAGCAGCAGCTGCCCTTCTGGGGTGCCATCTCCCTGTCCCTGCCCTGCGCCTGCCAGCCCT CCTGCTCTGGTGACTGAGGACCGCCAGGCAGGGGCTGGTGTGGGCGGGGGCGGGCCCTCCCGCAG TGCAAGCCGGCCTGGCGGGTGAGGTAGTAGGTTGTGTGGTTTCAGGGCAGTGATGTTGCCCTCGGA AGATAACTATAACAACCTACTGCCTTCCCTGAGGAGCCAGTGACACGACCCCATGGGAGGGCCGCCCCCT ACCTCAGTGACACGACCCACGGGAGGGCTGCCCCACCTCAGTGACCTGCAGGGGGCCTGAGCCGAAG CTGGGTGGGCATCTGGGAGCTAGATTCAATAAAGCTGTCTGACCATGAACTTGAAGTGGCCCTTTCA TTCTGGAAGCCAGGGACTGGGCTCTGGGCGCGTCTGCTGCGGGACAGTGTGGAGGCCCTGGCCACGA GCTGTTCCAGGTCATGAGACTTGTCTGCTC
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI000063
Synonyms:	LET7B let-7b MIRNLET7B hsa-let-7b
Locus ID:	406884



[View online »](#)