

Product datasheet for **SC400008**

MIRLET7f-1 Human MicroRNA Expression Plasmid (MI000067)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	UCAGAGUGAGGUAGUAGAUUGUAGUAGUUGUGGGUAGUGAUUUUACCCUGUUCAGGAGUAACUUAUCAAUCUUAUUGCCUCCUGA
Mature Sequence:	UGAGGUAGUAGAUUGUAGUU
Full insert:	>MI000067 TCTGCATCCAGCGGTCTGATAGAAAGTCAGTAACTAATTGTACAATATTTAAGATTAACCTTGCTTAA AGAGATGTAGTGCAGCATTTGTTTATGGCCTGGAAATAAATTAATTTAGAGATAAAGTCTGTAGCAAGTA CACTGGATGGGGTGGGAAACCTTTTGCTTCTGTCTTATTTCTCTGTGTCAGAATAAATGATTTTTT TATTTTGATTTATGCTGATAATTTTATGTTGAAATTTCTTTTCGAAAGAGATTGTACTTTCCATTCCAGA AGAAAACATTGCTCTATCAGAGTGAGGTAGTAGATTGTATAGTTGTGGGGTAGTGATTTTACCCTGTTCA GGAGATAACTATACAATCTATTGCCTTCCTGAGGAGTAGACTTGCTGCATTATTTCTTTTTATTTAGA TGATATTAAACTCAGAAGAATTAATTTTGACATTTTGTATTTACAGTTTATCAGTTAATTTCTCTGTT CAAGTAGTACAGTAGGCACAGATTAACATTTAAATTTTTCACATATGGTATATTTTCAGAAATTTGAAGTT AAGCAAAAATTTTAAATGAGTAGAGAAAGTAAGTAGCCTTCAGGAAATCTTCATAGAGGACCAGGCCCTTT TGGAATTGTGAATAGGTTTATTGCCTTACATCCTGGTACACATGTCCAAGG
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI000067
Synonyms:	LET7F1 let-7f-1 MIRNLET7F1
Locus ID:	406888



[View online »](#)