

Product datasheet for **SC400047**

MIR122 Human MicroRNA Expression Plasmid (MI0000442)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Human

Vector: pCMV-MIR (PCMVMIR)

Precursor: CCUUAGCAGAGCUGUGGAGUGGACAAUGGUGUUUGUGUCUAAACUAUCAAAACGCCAUUAUCACACUAAAUAGCUACUGCUAGGC

Mature Sequence: UGGAGUGUGACAAUGGUGUUUG

Full insert: >MI0000442

GCAAACAGAGTTCCTGTCCATTTTCACAGGAAGTTTTGCAGAAGGGGTTTCTGGGTGGCAGTGAACACAGGCAGTGGACCCCTCTACTACTCTACTGTTCCAA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The

E. coli Selection:

Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0000442](#)

Synonyms: MIR122A MIRN122 MIRN122A miRNA122 miRNA122A hsa-mir-122

Locus ID: 406906



[View online »](#)