

Product datasheet for **SC400053**

MIR1229 Human MicroRNA Expression Plasmid (MI0006319)

Product data:

| | |
|---------------------------|--|
| Product Type: | MicroRNA Expression Plasmids |
| Species: | Human |
| Vector: | pCMV-MIR (PCMVIR) |
| Precursor: | GUGGGUAGGGUUUGGGGAGAGCGUGGGCUGGGUUCAGGGACACCCUCACACCUGCCUCCACAG |
| Mature Sequence: | CUCUCACCACUGCCUCCACAG |
| Full insert: | >MI0006319 GTGCTGTGCTGGAACGCTGCCCTGCAGCAGCTCACAGGCTGCCTGGCTGCAGTTCTTCTCCGCAGT GGGAACATCGAGCACCCGAGGACAAGCTCTTCAACACGTCTGTGGAGGTGCTGCCCTTCGACGTAAGTT CGGTGGGTGTGGTGTGCACTGAAGCCGGGCTCCTGTGACCACCTGCCTCAGCCCTTGATTCTGCCCCCA GAACCCTCAGTCAGACAAGGAGGCCCTGCAGGAGGGCCGACCCGCCACCTCCGGTACCCTCGGAGCCCC GACGGTACCTCCAGATCGGTGGGTAGGGTTTGGGGGAGAGCGTGGGCTGGGTTTCAGGGACACCCTCTC ACCACTGCCCTCCCACAGGCTCCTTCTACAAGGGAGTGGCAGAGGGAGAGGTGGACCCAGCCTTCGGCCC TCTGGAAGCACTGCGCTCTCGATCCAGACGGACTCCCCTGTGTGGGTGATTCTGAGCGAGGTGAGGCTG GGCAGGGCCAGCCAGGCACGGGAGCAGCCAGGGAGACCCTGGCGGAGCTCAGAGGCCCTGCTGACCCT GGTCCGGCTCTTCCAGATCTTCTGAAAAAGCCGACTAAGCTGCGGGCTTCTGAGGGTACCCTGTGGC CAGCCCTGAAGCCCACATTTCTG |
| OTI Disclaimer: | All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA. |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Mammalian Cell Selection: | Neomycin |
| Components: | 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube |
| Also Contains: | Both mature and minor (star) miRNA can be over expressed from the construct. |
| RefSeq: | MI0006319 |



[View online »](#)