

## Product datasheet for **SC400083**

### MIR125B-1 Human MicroRNA Expression Plasmid (MI0000446)

#### Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	UGCGCUCCUCUCAGUCCUCGAGACCCUAACUUGUGAUGUUUACCGUUUAAAUCCACGGUUAGGCUCUUGGGAGCUGCGAGUCGUGCU
Mature Sequence:	UCCUGAGACCCUAACUUGUGA
Full insert:	>MI0000446 TCCGTTTTAGTGAAC TATTTGCCCAAATCTTGAAAGAGTTTTCTCTGATGTACTCGTGATCGTATGTGT ATGTGCGTGATTGTATATGCGCCCCAGATACTGCGTATGTGTGTATATGTATTAGGCTTAAACGGAA TCTCAATTTTGTGAAGGAAAGGAGCTTAGAGAAGAAATACCATACCACCTGTTTGTTCATCTTAGTTAT GAACCTCGAACAGAAATTGCCTGTCATTCTTGTGTTTGTCTTGTCTCAAGAAAGAAAACATTGTTG CGCTCCTCTCAGTCCCTGAGACCCTAACTTGTGATGTTTACCGTTTAAATCCACGGTTAGGCTCTTGGG AGCTGCGAGTCGTGCTTTGCATCCTGGAATTTGGTGAATTTTATTCTTTAAAGCAAAAACAAAAGAA AAGAAAGTTTGTCTGAGGTGATTGAGTATACCTCTGAGGTTTTTATTGTTAGATGGGATCAGGTGACCAG AGAGTGGCAGCTCCTGGATTTCTGTTAAAGTAGGTATATATTTGCTTTAATAGCGTGTGGCTGCTATG TGTATGTATGCATGTTATTGATGCCTGAGGCTTGTGGAACCTCATCTGCCTGGTCTATTCTTGAGAGCAA CACGCAGTAGT
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	<a href="#">MI0000446</a>
Synonyms:	MIRN125B1
Locus ID:	406911



[View online »](#)