

Product datasheet for **SC400110**

MIR128-1 Human MicroRNA Expression Plasmid (MI0000447)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	UGAGCUGUUGGAUUCGGGGCCGUAGCACUGUCUGAGAGUUUACAUUUCACAGUGAACCGGUCUCUUUUUCAGCUCUUC
Mature Sequence:	UCACAGUGAACCGGUCUCUUU
Full insert:	>MI0000447 CTGGAGTCAATGAAAGCAACTATTTCAAAGATCAGATTACTTACCAGTTTCACTAATAAAGATTTATTA CTTTAAACCTTTATCATAAAATGTATGCTTTGAATACTGTGAAGTACACTGCATATAAGGAGTGTGGTAT AGTATAAAGAACTTTCTGCAGGTAGTAATTATAGTGAAGATTTAGGTTTACAAAGCCCTAGCTGTTTT CTGTGTAGCTTTTATTATTCTTATGACTCTTGACAAGTTTGTAGCTTCACCATATACATTTAATATTTTG CAATAATTGGCCTTGTTCTGAGCTGTTGGATTGCGGGCCGTAGCACTGTCTGAGAGGTTTACATTTCTC ACAGTGAACCGGCTCTTTTTTCAGCTGCTTCTGGCTTCTTTTTACTCAGGTTTCCACTGCTTTTTTGTCT TTTTTAAATGCTGTATGAAGGTGTTAACATTTGTTTATTTTTTCATTAATTGTAATACCTTTAAATCAT GCATCATACTCAGAAATAGGGATTAGAATTTAAGTGACATCTTTGGCCTAATATAATTTACCTGTTAAAA ATTTGTGAAAGCTATTGCTTATTTCTTTTCAAAGTAGATTTGGATGACAAATTCAGTTTGGTTTTGAGC TTATCTAAAAAGCTTAATGTTTTAAGCATACCAAACCTCATTAGCAGAAG
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000447
Synonyms:	MIR128A MIRN128A MIRN128-1
Locus ID:	406915



[View online »](#)