

Product datasheet for **SC400121**

MIR1289-1 Human MicroRNA Expression Plasmid (MI0006350)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Human

Vector: pCMV-MIR (PCMVMIR)

Precursor: UUCUCAUUUUUAGUAGGAUUUAAAAACAAAACUGGUAAAUGCAGACUCUUGGUUUCCACCCCCAGAGAAUCCUAAACCGGGGUGGAGUCCAGGAUUCUGCAU

Mature UGGAGUCCAGGAUUCUGCAUUUU

Sequence:

Full insert: >MI0006350 CAAGTAATTCTCTGCCTCAGCCTCCCGAGTAGCTGGGACTACAGGCGCCGCCACCACGCGCCGGCTAAT
TTTTGTATTTTAGTAGAGACAGGATTTACCATGTTGGCCAGGATGGTCTCGATCTCTTGTACCTCGT
GATCCTCCCGCCTCGGCCTCCCAAAGGGCTGGGATTACAGGCGTGAGCCACCGCGCCCGGCCTTGTGTAT
AAATTGTTAATGCATTTCTTAGAGTTGTAACTCCATTTACCTGCCTCACTCTTGCCTCAGTGTTCT
CAATTTTGTAGTAGGAATTAACAACTGGTAAATGCAGACTCTTGGTTTCCACCCCAGAGAAATCCC
TAAACCGGGGTGGAGTCCAGGAATCTGCATTTAGAAAGTACCCAGGGTGATTCTGATAATTGGGAACA
CGTTTGTAGGAAGCGCTGCATGAGATTGTGAGCTTGAGCGTAGAGACTGTAATTTGGTCATCTGACTATAG
CGCTCAGCACAAAGGCCTGGTACAAAATAGGTGCTCAATAAGTGCTGAGTGAATGGGTGAATGAGTGGATA
GCTCGAGTGAAACAGCAAGGCTTCAGGTTAATCCACAGCCAAGGTAACCTGAGAAAGGAACGCTTAAAA
GCTATATTTTATCTTCTCCTGTAAGAACACAAAAAGATCTTAATAGCTAATTTGAGAAATAGAGAAAG
TGGAGAGTAAAAGGGGAAG

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The

Disclaimer: may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of

E. coli Kanamycin (25 ug/mL)

Selection:

Mammalian Neomycin

Cell

Selection:

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Both mature and minor (star) miRNA can be over expressed from the construct.

Contains:

RefSeq: [MI0006350](#)



[View online »](#)