

Product datasheet for **SC400146**

MIR1304 Human MicroRNA Expression Plasmid (MI0006371)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Human

Vector: pCMV-MIR (PCMVMIR)

Precursor: AAACACUUGAGCCCAGCGGUUUGAGGCCUACAGUGAGAUGUGAUCCUGCCACAUCUCACUGUAGCAUCGAACCCCUGGGCUCAAGUGAUUCA

Mature Sequence: UUUGAGGCUACAGUGAGAUGUG

Full insert: >MI0006371

TTAAATCTTTTTCTTTGGTTCATGTAAGGAGCTACTTGATGTTACTTTTTATTTATCATTCTCTGGCTCCTAAGCACACGTGATTGTGAAATTATTGTTTT

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The

E. coli Selection:

Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0006371](#)



[View online >](#)