

Product datasheet for **SC400168**

MIR138-2 Human MicroRNA Expression Plasmid (MI0000455)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	CGUUGCUGCAGCUGGUGUUGUAAUCAGGCCGACGAGCAGCGCAUCCUCUUACCCGGCUAUUUCACGACACCAGGGUUGCAUCA
Mature Sequence:	AGCUGGUGUUGUGAAUCAGGCCG
Full insert:	>MI0000455 GTTGCACTCCAGCCTCTGGCTTCCCAGGCCTTCCCAGTGCTCAGTACCTCCATTGTTTCTCTGCACCCCC TGCATTGTACAGGTTTTTTTTTATGTCTGTTCCCTCTTCCCATGAACATCAAGTCTGTGAGCCCTT GTGGGCAGGGACCCAGATTCCACCATATGCACCCTGAGCCCTGGTGCCATACTGGCCTGGGGGAGCCAT TGGTTGGGGAAGGTGTTGTATGGGTTCACTGTCTCAAACCGGAGGAGCCACAGCCAGAGGAGGAA GCCGGCGGAGTTCTGGTATCGTTGCTGCAGCTGGTGTGTGAATCAGGCCGACGAGCAGCGCATCCTCTT ACCCGGCTATTTACGACACCAGGGTTGCATCATACCCATCCTCTCCAGGCGAGCCTCGTGGGACCCGGC CGGACGCAAGGACCTGGGGATGGGGTGGGGTGGGAGTCTGTGTGTGAAAGCTTGTCTCTGTAGCTCAC TGGATCTCAGCCAGACTGCTCCTCAGAATCACCTGGCAGCTTTCAAAAACTCTGATGTCTGCCCTCA CCTCTCTTGCCGAGGTCCACAGAATCTCAAGGTGGGGCCCTGGAATCAGCATTGTTAGAGGTTCCAGT GTGATGCCGAAGTCTGTGAGTAGTCAAGG
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000455
Synonyms:	MIRN138-2
Locus ID:	406930



[View online »](#)