

## Product datasheet for **SC400179**

### **MIR146B Human MicroRNA Expression Plasmid (MI0003129)**

#### **Product data:**

<b>Product Type:</b>	MicroRNA Expression Plasmids
<b>Species:</b>	Human
<b>Vector:</b>	pCMV-MIR (PCMV MIR)
<b>Precursor:</b>	CCUGGCACUGAGAACUGAAUCCAUAAGGUGAGCUCUAGCAAUGCCCUGUGGACUCAGUUCUGGUGCCCGG
<b>Mature Sequence:</b>	UGAGAACUGAAUCCAUAAGGCU
<b>Full insert:</b>	>MI0003129 AACACATCTGGGCAGGAAGCGACAGATGGATCCCTGCCATGGCCCTGAGCCCCCTCCCCTCTTAACCTAC TCCAATTCCTGTTGGAATTTATGGAACCTCTGTCCATTTCTTAGCTGTGGGCCCCCATTTTCAGCTCCA AGCTCAGACCCTCCCTGGAATAGGAGTTCTCTTGGTATCACAAGTTCCTTATGTGACCCATCCTGGGCT CAACTTACTCATCCTGGGAACGGGAGACGATTACAGAAGAAAGCATGCAAGAGCAGCGTCCAGGCTGAA AGAAGTTTGGCCACCTGGCACTGAGAACTGAATCCATAGGCTGTGAGCTCTAGCAATGCCCTGTGGACT CAGTTCTGGTGGCCGGCAGTGTACAACATCAATGCCAAGGCCGTGGGGCAGCTGATGGTTTGGGCTCCC AACTTCCCAGCCAGGTGCTTCTGACGGCCACATCTTGCCCACTGGCCAAACCTTTAAATAACTTTGACT CGGGCTACTCTTATGCTCAAAGACGTCAGGGGCTCTCCAAATCTCTTACCCTGCCAGAAAGTCTTCTA TAGTACGGCCTCCACTTAGCTTTCACGCCTGATCTTCCATCGCATCCTGCTCATAACCTGCCACTAGTGA ACCCCTGCTGCTCCGGCTCCACACTTCTCATGCTG
<b>OTI Disclaimer:</b>	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
<b>E. coli Selection:</b>	Kanamycin (25 ug/mL)
<b>Mammalian Cell Selection:</b>	Neomycin
<b>Components:</b>	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
<b>Also Contains:</b>	Both mature and minor (star) miRNA can be over expressed from the construct.
<b>RefSeq:</b>	<a href="#">MI0003129</a>



[View online »](#)