

## Product datasheet for **SC400224**

### MIR191 Human MicroRNA Expression Plasmid (MI0000465)

#### Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVIR)
Precursor:	CGGCUGGACAGCGGGCAACGGAAUCCCAAAAGCAGCUGUUGUCUCCAGAGCAUCCAGCUGCGCUUGGAUUUCGUCCCCUGCUCUCCUGCCU
Mature Sequence:	CAACGGAAUCCCAAAAGCAGCUG
Full insert:	>MI0000465 TCCTCCTGACGGATTCCGAAGCCGCCACCACCGTCAACGCACCAGTGGCCAGGCTTCGTGGGACC AACAGGGCTAATCCTGCCCCACGGATCTAATCCGGATCGGGACCCGGTCCCAGGGGAGGCCCCACCTCC CCTTCTATTGTTCTGTGGCCAGGTGAGCGTCTTGTCCCTCTAGACTCCGTTTCAACCTACTCCCGG GTCTTCTCCTGGGACTGGGGATGGGCGGGTTCTGATGTGGCCCCAGGGCGAGTGACCTGGGGGAGGAG CTCCCCGCCCCCGCAACGGCTGGACAGCGGGCAACGGAATCCCAAAAGCAGCTGTTGTCTCCAGAGC ATTCCAGCTGCGCTTGGATTTCTGCTCCCTGCTCTCCTGCCTGAGCAGCGCCCTGGCCAGATGGGGTGCC CCTGACCCCCAGACATACTTTACTGAGCTGCTTGGGTCTCAGTTCCTCTCAGTTGCGCCCTCAGGCTGGA GGTGATGGGTGTAGACGTGGGAGAGCCGAGGCTGGTGGCCATCTCTGGAACCCTGGGGAGGATTGGCGAG GGAGGGTGGACCCAGGACCTCTGGGTAGGGCTGCAATGGTAGTACTCCCCAGGGCTGCTCAGCCCGC ATCTCGCTAAAGATAGCCCTCTTACTTCTGGGAAAAACCCCAACCGTCTTCTGT
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	<a href="#">MI0000465</a>
Synonyms:	MIRN191 miR-191
Locus ID:	406966



[View online »](#)