

## Product datasheet for **SC400229**

### MIR1914 Human MicroRNA Expression Plasmid (MI0008335)

#### Product data:

**Product Type:** MicroRNA Expression Plasmids  
**Species:** Human  
**Vector:** pCMV-MIR (PCMVMIR)  
**Precursor:** CGUGUGAGCCCGCCUGUGCCCGGCCACUUCUGCUUCCUUAAGCGCAGGAGGGGUCCGCACUGGGAGGGGCCUCAC  
**Mature Sequence:** CCCUGUGCCCGGCCACUUCUG  
**Full insert:** >MI0008335 TCCTCCCTGTCCACCTCTTTCTTGCCTCCACCCGCCCTTTCTGCCCCTCCTAGACCTGCTGCTGC  
CCCGGCCATCCCAGGCTTGGGCAGCGGGCGGGACGGGACAGCCAAGAGACTGAGCCAGCTTTGCC  
CCCACCCCTCTCTTGGGGCCGCGCACCACCTCGGGGAGGAGCCAGGGCGGGCCTGACCCACAAA  
GCTTTAGTTCCCCGGGGAGAGGCTGGGGACAAGCTACCCATTCTGGGTATGACACCTATCAGCACGTG  
TGAGCCCGCCCTGTGCCCGGCCACTTCTGCTTCTTAGCGCAGGAGGGGTCCCGCACTGGGAGGGG  
CCTCACAGCTACCCGCTAGCTTTCTGCCTGATGGTGAGCTTCCCCGTTGTTGCCCGCATCTTTGCG  
CTTTGCCTGACAGCTCAGGTAGGGCGACGTCTCCCCAGCTCAGCCCAGGCAGAAGTCAGGCCGCTC  
TGCACCTCGTGTGCTGGCCCGTGTGTTCCCTCACAGTGACGCACAGCCGCTCTGTGGGAGGGAG  
GGCTCTGGGAAAAGACCCCTTCTTGTGTTCCCATCGGTTTGTCCCAGGGCTGCGCTGGCCTCGGC  
ACACCAGTGCCACCCGCTGC

**OTI Disclaimer:** All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

**E. coli Selection:** Kanamycin (25 ug/mL)

**Mammalian Cell Selection:** Neomycin

**Components:** 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

**Also Contains:** Both mature and minor (star) miRNA can be over expressed from the construct.

**RefSeq:** [MI0008335](#)



[View online »](#)