

## Product datasheet for **SC400238**

### MIR196A2 Human MicroRNA Expression Plasmid (MI000279)

#### Product data:

<b>Product</b>	MicroRNA Expression Plasmids
<b>Type:</b>	
<b>Species:</b>	Human
<b>Vector:</b>	pCMV-MIR (PCMVMIR)
<b>Precursor:</b>	UGCUCGCUCAGCUGAUCUGUGGCUUAGGUAGUUUCAUGUUGUUGGGAUUGAGUUUUGAACUCGGCAACAAGAAACUGCCUGAGUUACAUCAGUCGGUUUUCGUCG
<b>Mature Sequence:</b>	UAGGUAGUUUCAUGUUGUUGGG
<b>Full insert:</b>	>MI000279 GCCAGATGCAGGGGAATAGCTCCACTGCGGTGCTGGGGTGGGGCAGGTATTTGGCCGCCAGTGGTCCCA TTTCACCAGATTTTCCTGAGCAATTTGCTTAGACCAGAACTAATGAAGGCAGTCAGACCCCTTACCCA CCCAGCAACCCAAAGTCTACTCTCTAGTCCTTAGGGAGGTTGTGGGGCGGAAAGGGGACGGGGCTGAA TTTCTTCCTTCCCAACCCCTTCCCTTCTCCTCCAGATAGATGCAAAGCTGAATCTCCCGCCCTGCTCG CTCAGCTGATCTGTGGCTTAGGTAGTTTCATGTTGTTGGGATTGAGTTTTGAACTCGGCAACAAGAACT GCCTGAGTTACATCAGTCGGTTTTTCGTGAGGGCCCCAACCCACCTCTCCACTCCTACCTCCCCAGTG GGACTGCCCCACTGCCCCCTCCAGATAGGGCAAAGTGGGTGCAGACCAAGGAGGACAAGCTGTGAGTGG GGTTGCAGACAAGTCTGGAGAACCCTGCTTTATGCCGTCTCTCAACCTATCCCAATCCCTGTGATTAA GGGAGAGGAGTCTTTCCTGCAGCTATTGTCTCCTGGGCCAGCTGATGGAGACAGACAAAGAGGGTCCAG CCAATTCCAAGGAAATTCAGGAGTTTGAGTAAGAGGTGGC
<b>OTI</b>	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence <a href="#">imiRBase</a> . The
<b>Disclaimer:</b>	flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
<b>E. coli Selection:</b>	Kanamycin (25 ug/mL)
<b>Mammalian Cell Selection:</b>	Neomycin
<b>Components:</b>	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
<b>Also Contains:</b>	Both mature and minor (star) miRNA can be over expressed from the construct.
<b>RefSeq:</b>	<a href="#">MI000279</a>
<b>Synonyms:</b>	MIRN196-2 MIRN196A2
<b>Locus ID:</b>	406973



[View online »](#)