

Product datasheet for **SC400250**

MIR199A1 Human MicroRNA Expression Plasmid (MI0000242)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMV MIR)
Precursor:	GCCAACCCAGUGUUCAGACUACCUUGUUCAGGAGGCUCUCA AUGUGUACAGUAGUCUGCACAUUGGUUAGGC
Mature Sequence:	CCCAGUGUUCAGACUACCUUGUUC
Full insert:	>MI0000242 GACCCCCAAAGAGTCAGACATTCTCCTGAGCCAGAAGCCACGATCCCAAACCTGCCTCCTGCTCCGC CTCCCCACTCTTTAGGATTTCTGAAAACCCAGAATTTCTCCAGATGCGAGCCGGCCCCAGCCTTGCC ACGTCAGAAGGGACAGAGCGGATCGTCTCGGGAAGAGTGGTGGTTTCTTGCTGCTCAGAGGTGCTGAG GCACGGCCTGGCCTGGTGGCCCCAGCGTCTGCCTGGGGGTTCTGCAGGATGGATAGCCGGCCCCGCCAA CCCAGTGTTCAGACTACTGTTCCAGGAGGCTCTCAATGTGTACAGTAGTCTGCACATTGGTTAGGCTGGG CTTGGGTGAGCGGCTCGTCGAGACAGGCCCCCCAAACTCGCCGGCAGGTGAGTGTCAATTTCCACCACCC CGTTTCCCACTGTGGCAGAGCCTCGCATAGAAGATTCGAGGGCCTGGGTGGGAGAAGAGGCACTGGGGAA AGGCAGAGGGAGCCCCGAGCCGGCCGAGGGGTTGGAGGCCCTGGCCCTATCAGTCTGCCATCCCCACCC CATGGCTGTAATGTCTTGTATTTTTAAATAAAGAGATATTGATGTCTTGTGCTCACTGAGGCATC TAAGAG
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000242
Synonyms:	MIRN199A1
Locus ID:	406976



[View online »](#)