

Product datasheet for **SC400257**

MIR200B Human MicroRNA Expression Plasmid (MI0000342)

Product data:

| | |
|---------------------------|--|
| Product Type: | MicroRNA Expression Plasmids |
| Species: | Human |
| Vector: | pCMV-MIR (PCMVMIR) |
| Precursor: | CCAGCUCGGGCAGCCGUGGCCAUCUUACUGGGCAGCAUUGGAUGGAGUCAGGUCUCUAAUACUGCCUGGUAUAUGAUGACGGCGGAGCCCGCAGC |
| Mature Sequence: | CAUCUUACUGGGCAGCAUUGGA |
| Full insert: | >MI0000342 AGGGAAGTCCCCCGCAGAGGGGTGGGGCAGAGGGCCGAGCGGGCGGGCAGAGGGCCCGTGTGAGCCC CACTCCGACCTAGTCTCGGCCGTCTGGCCAGGACTTTCGGCCCCCAGGTGCCACCCAGGACCCAA AGCTGGTGGCTGCTGGACTCGGCAGGGCTGGCGGGTGGGGCTACCCGGGCCCTGCCCTCCGGCGATGC TGTCCTCAGTGCCCCAGGAGGACGAGGCCCCAGCTACTGAGCTTCCAGCGAGTCCCATGCAACCCTC AGCCGGGCGGGCCCCGAGCCAGCTCGGGCAGCCGTGGCCATCTTACTGGGCAGCATTGGATGGAGTCAG GTCTCTAATACTGCCTGGTAATGATGACGGCGGAGCCCTGCACGCAGCGACCGGCCGACCCCGTCCC CCCCAGGGCTCCCGCGAGCCACACTCCCCTCCCACAGCCCGCTACCCGACCCACCCCGTCCG GGCCAGATGGCTCTGGGGAGAGTGACGGGGGACGTGACCTTCGGAAGGGGGCTGGGGCGCCCCACCA GCATCTCAGAGACCCCGGAATGCAAGGCCTCCTGCAGCTGCACCCAGGGCGCCCCACAGTGATCTGGGGA TTAGGACGCTCAGGTGTCCGGGCCCTGCCACAGCGCCTGGGCGGGAGCTGCTGTACACCCA |
| OTI Disclaimer: | All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA. |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Mammalian Cell Selection: | Neomycin |
| Components: | 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube |
| Also Contains: | Both mature and minor (star) miRNA can be over expressed from the construct. |
| RefSeq: | MI0000342 |
| Synonyms: | MIRN200B |
| Locus ID: | 406984 |



[View online »](#)