

Product datasheet for SC400273

MIR211 Human MicroRNA Expression Plasmid (MI0000287)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	UCACCUGCCAUGUGACUUGUGGGCUUCCCUUUGUCAUCCUUCGCCUAGGGCUCUGAGCAGGGCAGGGACAGCAAAGGGGUGCUCAGUUGUCACUCCACAGCA
Mature Sequence:	UUCCCUUUGUCAUCCUUCGCCU
Full insert:	>MI0000287 TACTCCAGTAATTCTAACGTCAAACAATTGTAAATGATCTGAAAGAGCAGTTTGTTCCTCCATTACT GGCAAGATTTTGTGGGAGACCATCTAGTCTGGCGAATGGTCATTGGGAATCATCATGTATGGGACTGT GATGCTGCAGAGTGGGTAGAAATGGGGATGGCCTACTAGGAAAAGTTGTAGATTTCCGAGAATCAGATCTC GAAGAATACATTGGTCGATGACTGAATGTTTTTGTCTCACATCTAATTTTTCTGGACTTGAAGTT TTTACAACACCCCATTTACCTGGCCATGTGACTTGTGGCTCCCTTTGTCATCCTTCGCCTAGGGCTC TGAGCAGGGCAGGGACAGCAAAGGGGTGCTCAGTTGTCACTCCACAGCAGGAGCTTCGAGGACTCCA GCCTAAGCCAAGAAGACCTGTTCTACTCTGTTGCTCGAAAAGAAGGTGATTTGGCCACCATATTCCTAA AGATTGAAAATATCAATCCTTGCTTGGGTTGGAAATCTCAAAGGCGGAGTCCAGGATATCTTAGTTCATT TTTCCAATATCTCCTGGTCAGGGCTTGCACAAAGCAGACCCTCAATGCATGTCTGCTGAAATGAACTTA GCTGTAATACTCTGGTGGACACAGCAGAGTTATTCAACATATGTAAGAGTCCCTGACTGTAGAAACACC TC
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence imiRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000287
Synonyms:	MIRN211 mir-211



[View online »](#)

Locus ID: 406993