

Product datasheet for **SC400285**

MIR219-2 Human MicroRNA Expression Plasmid (MI0000740)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Human

Vector: pCMV-MIR (PCMVIR)

Precursor: ACUCAGGGGCUUCGCCACUGAUUGUCCAAACGCAAUUCUUGUACGAGUCUGCGGCCAACCGAGAAUUGUGGCUGGACAUCUGUGGCUGAGCUCGGG

Mature Sequence: UGAUUGUCCAAACGCAAUUCU

Full insert: >MI0000740 CAGGCCAGACAGAGGAGATTCCCCTCGACCTGGGGCGGCGCTGAGACCAACGCGGGCCACTCTGCGGA
CTCCGGGGCGGATCGAAGGAGCCACGCCAGGGGCGATGGGGATCCCCGCACAGGGTCAGGATCCGCGTG
CTTGCTAGTCCACCTTACCGCCGATGGGCAGCCGCTGCTGAGCGCTGCTAGGTGGTCCGACGCGCAGC
TCCGAGCTGTCTGGACTGGCGCTTTTATGAGAGCGAAGGGGAACCCGCGAGGAGACCAGGGCCCTGAAC
CAGGGGCTTCGCCACTGATTGTCCAAACGCAATTCTTGTACGAGTCTGCGGCCAACCGAGAAATTGTGGCT
GGACATCTGTGGCTGAGCTCCGGGCGCAACAGGGGCGGGGGCCCCAGGGACAGGGCTCAGCGCGGGCGAG
ACCTCTCGGGGCGGCGGCTGGCTGGCCCCAGCGCGAGGATCGCGGTCCCGGCCGCGCGCACAGAGACGC
CGGTCCCTGCCACCAGCGCCGCCATCCGCATCCTCGTCTTCTTCTTCTCTGCGGGTGACCAGTCCCGG
CCGGCAGGGTCCTTCAGTCGATCGCGTCTTTTCTCTCTTCAAATTCGGCCGGAACGAAGCCGGCAC
GCAGCGCGGGGCCGAATGTGGTCCGCCCTCA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0000740](#)

Synonyms: MIRN219-2

Locus ID: 407003



[View online »](#)